

# Armstrong Medical

INDUSTRIES INC.

2017R



REVISED EDITION



#141R  
Small Business

800/323-4220 • [www.armstrongmedical.com](http://www.armstrongmedical.com)

# Keyless-Entry Auto-Locking Cart

With Optional Proximity Reader

## A-SMART®

*Auto-Locking Mini Cart w/optional Proximity Reader, pg. 103*

*Auto-Locking Cart w/optional Proximity Reader, pgs. 80-81*

*Auto-Locking Narrow Cart w/optional Proximity Reader, pg. 103*

**WIRELESS**  
Auto-Locking Cart  
w/optional  
Proximity Reader and  
**CONTROLLED  
SUBSTANCE  
DRAWER,**  
pgs. 77-79





# Table Of Contents

## Training, Pages 4-32

Actar CPR Manikins/AED Trainer. . .	4-7
CPR Manikins. . . . .	8-10
Laerdal CPR Manikins. . . . .	9-10
CPR Adjuncts. . . . .	11
Defib/Code Manikins. . . . .	12-16
Airway Management Trainers. . .	15-20
Vascular Access Manikins. . . . .	20-24
Procedure Simulators. . . . .	22-24
Trauma/Rescue Manikins. . . . .	25
X-Ray Positioning Manikin. . . . .	26
Hospital Training Manikins. . . . .	27-29
O.B. Manikins. . . . .	30-31
Exam/Procedure Trainers. . . . .	32

## Personal Protection, Pages 33-34

## Equipment, Pages 35-69

Anesthesia Chairs. . . . .	35-37
Glove Boxes and Security Box. . . . .	38
Suction Units. . . . .	39-40
Broselow® System. . . . .	41-45
Laryngoscopes. . . . .	46-48
Diagnostic Equipment. . . . .	49-53

Stethoscopes. . . . .	54
Privacy Screens. . . . .	55
Scales. . . . .	56
Hospital Equipment. . . . .	57
Resuscitators & Pressure Infusers. . .	58
Splints. . . . .	59-61
Collars, Boards & Immobilizers. . .	62-65
Eye Exam Treatment. . . . .	66
EMS Supplies & Cases. . . . .	67-69

## Cart Systems, Pages 2-3, 70-124

Cart Systems Table of Contents. . .	70
Benefits of A-SMART Carts. . . . .	71
Introduction/How To Order:	
A-SMART® Standard/Premier™. . .	72-75
Build Your Own Cart. . . . .	76
<b>Wireless Auto-Locking Carts. . . . .</b>	<b>77-79</b>
Auto-Locking Carts. . . . .	80-81
Bedside Carts. . . . .	82-83
Anesth./Mobile Workstations. . .	84-91
Emergency Carts. . . . .	92-99
Push-Button Carts. . . . .	100
Electronic Touchpad Carts. . .	101-102
Wide Carts. . . . .	102

ScopeSafe. . . . .	102
Mini and Narrow Carts. . . . .	103-105
Broselow Pediatric Emergency ColorCode Carts. . . . .	106-107
Premier Locking Bars. . . . .	108
Isolation Carts. . . . .	109
Cart Extender . . . . .	110
Glove Boxes . . . . .	110
Secure Storage. . . . .	111-112
Trellis Systems. . . . .	112
Cart Accessories. . . . .	113-121
Stainless Steel/Utility Carts. . .	122-123
Wire Shelving. . . . .	124

## Index, Pages 125-126

## Ordering Information, Pages 127-128

## National Sales Force Located To Serve You

with  
Pre-Sale Support • Delivery Support • Post-Sale Support



## Armstrong Medical INDUSTRIES INC.

★ Corporate Headquarters

575 Knightsbridge Parkway • P.O. Box 700  
Lincolnshire, IL 60069-0700

800/323-4220 • Local: 847/913-0101

FAX: 847/913-0138

[www.armstrongmedical.com](http://www.armstrongmedical.com)

[csr@armstrongmedical.com](mailto:csr@armstrongmedical.com)

Prices and products are subject to change without notice.

# Organize Your EMERGENCY DEPT

*With These Essentials*

*Auto-Locking  
Narrow Cart  
w/optional  
Proximity Reader,  
pg. 103*



**NEW** Aluminum  
Electronic Touchpad  
Supply Carts,  
see pages 101-102



Premier™ Aluminum Solid  
Color Emergency Carts,  
see page 97



*DUET Suction Unit  
see page 39*



*Broselow Pediatric Resuscitation System  
see pages 41-45*



*Actar D-fib, Actar AED Trainer,  
Actar Original, Actar Infantry  
see pages 4-7*



# A-SMART® CARTS + AREA SALES REP = SATISFIED CUSTOMER A-SMART ADVANTAGE

## CONSIDER THESE VALUE ADDED SERVICES

# Before You Make A Longterm Investment

**ASK YOUR LOCAL REPRESENTATIVE FOR A FREE DEMO OF OUR MOST POPULAR CARTS**

### *Armstrong Medical Offers:*

- Nationwide coverage (over 40 Sales Representatives)
- Pre-sale support
- Delivery support
- Post-sale support
- **The Best Overall Value Around**

### FEATURES OF A-SMART® PREMIER™ ALUMINUM CART

- **ALUMINUM** Carts are 45% lighter than Steel
- **Full Extension Ball Bearing Drawer Slides**  
*Ensures that your drawers will open and close smoothly every time.*
- **Quality Aluminum Construction**  
*One-piece outer shell construction providing you with added strength and stability. Locking components consist of zinc plated and galvanized parts.*
- **Double Side-Wall Construction**  
*Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.*
- **Two Interchangeable Locking Systems**
- **Customize Your Locking Configuration**  
*Lock individual drawers, double lock drawers, use multiple types of locks on the same cart.*
- **Slide-Thru Side Shelf**
- **Contour-Fit Plastic Top**
- **Four Vertical Aluminum Corner Mounting Tracks**  
*For easy accessory attachment (eliminates the need for holes).*
- **Stabilizing Frame With Bumper**  
*Casters are moved out from underneath to the outside of the cart. This extra wide base provides your cart with added stability.*
- **18 Colors Available**  
*Allows you flexibility in your carts appearance. Choose colors to specify department or function, or even use multiple colors on the same cart to differentiate drawers.*
- **Highest Quality Swivel Casters (two locking, one tracking)**  
*Quiet, durable, German-made casters. Smoother, easier maneuvering.*

**A-SMART® Carts Are Made In The USA**

**Armstrong Medical** INDUSTRIES INC.

575 Knightsbridge Parkway • P.O. Box 700 Lincolnshire, IL 60069-0700  
800/323-4220 • 847/913-0101 • FAX: 847/913-0138  
[www.armstrongmedical.com](http://www.armstrongmedical.com) [csr@armstrongmedical.com](mailto:csr@armstrongmedical.com)

# The Original **ONE-TO-ONE** CPR Manikin

Durability Tested For Over 500,000 Compressions

## CPR Teaching Features:

- **JAW THRUST, HEAD TILT, AND CLOSEABLE AIRWAY**
- All landmarks: Adam's apple, carotid groove, lower rib, nipples, and navel
- Chest compressions
- Visible chest rise and fall
- Adult and child CPR
- Abdominal thrusts
- No cleaning of internal parts necessary (lung/mouthshield)
- One manikin per student
  - More practice time maximizes psychomotor skills
  - Reduces the risk of cross contamination
  - Maximizes class time efficiency
- Latex-free
- Bag valve mask compatible

## AED Training Features:

- Electrodes are easy to remove
- Compatible with all brands of training electrodes
- All landmarks: lower rib, nipples, and collar bone
- Electrodes don't interfere with chest rise or chest compressions

AA-1305 Actar D-fib™ 5 Pack with Nylon Back Pack and 50 Lung/Mouthshields, \$375.00

AA-1310 Actar D-fib™ 10 Pack with Nylon Carry Bag and 50 Lung/Mouthshields, \$675.00

AA-1325 Actar D-fib™ 25 Pack with Rolling Cart, Nylon Cover, and 100 Lung/Mouthshields, \$1,595.00

AA-1312 Actar D-fib™ Lung/Mouthshields, 50/pkg., \$22.00



30 + 2 CD, sets the pace for at least 100 beats per minute.



Ventilation



Immobilization and Jaw Thrust



Bag Mask Compatible



Head Tilt



AED Training



Chest Compression



AA-1312

AA-1305 Actar D-fib™ 5 Pack



AA-1310 Actar D-fib™ 10 Pack

**Great for Schools!**

AA-1325 Actar D-fib™ 25 Pack







Actar D-fib™ Manikin sold separately. See page 4.

## Realistic AED Sequencing

## CPR/AED GUIDELINES ON CD



Actar D-fib™ Manikin sold separately. See page 4.

## More Hands-On Training

# Actar AED Training System

The Actar AED Trainer™ 5 and 10 Packs allow multiple students to practice AED training scenarios simultaneously. This not only reduces class time but also gives the students more hands-on training, which leads to greater skill proficiency and confidence. The students simply follow the prompts from the audio CDs or a conventional AED trainer which gives the instructor more opportunity to monitor and assist the students.

The Actar AED Trainer provides realistic AED sequencing which includes electrode pad placement, power on button, analyze button, and shock button. The Actar AED Trainer comes with two audio CDs. One CD contains scenarios for lay responders, the other contains scenarios for healthcare providers. Actar AED Training Electrodes are reusable and won't leave a residue on your manikin. Each pair of electrodes comes with a simulated "nitro" patch.

### 1. AED TRAINER (set of 5 or 10)

Each trainer includes three buttons that simulate "on/off" (green), "analyze" (yellow), and "shock" (red), and a socket for the electrode cable. Units stack for efficient storage.

### 2. TRAINING ELECTRODES

Trainer 5 Pack includes twenty-five pairs of training electrodes, 10 Pack includes fifty pairs. Each pair includes a "nitro" patch (simulation) for use with scenario 5 (CD track 5). The training electrodes will not function on an actual AED.

### 3. CABLES (set of 5 or 10)

Trainer 5 Pack includes five cables, 10 Pack includes ten cables. Each cable includes two connectors that snap securely to the electrode pads and a plug that connects to the AED trainer. Velcro ties and zippered pouch are provided for tidy storage.

### 4. AUDIO CDs - AED Training Scenarios

Trainer 5 and 10 Packs include two audio CDs. One CD contains scenarios for lay responders, the other contains scenarios for healthcare providers. All scenarios are produced in accordance with the 2010 Guidelines for CPR and ECC.

### 5. PADDED CASE (for 5 or 10 packs)

Stores components in a protective case which can be attached to the Actar D-fib™ 5 or 10 Pack carrying cases.

**AA-1355 Actar AED Trainer 5 Pack, 1 Electrode Package (25/pkg.), 5 Electrode Cables, 2 Audio CDs, and Padded Case, \$350.00**

**AA-1360 Actar AED Trainer 10 Pack, 2 Electrode Packages (25/pkg.), 10 Electrode Cables, 2 Audio CDs, and Padded Case, \$595.00**

**AA-1380 Electrodes, 4 packages, 25 pairs/pkg., \$180.00**



Visit Our Website At  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

# Actar 911™

## The Original **ONE-TO-ONE** CPR Manikin Introduced In 1990 (Millions of Students Taught)

### *Patrol Plus*



#### AA-1840 Actar 911 Patrol™ Plus (5 Pack), \$465.00

Cadet piston converts an adult Actar manikin into a child manikin and it doubles to assist abdominal thrust practice.

Patrol Plus includes five adult manikins, five cadet pistons, one hundred disposable lungs, five triple-action airways, instructions, video, nylon Patrol

bag with shoulder strap, and Cadet Pac bag. Weighs only 13 lbs.

#### AA-1825 Actar 911 Patrol™ (5 Pack), \$365.00 (not pictured)

Includes five adult manikins, one hundred disposable lungs, five triple-action airways, instructions, video, and nylon Patrol bag with shoulder strap. Weighs only 11 lbs.



### *Squadron Plus*



#### AA-1845 Actar 911 Squadron™ Plus (10 Pack), \$775.00

Each cadet piston converts an adult Actar manikin to a child manikin and doubles to assist abdominal thrusts. Includes ten adult manikins, ten cadet pistons, one hundred disposable lungs, ten triple-action airways, video, instructions, nylon Squadron bag with shoulder strap, and two Cadet Pac bags. Weighs 21 lbs.

#### AA-1830 Actar 911 Squadron™ (10 Pack), \$575.00 (not pictured)

Includes ten adult manikins, one hundred disposable lungs, ten triple-action airways, video, instructions, and nylon Squadron bag with shoulder strap. Weighs only 17 lbs.



#### AA-1810 Actar™ Heads, 10/pkg., \$79.50



### Actar 911 Cadet Pac™ and Accessories



#### AA-1835 Actar Cadet Pac, \$100.00

The innovative Cadet Compression Piston converts an adult Actar manikin to a child version in seconds. Simulates child chest resistance with 1.5" max. compression depth. Piston doubles to assist abdominal thrusts. Cadet Pac includes five cadet pistons, instructions, bag, and shoulder strap. Weighs 2 lbs.


AA-1812 Disposable Lungs for Actar™ Manikins, 100/pkg., \$35.00; 2 or more packages, \$25.00/pkg.

AA-1813 Actar™ Adult Chest Plate, \$25.50

AA-1814 Actar™ Adult Piston, \$25.50

AA-1817 Actar™ Triple-Action Airways, 10/pkg., \$17.50





## Infantry

Infant CPR Training System

### AA-1860 Actar 911 Infantry™ 10 Pack, \$620.00

Includes ten infant manikins, one hundred disposable lungs, ten adapter caps, video, instructions, and nylon Infantry bag. Weighs just 12 lbs.

### AA-1855 Actar 911 Infantry™ 5 Pack, \$335.00 (not pictured)

Includes five infant manikins, one hundred disposable lungs, five adapter caps, video, instructions, and nylon Infantry bag. Weighs only 7 lbs.

### AA-1871 Actar 911 Infantry™ Chest Plate, \$13.50

AA-1862 Infantry Lungs, 100/pkg.,  
\$35.00; 2 or more, \$25.00/pkg.

AA-1870 Infantry Heads, 10/pkg.,  
\$76.00

AA-1865 Infantry Adapter Caps,  
10/pkg., \$19.00



AA-1860



AA-1862

Visit Our Website At  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

## Little Joe™

### Adult Manikin

#### Extremely Sanitary

Prevents cross contamination with an individual face piece for each student and a disposable, non-rebreathing airway.

#### Lifelike Features

Head tilt/chin lift mobility. Plus, all essential landmarks, including navel for obstructed airway practice. Longer torso for abdominal thrusts. Realistic rib cage overlay for landmarking.

#### Quick Clean Up

After class, the inexpensive face pieces can be discarded, or sanitized and reused. Airways are disposable.

#### Durable and Portable

Built to withstand high-volume usage. Easy to transport.

#### Double Your Training Capacity

AA-5872 Little Joe™ 2 Pack, Caucasian, \$430.00

AA-5872B Little Joe 2 Pack, African-American, \$430.00

AA-5872E Little Joe 2 Pack, Hispanic, \$430.00

Includes two manikins, six face pieces, six airways, and dual carrying bag.

AA-5875 Little Joe 1 Pack, Caucasian, \$222.00

AA-5875B Little Joe 1 Pack, African-American, \$222.00

AA-5875E Little Joe 1 Pack, Hispanic, \$222.00

Includes manikin, three face pieces, three airways, and single carrying case.

AA-5615 Face Pieces, Caucasian, 10/pkg., \$52.10


AA-5615B Face Pieces, African-American, 10/pkg., \$52.10

AA-5615E Face Pieces, Hispanic, 10/pkg., \$52.10

AA-5880 Airways, 24/pkg., \$65.30

AA-5803 Dual Carrying Bag, \$40.00  
2526 Single Carrying Case, \$59.20

AA-5872



AA-4200

AA-4200B

### CPR Kyle™ Three-Year-Old Child Manikin

AA-4200 CPR Kyle, Caucasian, \$344.00

AA-4200B CPR Kyle, African-American, \$344.00

AA-4205 CPR Kyle with **Jaw Thrust**, Caucasian, \$372.00

AA-4210 Airway Systems, 24/pkg., \$65.30

AA-5215 Face Pieces, Caucasian, 10/pkg., \$52.10

AA-5215B Face Pieces, African-American, 10/pkg., \$52.10



Manikin Face Shields, see pg. 11



## Chris Baby™ CPR Manikin

AA-1100 Chris Baby™, \$425.00

Includes five face pieces, one hundred air head bags, and soft carrying case.

*Eliminate cross contamination.*



AA-1100

### Realistic Anatomy

- Ventilation requires proper head tilt to open the airway
- Exhaled air from the mouth and nose creates realism in training
- Chest rises with proper ventilations
- Stomach clearly shows inflation
- CPR or bag mask ventilation is possible
- Rib cage has a realistic feel
- Lever allows complete airway obstruction

### Minimal Time Spent On Clean Up and Disinfection

- All parts in contact with the trainee's exhalations are for individual use only
- Face pieces are easily disinfected
- Air head bags are single-use and disposable
- After class, Chris Baby cleans up in a few minutes – no disassembly or inside cleaning required

AA-1160 Face Pieces, 5/pkg., \$55.00

AA-1170 Air Head Bags, 100/pkg., \$32.50

## Baby Anne™ 4 Pack

AA-2410 Baby Anne 4 Pack, \$445.00

The Baby Anne 4 Pack is a low-cost alternative for infant CPR training from Laerdal.

Includes twenty-four airways, forty foreign body practice objects, case, and instructions.

- Weighs only 12 lbs. total – easy to transport
- Manikin is 22" long
- Includes unique airway obstruction feature

- Requires head tilt/chin lift to ventilate
- Provides realistic chest compression and lung inflation

AA-2406 Baby Anne Face Masks, 6/pkg., \$61.20

AA-2424 Baby Anne Airways, 24/pkg., \$18.80

05 04 00 Replacement Jaw, \$8.40

05 05 00 Mask Connector, \$7.20



AA-2410

## Laerdal™ Little Junior™ Child CPR Training Manikin

AA-2600 Laerdal Little Junior Child CPR Training Manikin with Carrying Bag, \$215.00

Includes two faces, two airways, six manikin wipes, soft pack, and

directions. Soft pack doubles as training mat.

AA-2610 Little Junior 4 Pack, \$776.00

4 Pack includes eight faces, eight airways, twenty-four manikin wipes, soft carrying bag, and directions.

Continuing in the tradition of Little Anne™ and Baby Anne™, Little Junior offers:

- Anatomical landmarks and realistic features, like natural airway obstruction and articulating jaw

- Realistic compliance for compressions and ventilations
- Realistic chest rise and fall
- Quick clean up with low-cost airways
- Airways with one-way valve prevents cross contamination

AA-2606 Face Pieces, 6/pkg., \$61.20

AA-2625 Airways, 25/pkg., \$55.10

AA-2630 Airways, 100/pkg., \$204.00



AA-2600

# Laerdal® CPR Manikins

## Laerdal Little Family Pack

AA-1700 Little Family Pack, Light Skin, \$556.00

AA-1700B Little Family Pack, Dark Skin, \$556.00

Little Family Packs come with one Little Anne™ CPR Manikin, one Little Junior™ CPR Manikin, one Baby Anne™ CPR Manikin, and a 3-pack carrying case.

### Little Anne features:

- Oral and nasal passages
- Natural obstruction of airway
- Head tilt/chin lift and jaw thrust
- Anatomically correct landmarks and sternal notch
- Realistic resistance for chest compressions
- Disposable non-rebreathing airway
- Audible feedback for chest compressions

### Little Junior features:

- Oral and nasal passages
- Natural obstruction of airway
- Head tilt/chin lift and jaw thrust
- Realistic airway
- Anatomically correct landmarks and sternal notch
- Realistic resistance for chest compressions
- Audible feedback for chest compressions
- Hygienic convenience; individual face mask for each student, face can be easily cleaned between students
- Non-rebreathing airway
- Disposable airway changed after class



AA-1700

### Baby Anne features:

- Natural obstruction of airway
- Head tilt/chin lift and jaw thrust
- Foreign body airway
- Realistic resistance for chest compression

**Child CPR, pages 7-9**

## Resusci Anne® QCPR

AA-1650 Resusci Anne® QCPR Torso, \$1,580.00

AA-1660 Resusci Anne® QCPR Full Body, \$2,061.00

AA-2050 Resusci Anne® Baby QCPR, \$939.00

AA-1670 SkillGuide, \$64.25



AA-1006 Face Pieces, 6/pkg., \$61.20

AA-1524 Airways, 24/pkg., \$116.00

Laerdal's QCPR Technology enables this manikin and feedback devices to measure the quality of CPR, providing real-time and summative feedback on compression rate, depth, release, hands-off time, and other components of high-quality CPR as defined in the American Heart Association Consensus Statement.

- Realistic anatomy including head tilt, chin lift, compression depth, compression force, and chest rise
- Sensor indicates correct hand placement
- Ventilation system provides appropriate chest rise with BVM (Bag Valve Mask) and MTM (Mouth to Mouth)
- Enhanced feedback with *optional* SkillGuide





## Manikin Face Shields

AA-1880 CPR Manikin Face Shields, 6 packs of 36, \$65.00



## Resusci® Manikin Face Shields

15 12 01 Resusci Manikin Face Shields, 6 Rolls of 36, \$67.30



## AE-5307 CPR Training Mouthpieces, 10/pkg., \$15.00

Includes filter. Standard diameter fitting allows use with most CPR masks. Valves are single use and for manikin training only.



## AE-550 CPR Backboard, \$100.00

High-density, polyethylene backboard offers firm surface for administering CPR. 16"W x 23 1/2"L.



## AA-8300 Manikin Storage Soft Case, \$155.00

Heavy-gauge, soft-lined Cordura fabric. Dual inside straps hold manikin secure for transport or storage. 69"L x 9 1/2"W x 20"D.

*Manikin not included.*



## AA-6090 Kneeling Pad, 6/box, \$100.00

Poly foam pad for CPR practice. 12" x 24" x 1 1/2".



## Anti Choking Trainer

AA-6080 ActFast™ Anti Choking Trainer, Red \$139.00

AA-6082 ActFast™ Anti Choking Trainer, Red, 4/pkg., \$429.00

AA-6085 ActFast™ Anti Choking Trainer, Blue (no back pad), \$129.00

AA-6087 ActFast™ Anti Choking Trainer, Blue (no back pad), 4/pkg., \$429.00

The training vest can be used in a variety of positions to create realistic scenarios, including: standing, sitting and laying on the ground. Also, it allows the student to practice self rescue, using a counter top or chair.

- Transparent PVC airway allows user to see release of obstruction.
- Large easy snap buckles with adjustable straps. One size fits all.
- Back Pad (Red only) crumple zone for practicing back blows.
- Navel aids in correct hand placement.

# Defib Simulators and Crisis Manikin

## 12-Lead Arrhythmia Simulator

AA-601 12-Lead Simulator, Zoll, \$675.00

AA-602 12-Lead Simulator, Physio "Quick Combo", \$675.00

AA-604 12-Lead Simulator, Philips MRx, \$675.00

Simulate electrical capture with your external pacer. Press simulator capture key to select one of four preset pacing capture levels: 70, 80, 90, or 100 mA.



AA-602

When the pacer current is greater than the selected capture level, paced beats appear on your monitor. Waveforms for pacing include sinus brady (two), 1st degree AV block, 2nd degree type I AV block, 2nd degree type II AV block (with and without PVCs), and 3rd degree AV block.

Simulate cardioversion with your manual, semi-automatic, or automatic defibrillator. Before defib discharge, activate "convert" feature to select another rhythm to run immediately after defib discharge. Waveforms for defibrillator training include V-Fib, V-Tach (high rate and low rate), Torsade, A-Fib, A-Flutter, PSVT, Sinus Tach, Sinus

Rhythm, Sinus Rhythm with PVCs, Asystole, and NSR.

Generate realistic 12-lead ECG rhythms. Connect your ECG cable to the simulator's ten ECG snaps. Independent chest lead and limb lead ECG signals create realistic 12-lead ECGs for each rhythm. Generate ST segment and T wave abnormalities, including Anterior MI, Inferior MI, Antero-Septal ST Elevation, Anterior ST Depression, Lateral ST Elevation, and Inferior ST Elevation. LED indicators include pacer pulse detection, defib discharge detection, and low battery.

Three-year warranty. Weighs 2 lbs.

## "Hands-Free" Defib Pad Adapters for Defib and CRiSis™ Manikins

For use with CRiSis Manikins and the Complete Defib Training Package.

AA-6030 Physio Control "Quik Combo"™, \$225.00

AA-6021 Zoll, \$225.00

## Interactive ECG/Pacing Simulator

AA-550 Interactive ECG/Pacing Simulator, \$750.00

This simulator combines defibrillator and external pacer training features with adult and pediatric ECG rhythms. Built-in circuitry allows users to practice defibrillation and external pacing with or without a manikin\*. The simulator features seventeen adult and pediatric rhythms for pacing on any defibrillator. It has the ability to automatically convert

from one rhythm to another when defibrillated. Uses one 9V battery.

\* Extra adapter may be required.


Z02 Zoll Defib Pad Adapter, \$150.00

Z10 Physio Quick Combo Defib Pad Adapter, \$150.00

Z13 Philips MRx Defib Pad Adapter, \$150.00



AA-550



## CRiSis™ Patient Manikin

AB-1700 CRiSis Patient Manikin, Complete System, \$3,647.00

- CRiSis Manikin with Airway Larry® airway management head
- I.V. training arm
- New and improved blood pressure arm
- Defibrillation chest skin

AB-1730 CRiSis Manikin with Heart and Lung Sounds,

ECG/Pacing Simulator, and Case, \$6,262.00

**Heart sounds** include normal, holosystolic, systolic click, aortic regurgitation, midsystolic, atrial septal defect, mitral stenosis, S3 gallup, PDA, pulmonary stenosis, S4 gallup, and VSD. **Lung sounds** include normal,

rhonchi crackle, bronchial, egophony, wheeze, stridor, pulmonary edema, mono wheeze, coarse crackle, friction rub, fine crackle, cavernous, and pectoriloquy.





AB-1930

Four-Lead ECG Monitoring

## Infant CRiSis™ Manikin

AB-1930 Infant CRiSis Manikin with ECG Chest Cover, \$1,704.00

The most comprehensive PALS trainer available, the Infant CRiSis Manikin is designed for practice of I.V. therapy, umbilical catheterization, airway management, intraosseous infusion, and CPR. Comes equipped with four-lead ECG monitoring. Teach infant life saving techniques to individuals and groups. Tactile and visual realism provides students with excellent training.

Manikin includes:

- Airway management trainer head with intubation
- Intraosseous infusion simulator
- I.V. arm
- I.V. leg
- ECG chest cover
- Umbilical cannulation chest skin
- Hard carrying case

AB-1925 Infant Crisis Manikin (w/o ECG chest cover), \$1,428.00

## Pediatric ALS Trainer

AA-9300 Pediatric ALS Trainer, \$912.00

AA-9310 Pediatric ALS Trainer with ECG/Pacing Simulator (AA-550, see pg. 14), \$1,647.00

AA-9301 Replacement Bones, 12/pkg., \$109.20

AA-9302 Replacement I/O Leg Skins, 4/pkg., \$54.10

- Bag valve mask ventilation
- Oral and nasal intubation
- Three-lead monitoring – can be used with most ECG simulators
- Brachial pulse

- External chest compressions
- I.V. arm and hand
- Intraosseous access (both legs)
- Scalp vein for palpation

- Cranial sutures and fontanel
- NG feeding tube
- Sellick maneuver
- Jaw thrust



AA-9300

## Laerdal™ ALS Baby™

AA-2800 Laerdal ALS Baby without Heartsim 200, \$1,795.00

Includes carrying case, five replacement leg pads, pulse bulb, and directions. Represents 1-3 month old.

AA-2805 Laerdal ALS Baby with Heartsim 200, \$3,135.00

Includes Heartsim 200, carrying case, five replacement leg pads, pulse bulb, and directions. Represents 1-3 month old.

AA-2820 Replacement Leg Pads, 5/pkg., \$237.00

08 33 10 Replacement Face Skin, \$170.00



AA-2800

### Professional Rescuer CPR

- Practice bag valve mask ventilation
- Practice insertion of oral/nasal airways
  - Practice external chest compressions
  - Simulate brachial pulse

### Vascular Access

- Practice intraosseous needle insertion with aspiration of simulated bone marrow

### Airway Management

- Teach and practice bag-valve mask ventilation oral and nasal intubation, gastric tube insertion, Sellick maneuver, and LMA insertion

### ECG Monitoring

- Three-lead monitoring when used with optional simulator

# STAT Manikin



## STAT Manikin

AA-7300 STAT Manikin with Carrying Case, \$5,415.00

The STAT Manikin's body demonstrates a giant step forward in manikin design. The joints' range of motion is near human. Rib cage motion mimics lifelike respirations. Multiple pulse points demonstrate blood pressure variation.

### Features:

- Twelve pulse points
- Accurate human anatomy
- CPR
- Multiple injection sites
- Advanced I.V. arm
- Oral and nasal intubation
- Four-lead monitoring
- Tongue swells
- Laryngospasms
- Cricothyrotomy
- Bilateral chest tube placement with simulated fluid discharge
- Bilateral pneumothorax relief
- Stomach distention
- I.V. arm therapy
- Defibrillation
- Arrhythmia training

The STAT Manikin comes with the AA-550 Interactive ECG/Pacing Simulator (see page 12).


Shipping size: 57" x 34" x 12"  
Shipping weight: 100 lbs.

- Z02 Zoll Defib Pad Adapter, \$150.00
- Z10 Physio Quick Combo Defib Pad Adapter, \$150.00
- Z13 Philips Defib Pad Adapter, \$150.00




AA-7315 Optional STAT Manikin Intraosseous Leg, \$246.90

For more details, see AA-7310 Intraosseous Leg with Stand on page 23.




Accepts All Airways



Advanced I.V. Arm



Cardiac Care



Bilateral Chest Tube Placement



Bilateral Pneumothorax Relief



# AED and Airway Trainers

## Laerdal™ AED Trainer 2 Manikin Package

- AA-1942N AED Trainer 2 with Manikin, \$560.00
- AA-1940N AED Trainer 2 with Case only, \$465.00
- AA-5910 Lungs, 100/pkg., \$40.00

The AED (Automated External Defibrillator) Trainer 2 makes complex defib training simple. The unit mimics the Forerunner® 2 and has ten preset comprehensive training scenarios available at the touch of a button.

AED Trainer 2 features:

- User interface
- Rescue voice prompt

- AA-1941 Replacement Training Electrodes, \$26.85
- AA-1945 Pediatric Training Electrodes, \$52.55
- AA-1943 Optional Remote Control, \$53.60

- Adult training electrode pads

Manikin features:

- AED pad placement practice
- CPR practice
- Single-use airway/lung/face shield system
- Anatomical landmarks: sternum, rib cage, and substernal notch



AA-1942N

## Laerdal™ AED Trainer 3


- AA-1960 Laerdal AED Trainer 3, \$363.00
- AA-1961 Replacement Pads, \$31.10
- AA-1963 Optional Remote Control, \$41.80

The AED Trainer 3 reflects the unique functionality of the Philips HeartStart FR3 AED. It is configured with eight real-world scenarios that have been developed in accordance with internationally recognized responder programs. They are AHA and ECC 2010 Guidelines

compliant. Many parameters on each scenario can be easily configured to conform to local protocols or future Guidelines updates.

The instructor can change scenarios by pressing the buttons on the device or by using an optical remote control.

AA-1960  
(shown with optional accessories)



AA-3900

## Critical Care Airway Trainer

- AA-3900 Critical Care Airway Trainer, \$1,310.00
- AA-3910 Replacement Lungs/Stomach, 3/pkg., \$16.40

This very difficult intubation head, coupled with the CPR torso, supplies the challenges for initiating respiratory care to the critical patient. By using the features of tongue swelling and laryngospasms, the instructor creates a

critical airway management scenario that requires the student to perform cricothyrotomy. In addition, two external and one internal bleeding wounds demand attention with suctioning and bleeding control.

## Patient Education Tracheostomy Care Set

- AB-1136 Patient Education Tracheostomy Care Set, \$878.00
- AB-1137 Adult Manikin only with Carrying Case, \$648.00
- AB-1138 Infant Manikin only with Carrying Case, \$295.00

The Adult Manikin includes oral and nasal passages, the pharynx, epiglottis, trachea, esophagus, tracheostomy site, and cricoid cartilage. Students can practice suctioning techniques, cuff inflation, dressing changes, and other techniques. The Infant Manikin can be


used to teach parents and caregivers basic tracheostomy care skills.

AB-1136 includes one adult and one infant manikin, one adult and one infant tracheostomy tube, an instruction guide, and a hard carrying case.



AB-1136

# ALS Trainers



**Adult ALS Trainer**

AA-4400

**AA-4400 Adult ALS Trainer Full Body with Simulator**

(see AA-550, pg. 12), \$3,020.00

**AA-4405 Adult ALS Trainer Full Body only, \$2,720.00**

**AA-4410 Adult ALS Trainer Torso with Simulator**

(see AA-550, pg. 12), \$2,615.00

The Adult ALS Trainer is a comprehensive, megacode manikin. It has an airway management trainer for both oral and nasal intubation. Practice airway management procedures for bag valve mask, Combitube®, and LMA. Laryngospasms can also be simulated.

All airway landmarks are realistically featured.

When the ALS Trainer is used with a simulator, standard ECG monitoring can be performed, as well as defibrillation, cardioversion, and pacing. The manikin also includes an I.V. injection arm and

several I.M. injection sites located on the deltoid and thighs.

See AA-550 on page 12 for Simulator features.


**AA-4415 Adult ALS Trainer Torso only, \$2,405.00**

**Z02 Zoll Defib Pad Adapter, \$150.00**

**Z10 Physio Quick Combo Defib Pad Adapter, \$150.00**

**Z13 Philips MRx Defib Pad Adapter, \$150.00**

**Use with the 12-Lead Arrhythmia Simulator, see pg. 12**



AA-3400

ECG Simulator not included

## Twelve-Lead ECG Placement Trainer

**AA-3400 Twelve-Lead ECG Placement Trainer, \$1,050.00**

A replica of a human thorax, this manikin has all appropriate landmarks necessary to train students in proper placement of ECG leads. Using reusable, color-coded magnets for leads, students learn correct placement as the appropriate LED lights up. A lighted signal also indicates when all leads are placed correctly.

The AA-3400 can be used with a rhythm simulator to enhance training by adding ECG readings that can be viewed on your monitor or printed on your twelve-lead ECG machine.

**Visit Our Website At**  
[www.armstrongmedical.com](http://www.armstrongmedical.com)



## Armstrong Intubation Manikin

AA-3000 Armstrong Adult Intubation Manikin with Case, \$995.00

The Adult Intubation Manikin's head extends for correct positioning, the lungs have a realistic compliance, and the skin is soft.

Practice bag valve mask ventilations, oral or nasal intubation, EOA, Combitube®, and other airway management techniques. Works well with curved or straight laryngoscope blades. Comes with realistic cricoid cartilage (for cricoid

pressure) and realistic teeth. Laryngospasms can be simulated.

Manikin features oral pharynx, tongue, epiglottis, vocal cords, stomach, trachea, right and left bronchus, esophagus with stomach, articulating jaw, upper and lower teeth, bilateral nasal passages, and soft face skin.

AA-3000L Lungs and Stomach, 3/pkg., \$50.00



AA-3000



AB-1875

## Larry Intubation Trainer

AB-1875 Larry Head with Mounting Board, \$893.00

Includes mounting board, case, and lubricant.

AB-1830 Larry Torso Manikin (not pictured), \$1,224.00

Disposable lower airway, case, and lubricant.

AB-1870 Larry Head only (not pictured), \$623.00

- Find more tracheal space for improved maneuvering
- Apply cricoid pressure, change the position of the trachea, and close the

esophagus in order to practice the Sellick maneuver

Attaches to either CPARLENE® or Resusci® Anne™ CPR manikins.


## Advanced Larry Intubation Trainer

AB-1650 Advanced Larry, \$1,275.00


AB-1655 Advanced Larry Head with Stand, \$1,122.00

These trainers offer tongue swelling and laryngospasms in addition to all the features of the standard models (see above). Advanced Larry simulates a non-anesthetized patient for practicing intubation, ventilation, suction, and CPR techniques. Practice oral, nasal, and digital intubation, as well as ET, EOA, PTL, LMA, EGTA, and Combitube®

insertion. Applying pressure to the cricoid cartilage changes the position of the trachea and closes the esophagus to allow realistic practice of the Sellick maneuver. With its slightly anterior position, swelling tongue, and vocal cords, Advanced Larry is perfect for any level of training.



AB-1650



AA-3100

## Laerdal™ Airway Management Trainer

AA-3100 Laerdal Airway Management Trainer, \$2,215.00

Includes mounting board, separate larynx model, spray lubricant, sanitation kit, and hard case. Simulates a non-anesthetized patient for teaching intubation, ventilation, and suction.

**Intubation:** Practice esophageal, bronchial, oro- and naso-tracheal, and oro- and naso-pharyngeal intubation. Features full-range head motion and a complete realistic airway anatomy. Tube

placement can be checked by visual inspection of lung expansion during ventilation, auscultation of breath sounds, or bronchoscopic evaluation of tip position.

**Ventilation:** Allows practice for correctly establishing an open airway. Manual or mouth-to-mask resuscitators can be used.

**Foreign Body Removal:** Practice oropharyngeal and tracheal suctioning.

# Intubation Trainers



AA-3050

## Infant Intubation

### AA-3050 Infant Airway Intubation Head with Carrying Bag, \$695.00

AA-3050 features both realistic size and anatomy for training and practicing airway management skills on an infant. Represents newborn - 1-month-old child.

Manikin head comes with lubricant and a carrying bag.

Quality and affordability make the AA-3050 an ideal training tool.

Equipment Sizes:

- Laryngoscope Blade: #0

- Endotracheal Tube: 3.0 mm
- Nasotracheal Tube: 3.0 mm

Features:

- Landmarks include uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, esophagus, inflatable lungs, and stomach
- Vocal cords, highlighted in white, provide easy laryngoscope viewing

## Child Intubation

### AA-3025 Three-Year-Old Airway Intubation Head with Carrying Bag, \$695.00

AA-3025 features both realistic size and anatomy for training and practicing airway management skills on a child. Quality and affordability make the AA-3025 indispensable for training.

Manikin head comes with lubricant and a carrying bag.

Equipment Sizes:

- Laryngoscope Blade: #2

- Endotracheal Tube: 3.5 mm
- Nasotracheal Tube: 2.5 mm

Features:

- Landmarks include uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, esophagus, inflatable lungs, and stomach
- Vocal cords, highlighted in white, provide easy laryngoscope viewing



AA-3025



AB-1900

## Child Airway Management Trainer

### AB-1900 Eight-Year-Old Child Airway Management Trainer with Hard Case, \$699.00

This anatomically correct, eight-year-old child manikin is perfect for practicing pediatric airway management skills, including oral, digital, and nasal intubation, and suction. Landmarks include nostrils, teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false and true cords, cricoid

ring, tracheal rings, trachea, esophagus, lungs, and stomach.

Equipment Sizes:

- Laryngoscope Blade: #2
- Endotracheal Tube: 4.5 mm
- Nasotracheal Tube: 3.0-3.5 mm

## Infant Airway Management Trainer

### AB-2000 Infant Airway Management Trainer, \$668.00

Mounted on a stand. Includes stand, hard case, spray lubricant, lungs, and stomach. Represents a 3-month-old child.

Equipment Sizes:

- Laryngoscope Blade: #1
- Endotracheal Tube: 4 mm

### AB-2010 Infant Airway Management, Head Only, \$469.00

Can be attached to Resusci® Baby Infant CPR Manikin.

Infant Airway Management Trainer features realistic size and anatomy for practicing airway management skills on an infant.



AB-2000



## Laerdal™ Infant Intubation Model

AA-3200 Laerdal Infant Intubation Model with Spray Lubricant and Protective Cover, \$745.00

Features realistic oropharynx, tongue, epiglottis, larynx, vocal cords, trachea, and a detachable tongue and airway set that allows for simple replacement of these parts. Correct tube placement can be visually checked via the transparent

trachea or the practical inflation test. Represents a three month old.

### Equipment Sizes:

Laryngoscope Blade: #1  
Endotracheal Tube: up to 4.0 mm



AA-3200



AA-3250

## Laerdal™ Neonatal Intubation Trainer

AA-3250 Neonatal Intubation Trainer, \$607.00

The Neonatal Intubation Trainer provides practice of intubation skills on a newborn baby. Robust and realistic, this model allows students to undertake training that is directly transferable to the clinical setting.

- Realistic anatomy
- Oral and nasal intubation
- Bag valve mask ventilation
- Correct tube placement can be checked by practical inflation tests


## Complete Auscultation Training System

AB-2710 Complete Auscultation Training System (Posterior Board not included), \$3,615.00

The Complete Auscultation System includes an Auscultation Trainer (see below), a Smartscope™, an Auscultation Board with Case, a remote control, and an amplifier/speaker system. Coupled with the Auscultation Trainer and Smartscope, the Auscultation Board guides students through the transition from anatomical drawings to a three-dimensional torso. Using the same technology as the Auscultation Trainer, the Auscultation Board will cause the Smartscope to emit various sounds dependent on the location of the stethoscope.



AB-2710



AB-2700

## Auscultation Trainer and Smartscope™

AB-2700 Auscultation Trainer and Smartscope, \$2,948.00

AB-2705 Additional Smartscope, \$1,010.00

Lung sounds include normal, wheeze, mono wheeze, fine crackle, coarse crackle, rhonchi crackle, stridor, cavernous, bronchial, pulmonary edema, friction rub, egophony, and pectoriloquy.

Heart sounds include normal, aortic regurgitation, mitral stenosis, pulmonary stenosis, holosystolic, midsystolic, S3 gallow, S4 gallow, systolic click, atrial septic defect, PDA, and VSD.

### Additional features:

- Invisible auscultation sites
- Palpable landmarks
- Completely wireless
- Remote controlled
- No internal speakers
- Sounds vary at different sites
- Color-coded heart and lung sound sites in manual
- Includes a hard storage case

# Airway Management Head and I.V. Models



AA-3310

## Adult Deluxe Airway Management Trainer

AA-3310 Adult Deluxe Airway Management Trainer, \$1,539.00

This adult airway trainer utilizes lifelike skin that stretches more than human skin and will not be damaged by aggressive intubation. Practice cricothyrotomy and other surgical techniques. The feel and anatomy emphasizing accurate visualization features makes this trainer the choice for advanced airway management with the convenience of board mounting. Features include laryngospasm and tongue edema. Both lungs and stomach

are visible for easy determination of successful intubation. All the anatomical landmarks are shown and the vocal cords are in a neutral position and highlighted for easy viewing. The upper teeth break out if inappropriate technique is applied during intubations. The cricoid area replacement skin allows for an economical way to teach and practice surgical techniques for cricothyrotomy.

## Intradermal Injection Simulator

AB-1010 Intradermal Injection Simulator, \$215.00

Practice intradermal injections safely and conveniently with the Intradermal Injection Simulator.

This unit consists of a molded arm section from the wrist to just below the elbow. Vinyl skin provides a realistic feel and appearance for a successful training experience.

Unit features:

- Eight sights for intracutaneous injections
- Skin welts which form during proper fluid injection
- Reusable practice sights

The Injection Simulator includes sealant, a syringe, a teaching guide, and a portable storage box.



AB-1010



REALISTIC  
DERMALIKE™  
SKIN

AA-9675

## Advanced 4-Vein Venipuncture Arm Pad

AA-9675 Advanced 4-Vein Venipuncture Arm Pad, \$341.00

Practice venipuncture techniques and "no stick" I.V. catheter placements.

The latex-free Arm Pad features realistic blue veins below a tissue-like pad. The Arm Pad's skin rolls when the vein is palpated. Veins are present in three different sizes and occur at various depths.

Increase difficulty of procedure and thus improve skills.

This lightweight, portable, economical, non-staining training model is an excellent alternative to larger, heavier, more expensive models.

## NEW Ultra Sound Phantom

AA-9825 Ultra Sound Phantom, \$468.00

AA-9826 Simulated Blood, 8 oz., \$29.00

The phantom provides unmatched training opportunities for ultrasound guided vascular accessing. Prefilled with stain resistant blood, this model also contains three different sized vessels for additional learning.

- Can be accessed numerous times and is self healing.
- Realistic that movement of soft tissue can be observed as needle tip enters the vessel.
- Flesh toned, curved contour helps transducer positioning





## Chester Chest™ *REALISTIC DERMALIKE™ SKIN*

AA-9800 Chester Chest, \$798.00

Life-size torso model for teaching and practicing the most common, long-term vascular access routes.

- Implanted ports
- External central catheters
- Peripheral central catheters

Chester Chest has three accessing inserts.

- Tilting/tipping ports
- Wandering/shifting ports
- Deeply placed ports

Practice fluid infusion, blood withdrawal, subclavian and external jugular line placement, and heparinization.

AA-9801 Case for Chester Chest, \$85.00

405 Replacement Chest Skin for Chester Chest, \$113.00

AA-9806 Replacement Arm with Single Lumen Port, \$273.00 (shown)

410 Optional Triple Lumen Port, \$68.00 (shown)

454 Optional Connector, \$25.00 (shown) \*Required if ordering 410.

417 Optional Implanted Peripheral Port, \$69.00 (shown)



AA-9800



AB-1145

## Central Lines Injection Trainer

AB-1145 Central Venous Cannulation Simulator, \$878.00

Designed for ACLS and ATLS testing stations. Landmarks include clavicle, sternal notch, and sternocleidomastoid. The trainer shows muscle and bone section, subclavian, internal jugular, external jugular, and carotid.

Includes replaceable skin, tubing sealant, one quart of imitation blood, needles, and a case.

- Subclavian vein approaches
- Internal jugular vein approaches
- External jugular vein approaches
- Swan-Ganz catheter insertion

AB-1146 Cannulation Skin and Vein Replacement Kit, \$195.00

AB-1147 Muscle and Bone Replacement, \$203.00

## Peter PICC™ *REALISTIC DERMALIKE™ SKIN*

AA-9815 Peter PICC™ Central Catheter Line Insertion Simulator with Soft Case, \$728.00

The only teaching model that allows the practice of placing peripherally inserted central catheter lines (PICC lines).


Model features anatomically correct cephalic, basilic, median basilic, jugular, subclavian, and superior vena cava. Also features a movable chin to simulate occlusion of the jugular to prevent the PICC line from taking this route.

Standard I.V. catheter placement is also possible with Peter PICC.

Measure proper catheter length using palpable ribs or measure catheter length outside the body. Confirm proper placement by viewing distal tip of the catheter in the viewable superior vena cava. Identify underlying veins below the replacement skin.



AA-9815



AB-1005


## Arterial Puncture Arm

AB-1005 Adult Arterial Puncture Arm, \$561.00

This arm is ideal for practicing to draw arterial blood samples at either radial or brachial sites. The instructor can control pressure to produce a lifelike backflow of blood in the syringe.

AB-1007 Arterial Arm Replacement Skin with One Set of Arteries, \$151.00

# I.V. Arm Trainers and Pneumothorax Simulator



AB-1000

## I.V. Arm

AB-1000 I.V. Arm, Caucasian, \$547.00

AB-1001 I.V. Arm, African-American, \$547.00

Realistic skin rolls as you palpate the vein and the characteristic "pop" can be felt as the needle penetrates the vein. Multiple access points allow practice for drawing blood, intramuscular injections, I.V. catheter insertions, and infusions. Comes with a hard carrying case.

AB-1002 Replacement Skin and Vein Kit, Caucasian, \$125.00

AB-1003 Replacement Skin and Vein Kit, African-American, \$125.00

AB-1004 Replacement Vein, \$32.00

## I.V. Arm and Hand

AA-4340 I.V. Training Arm and I.V. Training Hand, \$456.00

This is a fully functional I.V. training arm and hand. The student can palpate the arm and hand, which are made of a supple, resilient, lifelike skin. The vein can be located, the needle introduced, and then blood withdrawn or fluid injected. Simulated blood and a carrying tray are included.

AA-4346 I.V. Arm Only, \$265.00

AA-4347 I.V. Hand Only, \$265.00

141R Skin for Right Hand, \$48.00

142 Skin for Arm, \$48.00


143 5-ft. Vein for Hand and Arm, \$22.50

144 I.V. Reservoir, \$13.30

225 Blood Powder (for 1 gallon), \$17.40



AA-4340



AB-950

## Advanced Injection Training Arm

AB-950 Advanced Injection Arm, Caucasian, \$635.00

AB-951 Advanced Injection Arm, African-American, \$635.00

Includes simulated blood, two syringes, two fluid bags, a manual, and a case.

Complete venous access includes basilic, cephalic, accessory cephalic, median cubital, dorsal metacarpal, digital, and more. Blood flows to all injection points simultaneously from a single fluid bag. Soft, flexible fingers, wrist, and skin make for excellent realism.

LF-1124U Replacement Veins, \$114.00

LF-1122U Replacement Skin and Vein Kit, Caucasian, \$175.00

LF-1123U Replacement Skin and Vein Kit, African-American, \$175.00

## Pneumothorax Simulator

AA-9500 Pneumothorax Simulator, \$595.70

This simulated torso helps emergency personnel learn the life-saving techniques associated with pneumothorax injuries.

- Life size
- Features rib cage, liner, and vinyl skin
- Ideal for PHTLS and EMT-P classes
- Excellent review for ER, OR, or ICU

The pneumothorax simulator includes a pump assembly for lung inflation, a dual

lung, a rib cage, replaceable skin, and a carrying bag.

The self-sealing dual lungs may be punctured approximately 75 to 100 times each, allowing double usage from original and replacement lungs. Procedural needle is not included.

### Replacement Parts:

1805 Replacement Lung, \$39.80

1807 Outer Skin, \$116.30



AA-9500





## Intraosseous Leg

**AA-7310 Intraosseous Leg with Stand, \$274.40**

Intraosseous access to the proximal tibia is offered with the I/O Leg Adult STAT Simulator. Use the I/O Leg to enhance training and understanding of procedures associated with B.I.G.™ and EZ-IO™ devices. This I/O device is not fluid capable and no blood source is accessible through the injection site.

*Use in conjunction with STAT Man (see page 14).*



## Intraosseous Infusion Simulator

**AB-1160 Intraosseous Infusion Simulator, \$663.00**

The Intraosseous Infusion Simulator enables students to practice intraosseous infusion with incredible accuracy and realism.

Provides a very accurate and consistent simulation of procedure. Due to the innovative bone design, students can palpate landmarks and puncture each bone. To decrease replacement costs, bones have been designed so all four sides can be punctured multiple times.

Includes a case, ten bones, two pairs of replacement leg skins, and accessories.

**AB-1161 Replacement Bones, 10/pkg., \$133.00**

Each side of the four-sided bone can be punctured approximately 4-6 times.

**AB-1162 Replacement Skins, 4/pkg. (2 pair), \$80.00**

## Aspiration Needles

**AE-1715 JAMSHIDI® Disposable Illinois Sternal/Iliac Aspiration Needle, 15 Gauge, \$25.00**  
15 gauge adjustable length: 3/8" - 1 7/8"

**AE-1718 JAMSHIDI® Disposable Illinois Sternal/Iliac Aspiration Needle, 18 Gauge, \$25.00**  
18 gauge adjustable length: 9/16" - 1 1/2"



## NEW Humerus Intraosseous (IO) Trainer


**AA-3600 Humerus Intraosseous (IO) Trainer, \$362.10**

**AA-3603 Replacement Injection Pads, 3/pkg., \$83.70**

The Humerus IO Trainer is the perfect tool to train the access procedure for fluid management. Humerus access is now accepted as the best alternative when IV access cannot be obtained. This trainer presents three orientations that can be demonstrated to provide the correct simulation for application. The

realistic trainer includes the correct hallmarks, which are needed to provide accurate and safe access for fluid administration. The trainer has an easy-to-access replacement bone section supporting rapid-fire training, while presenting the realism needed to make training effective.

# Procedure Simulators and Casualty Simulation Kits



## Intramuscular Injection Aid

**AB-1025 Intramuscular Injection Simulator, \$618.00**

This is the perfect volunteer for training health care students in the proper procedure for administering intramuscular injections. One buttock is shown dissected to reveal the major muscles, nerves, and vessels. The other side offers all the bony landmarks used to identify the desired gluteal injection area, including the superior end of the femur, the posterior and anterior

superior iliac spines, and the sacrum. Lateral injections in the thigh can also be demonstrated and practiced.

Skin texture, muscle texture, bone shape, and position all realistically resemble a live patient. Perforation of the tissue with a needle replicates the feeling of administering an actual injection.

Includes syringes, a teaching guide, and a carrying case.

## Wound Care Simulator

**AA-9680 Seymour II Wound Care Simulator, 538.00**

**AA-9681 Optional Carrying Case, \$88.00**

All four stages of pressure ulcers are displayed, some with the additional complications of undermining, tunneling, sinus tract, slough, infection, and exposed bone. Also displayed are a suspected deep tissue injury, unstageable full eschar/slough wound, and a 5 1/2" dehisced wound. The Stage III and Stage IV are positioned so that a "bridging"

dressings for use with a vacuum assisted closure and negative pressure wound therapy devices can be demonstrated and practiced. You are also able to demonstrate and practice wound cleansing, classification, staging, and assessment, as well as the measuring of wound length, depth, undermining, and tunneling.



AA-9680



AA-4615

## Deluxe Casualty Simulation Kit

**AA-4615 Deluxe Casualty Simulation Kit, \$815.00**

The Deluxe Casualty Simulation Kit contains seven bleeding strap-on wounds and fifteen never before seen moulage stick on wounds. The full line of accessories allows you to create seventy individual wounds. Wounds include

lacerations, fractures, burns, projectile entry, and much more. You can then use the included make up to create realism in your evolutions. The kit comes packaged in a carry case and contains instructions.

## Casualty Simulation Kits

**AA-4600 EMT Casualty Simulation Kit, \$1,076.10**

This kit contains everything necessary to stage incredibly realistic casualty scenarios. This 95-piece kit has burns, lacerations, compound fractures, gunshot wounds, eviscerations, and amputations. It also includes several varieties of make-up, simulated blood, and adhesives, so any injury can be represented. Packed in a rugged carrying case.

**AA-4625 Multiple Casualty Simulation Kit, \$383.60**

This budget-conscious alternative is a toned down version of the EMT Casualty Simulation Kit. It is designed for intermediate and advanced first aid courses. It features compound fractures, an open amputation, a sucking chest wound, a gunshot wound, lacerations, simulated blood, adhesives, grease paint, reservoirs and pumps, sprayer, and an instruction guide.



AA-4600



## Prepare For The Worst

### Terry Trauma Manikin™

AA-7200 Terry Trauma Adult Trauma and CPR Manikin, \$1,495.00

AA-8300 Terry Trauma Carrying Case, \$155.00

Complete with burns, lacerations, and broken bones, Terry Trauma offers realistic experience in identification, assessment, treatment, and transport of disaster victims. His flexible joints make him perfect for any kind of rescue, extrication, immobilization, or restraint.

Terry comes with two heads, a rigid rescue head and a second head for sanitary CPR training. He also includes twelve moulage pumps, simulated blood, make-up blood, CPR faces, airways and lungs, a chest/skin overlay, and clothing. Weighs 35 lbs.

Use any of the following wounds for a custom disaster scenario:



- ✓ Sucking chest wound
- ✓ Forehead and forearm lacerations
- ✓ Broken clavicle with seat belt contusion
- ✓ 2nd and 3rd degree burns of the face, hand, forearm, chest, and back
- ✓ Protruding intestines
- ✓ Compound fractures
- ✓ Jaw wound
- ✓ Projectile entry and exit on arm
- ✓ Impalement of the lower leg
- ✓ Crushed foot

### Rescue Randy™

AA-7500 Rescue Randy, 5'5" Ht., 55 lbs., \$870.10

AA-7505 Combat Challenge Randy, 6' Ht., 55 lbs., \$970.10

AA-7550 Rescue Randy, 5'5" Ht., 105 lbs., \$970.10

AA-7555 Combat Challenge Randy, 6' Ht., 105 lbs., \$1,079.20

AA-7575 I.A.F.F. Rescue Randy with Additional Reinforcements, 6'1", 165 lbs. (Does not include pants), \$1,376.00

AA-8300 Soft Carrying Case for 5'5" Manikin, \$155.00

AA-7570 Soft Carrying Case for 6' Manikin, \$170.40

Rescue Randy is a rugged, life-size adult manikin designed to simulate victims for emergency training programs. Randy's extra-tough, steel-framed body and flexible major joints can be flexed to fit any rescue scenario.

Rescue Randy is filled with a weighted foam that emulates both the weight and resiliency of a human body.

With the help of the I.A.F.F., a Rescue



Randy Manikin was developed which includes additional reinforcement for increased durability.

Ideal for training exercises:

- Extrication of victims from vehicles
- Fire rescue from smoke-filled rooms
- Removal of victims from accidents
- Application of anti-shock trousers

- Immobilization/restraint techniques
- Suspended victim scenarios

✓ EXTRA WEIGHT can be added to Randy at \$6.50/lb.

1112 Weighted Vest, 20 lbs., \$234.60

1114 Weighted Vest, 30 lbs., \$247.90

1115 Weighted Vest, 40 lbs., \$260.10

# X-ray Positioning Manikin

## X-Ray Positioning Manikin

AA-7800 X-ray Positioning Manikin, \$13,495.00

AA-8300 Manikin Storage Soft Case, \$155.00 (see page 11)

The X-ray Positioning Manikin allows students of radiography to repeatedly practice any standard projection without discomfort or danger to patients or themselves.

The manikin lies naturally in the neutral, prone, or supine positions. With the use of positioning aids, it is possible to place the manikin in other postures for all standard projections.

This life-size, human figure contains a fully-flexible, articulated skeleton made of a specially formulated polyurethane for realistic x-rays. No metal parts have been used in the articulation of the skeleton. The skin surface and internal organs (heart, larynx, lungs, and kidneys) are constructed of radiolucent materials.

The manikin can be positioned and for a whole body CT investigation. Good

images in terms of body outline can be obtained. Certain areas of the skeleton structure e.g. thoracic vertebrae and the head of the femur gave similar values to the typical patient. **NOTE: As the manikin is very different from the human body in terms of attenuation, it cannot be used for training in or practice of exposure control.**

**Manikin weighs approx. 30 lbs.**

## Life-Size Human Figure





## Advanced Geri™ and Keri™

AB-2425 Advanced Geri Geriatric Manikin, \$2,831.00

AB-2525 Advanced Keri Age-Generic Manikin, \$2,831.00

All the features of the Complete Geri and Keri Manikins, listed below, plus:

- I.V. Training Arm (AB-1000M) with replaceable skin, palpable veins, realistic flashback, and intramuscular injection capability.
- Blood Pressure Training Arm (AB-1035M) reproduces Korotkoff sounds, variable systolic and diastolic levels, pulse rate, volume, and auscultatory gap.
- Keri: 5'; Geri: 4'10"; each manikin weighs 36 lbs.



AB-2425 Advanced Geri Manikin Shown

## Realistic, Full-Featured Geriatric and Nursing Manikins

*More Than Thirty-Five Procedures Can Be Practiced*

## Complete Geri™ and Keri™

AB-2400 Complete Geri Geriatric Manikin, \$1,734.00

AB-2500 Complete Keri Age-Generic Manikin, \$1,734.00

### • Procedures

Urinary catheterization, enema administration, douching, eye and ear irrigation, intramuscular injection, oral and nasal lavage, gavage, and suctioning.

### • Care

Ostomy care, tracheostomy care, hair

washing, patient transfer and positioning, denture removal, hearing aid removal and insertion, bed baths, and clothing changes.

### • Superior range of motion

Trunk, shoulders, elbows, wrists, fingers, neck, hip, knees, ankles, and toes all

bend and flex realistically. Each manikin weighs 27 lbs.

### • Exams

Prostate exam, pap smear, visual inspection of ulcers and moles, and dilated pupil comparison.

AB-2405 Complete Geri Manikin with **Heart and Lung Sounds**, \$4,233.00

AB-2505 Complete Keri Manikin with **Heart and Lung Sounds**, \$4,233.00

Same features as the Complete Geri and Keri, plus these great heart and lung sounds.

Heart sounds include normal, holosystolic, systolic click, aortic

regurgitation, mid systolic, atrial septic defect, mitral stenosis, S3 gallup, PDA, pulmonary stenosis, S4 gallup, and VSD.

Lung sounds include normal, rhonchi crackle, bronchial, egophony, wheeze,

stridor, pulmonary edema, mono wheeze, coarse crackle, friction rub, fine crackle, cavernous, and pectoriloquy.

## Basic Geri™ and Keri™


AB-2450 Basic Geri Geriatric Manikin, \$1,377.00

AB-2550 Basic Keri Age-Generic Manikin, \$1,377.00

These manikins include many of the features of the Complete Geri and Keri Manikins except internal fluid reservoirs, ostomy care, enema administration,

oral/nasal lavage, gavage, and suctioning. Please note: Urinary catheterization can not be performed on the Basic Geri™ Manikin. Each manikin weighs 27 lbs.

AB-2550 Basic Keri Manikin





Please note that the walker is not included and the manikin cannot stand without support from personnel.

Grandma Chase™ was first developed by Dr. Chase around the turn of the century when a group of nurses asked him if he would create a training aid that they could use to teach patient care.

Today Grandma Chase™ has evolved into one of the most comprehensive patient care manikins and helps meet the ever changing challenges in training nursing students.

*Practice stoma and ulcer care, hygiene, patient transfer, and other basic nursing skills on this versatile manikin.*

## Grandma Chase Patient Care Manikin

AA-8500 Grandma Chase, \$1,145.00

Includes sweat suit and interchangeable genitalia in a carrying case.

### Training Applications:

- **Patient transfer/handling**

Trains for changing occupied beds, changing patients' clothing, positioning in chairs and wheelchairs, transfer from bed to chair, assisting patients to bedside commode, prepping and positioning a patient for surgery, and assisted bedpan use.

- **Patient hygiene**

Trains for bed baths, flossing and brushing teeth, decubitus ulcer care, washing and combing hair, shaving, and fingernail care.

- **Other applications**

Trains for bandaging fingers and toes, administering liquid feedings, catheter insertion, and ostomy care.

### Manikin Features:

- **Head and hair**

The neck is flexible. The grey hair can be shampooed and combed.

- **Nose and mouth**

Either can be used for stomach tube insertion. Teeth can be flossed and brushed.

- **Eyes and ears**

One ear accepts a hearing aid. Ears can be cleaned with dry cloth or cotton swab. Both eyes can be irrigated with eye drops and are removable.

- **Hands and toes**

Individual fingers and toes make Grandma Chase perfect for bandaging practice.

- **Male and female genitalia**

Both are provided and held in place with Velcro®. Catheterization or a vaginal douche can be simulated without passing fluid. Perineal care can be practiced.

- **Joints**

Manikin bends at the elbow, knee, shoulder, and hip.

- **Decubitus ulcer**

A foam decubitus ulcer is included and can be placed anywhere on the manikin to simulate a bed sore.

- **Other features**

Grandma Chase also comes with a tracheostomy, a rectum that accepts rectal thermometers or water enemas, and an ostomy stoma. Manikin weighs 38 lbs.

***GREAT FOR LONG-TERM, GERIATRIC,  
OR CONVALESCENT CARE.***





## Advanced Training Manikin

AA-8501 Advanced Grandma Chase with Case, Gown, and Interchangeable Male/Female Genitalia, \$2,195.00

Advanced Grandma Chase is the most comprehensive patient care manikin available. Designed specifically for nurse training programs, this lifelike, durable manikin allows high-quality practice and instruction. Advanced Grandma Chase offers all the realism of the Basic Grandma Chase manikin except tracheostomy, plus the added features below.

### New Hands-on Training Features:

- Carotid pulse
- Full-functioning airway with tongue, epiglottis, esophagus, and trachea with vocal cords
- Right and left bronchus for assessing intubation skills
- Articulating jaw
- Chest compresses for CPR
- Intramuscular injection sites at the deltoids, right gluteal muscle, and at the thighs
- I.V. practice capability on right forearm and dorsal side of right hand
- Interchangeable male/female genitalia with full-functioning catheterization – voids fluids
- Enema capability
- Stomach reservoir
- Lightweight and portable with detachable arms and legs

AA-8505 Replacement I.V. Arm and Hand, \$251.00

4024 Replacement Skin for I.V. Arm, \$75.00

4026 Replacement Veins for I.V. Arm, \$33.00

**INCLUDES FREE  
MANIKIN SOFT  
STORAGE CASE!**




*Wheelchair and glasses are not included.*

## Ideal For Nurse Training Programs



*I.V. Practice*

*Airway Management*



*Multiple Injection Sites*

# Pediatric IV Arm and Obstetric Manikins



AB-1015

## Pediatric I.V. Arm

AB-1015 Pediatric I.V. Arm, \$408.00

Practice venipuncture and intramuscular injection techniques used with children. Soft foam is used to simulate the deltoid muscle. Simulated bone in the shoulder defines and limits the injection area. Includes a blood donor set, an infusion set, a 3cc syringe, 1 pint of

simulated blood, a teaching guide, and a carrying case.

AB-1016 Replacement Skin and Vein Set, \$125.00

## O.B. Manikins

AA-9505 Deluxe O.B. Manikin, Caucasian, \$660.00

AA-9505B Deluxe O.B. Manikin, African-American, \$660.00

AA-9560 Basic O.B. Manikin, Caucasian, \$625.30

AA-9560B Basic O.B. Manikin, African-American, \$625.30

Realistically practice emergency childbirth techniques with this lifelike O.B. manikin.

- Opaque and transparent overlays
- Vertex, normal, and breech delivery
- Umbilical cords to clamp and cut

Includes disposable umbilical cords, clamps, an extra vulva, a transparent overlay, a skin-tone overlay, simulated blood powder, and a carrying case.

The Deluxe Model also includes a pregnant belly overlay for palpating the fetus (Leopold Maneuver).



AA-9505



AA-9550

## Forceps OB Delivery Manikin

AA-9550 Forceps OB Delivery Manikin, \$812.00

Anatomically correct pelvic model with premature and full term fetuses for practicing forceps and vacuum assist delivery. Includes simulated blood powder, an extra vulva, and a soft carrying bag. Pregnant belly overlay and transparent overlay can be added to this manikin.

- Realistic practice of delivery techniques (forceps, vacuum assist)
- Lifelike pelvic cavity
- Removable abdominal overlay
- Soft skin vulva (anatomically correct)
- Two infant fetuses (premature and full term with fontanel and cranial sutures)

## NEW Newborn Stat Patient Simulator

AA-7375 Newborn Stat Patient Simulator, \$2,850.00

- Advanced Intubation Head
- Tongue Edema
- Laryngospasm
- Accepts LMA and ET, OPA, NPA, LT
- Pulse Locations (carotid, brachial, femoral, plus L radial)
- PICC Site — Right Arm
- IV sites (umbilical vein, right side antecubital, dorsal hand, left saphenous pedal)
- I/O Leg — Right Side
- ECG Monitoring — 4 Lead
- Defibrillation (20 joules to 360 joules capable)
- Chest Decompression (bilateral sites no manikin damage, pad replacement only)
- Chest Tube Placement (left midaxillary site- fluid with insertion)



- Urinary catheterization (fluid return with insertion)
- Rectal Medication



## NEW Neonate Simulator with ECG

AB-2100 Neonate Simulator with ECG, \$2,150.00

- Defibrillation shock can be delivered through manikin or simulator
- Connect defibrillator/external pacer to simulator using adapters
- Built-in circuitry allows you to defibrillate and pace directly into the ECG simulator and observe ECG rhythms through the PADS connector
- Convert feature
- Select another rhythm to run immediately after defib discharge
- Pacing can be done on any manufacturer's defibrillator
- Battery saver feature powers-off simulator automatically when not in use
- 6 waveforms
- 17 adult/pediatric rhythms



## NEW Micro-Preemie Simulator

AB-1975 Micro-Preemie Simulator, \$310.00

- Simulates a 25-week neonate - the smallest, most realistically proportioned pre-term infant available
- Entirely flexible and offers training opportunities for a number of commonly performed procedures
- Extraordinary realism is achieved with the preemie's pose and lack of muscle tone
- Perfect for NICU, emergency medical

response teams, flight medics, mock code scenarios, nursing schools, medical schools, hospitals, and family education

- Includes Micro-Preemie Simulator, diaper, hat (color may vary), umbilicus, omphalocele, neural tube defect, bilateral chest block, unilateral chest block, airway, 3 cc syringe, 25-gauge butterfly needle, lubricant, blood powder, and bulb and tube assembly
- 5-year warranty

## Nita Newborn™ Training Aid

AA-9850 Nita Newborn with Soft Carrying Case, \$653.00

4-lbs., 16" female manikin for practicing vascular access in infants and newborns. Permits accessing, securing, dressing, site care, and maintenance for venipuncture, central catheters, PICC lines, and umbilical catheters. Nasal and oral openings allow placement of nasal cannulas, NG tubes, and feeding tubes.

Skin rolls as you palpate the vein and backflow confirms needle location. Self-sealing veins permit withdrawal, infusion, and heparinization. Numerous veins are accessible.

Includes a fluid bag with attached tubing, simulated blood, a diaper, and a case.



# Procedure Simulators



## Eva™ Pelvic Exam Model

AA-9540 Eva™ Pelvic Exam Set with Case, \$724.20

Eva can be used for anatomical instruction, abdominal palpation, speculum instruction, or advanced tactile comparisons.

Covered in soft, easy-to-clean vinyl, Eva offers detailed and lifelike features for realistic simulation.

Eva includes nine interchangeable inserts:

- Normal cervix
- Normal parous cervix
- Cervix with endocervical polyp
- Cervix with ectropion characteristics
- Cervix with neoplasia (carcinoma)
- Normal uterus
- Pregnant uterus (ten weeks)
- Two adnexal masses

## Female Catheter Simulator



AB-1045 Female Catheter Simulator in Carrying Case, \$628.00

Simulates the torso of a middle aged woman properly positioned for catheterization. External genitalia and perineum are lifelike and the labia minora can be spread to reveal the clitoris, urethral opening, and vaginal introitus.

Inside there is a bladder reservoir, urethra, and a simulated urethral sphincter.

Normal resistance and pressure are experienced as the catheter passes through the urethra, past the sphincter, and into the bladder. Fluid then passes through the catheter as it does with a real patient.

This Female Catheterization Simulator comes with lubricant, catheters, a teaching guide, and a carrying case.



## Male Catheterization Simulator

AB-1050 Male Catheterization Simulator, \$628.00

A lubricated catheter can be inserted in the urethral orifice, passed through the urethra, and entered into the bladder. Fluid flows from the catheter when the bladder is successfully entered. The normal restrictions caused by the mucosal folds, the bulbous urethra, and the internal urethral sphincter are present.

The simulator teaches proper positioning and movement of the penis for catheterization with a minimum of patient discomfort.

The Male Catheterization Simulator comes with lubricant, six catheters, a teaching guide, and a carrying case.

## Cath Trainers

AA-4700 Male Catheterization Trainer, \$258.10

AA-4725 Female Catheterization Trainer, \$258.10

Trainers have realistic feel and flexibility. 16 FR catheter can be used to demonstrate balloon inflation procedures. A fluid reservoir allows for urine output upon successful completion of procedure. Instructors are able to teach sterile technique for catheter insertion/removal, as well as patient prep for the procedure.





# CPR Microshield, Microkey, and Micromask



AE-5210

**AE-5210 CPR Microshield™ in Sealed Pouch, 10/box, \$70.00**


**AE-5215 CPR Microshield Protection Pack with Gloves, 10/box, \$80.00**

Waterproof pouch with Microshield and light blue, latex-free Nitrile gloves.



AE-5200

**AE-5200 CPR Microshield™ in Vinyl Case, 10/box, \$60.00**



AE-5510

**AE-5510 Microshield-plus CPR Barrier Device, 10/pkg., \$72.50**

- CPR compatibility; does not interfere with CPR technique
- 36 square inches of face coverage
- 2 cm. tunnel with a one-way valve and a 3M Filtrete™ filter
- For use with adults over 40 Kg. and children 10 Kg. - 40 Kg.

## Your Key To Safe CPR



### CPR Microkey™

**AE-5225 CPR Microkey, Blue Case, \$8.95**

Carry CPR preparedness on your key chain with CPR Microkey. Each key ring comes with the popular Microshield (shown above). Compact and easy to use, the Microshield leaves CPR technique unchanged once the unit is in position. The mouthpiece/airway's one-way valve eliminates contact with the victim's secretions.

### CPR Micromask

**AE-5300 CPR Micromask, \$11.00**

**AE-5375 CPR Micromask with O<sub>2</sub> Inlet and Head Strap, \$12.50**

**AE-5305 Replacement Non-rebreathing Valve and Filter, \$4.30**

- One-way, non-rebreathing valve
- 3M Filtrete™ Two-Way Filter – helps to protect both patient and rescuer
- Large valve opening for maximum ventilation volume with minimal effort
- Latex-free, soft, blue nitrile gloves included
- Can be utilized with bag mask resuscitators
- Soft mask conforms to various facial features and provides an excellent seal
- Suitable for adults, children, or infants
- AE-5300 and AE-5375 both have adjustable cuff inflation for even better facial conformity



AE-5300



# Resuscitation Masks and Shields



## Laerdal™ Pocket Mask™

AE-5000 Pocket Mask with Filter and One-way Valve, \$15.00  
 AE-5075 Pocket Mask with Oxygen Inlet, Head Strap, Filter, and One-way Valve, \$21.00

- Includes blue, latex-free gloves and anti-microbial wipe
- Prevents mouth-to-mouth contact with the victim
- Transparent so rescuer can check victim's lip color and secretions
- Replaceable one-way valve and filter
- Oxygen inlet version facilitates delivery of supplemental oxygen

82 04 10 One-way Valve, \$7.65  
 82 06 10 Filter, \$2.50

## Laerdal™ Face Shield Key Ring

AE-5700 Face Shield Key Ring, \$5.20

A Resusci Patient Face Shield contained in an easy-to-carry key ring. The Face Shield can be used for mouth-to-mouth resuscitation on adults or children, or mouth-to-mouth/nose on infants.

- Helps overcome hesitation to start resuscitation
- Prevents direct contact with victim's mouth, nose, and face
- Plastic sheet and hydrophobic filter help prevent possible contact with victim's saliva or blood



AE-5700



46 00 00

## Laerdal™ Face Shield

46 00 00 Laerdal Face Shield, case is five boxes of 10 sealed wallet packs (50 total), \$102.00/case

Laerdal Face Shields are designed for protection in real-life rescue applications. They can be used for mouth-to-mouth resuscitation on adults or children, or mouth-to-mouth/nose on infants.

- Helps overcome hesitation to start resuscitation
- Prevents direct contact with victim's mouth, nose, and face
- Plastic sheet and hydrophobic filter help prevent possible contact with victim's saliva or blood
- Fits into wallets, pockets, handbags, and first aid kits
- Only 3<sup>3</sup>/<sub>4</sub>" x 2<sup>3</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>16</sub>"

## Laerdal™ Pediatric Pocket Mask

AE-5100 Pediatric Pocket Mask with One-way Valve and Blue/Yellow Soft Pack, \$15.00

82 04 10 One-way Valve, \$7.65

Utilizing the same low-resistance one-way valve as the traditional Laerdal Pocket Mask (see above), in combination with a transparent soft silicone mask, The Laerdal Pediatric Pocket Mask is an excellent companion

for the professional rescuer and first responder who have contact with infants and children. Includes blue, latex-free gloves and anti-microbial wipe.



AE-5100



## EXTRAORDINARY ERGONOMIC DESIGN

Meets The 10 Most Commonly Recommended Features For An Ergonomic Chair.

AC-979 Anesthesia Chair, Greystone, Non-conductive, Naugahyde, \$720.00

AC-979A Anesthesia Chair with Adjustable Arms, Greystone, Non-conductive, Naugahyde, \$810.00

- Wide, non-tilt, contoured seat (21" x 19" x 3")
- Wide seat back (18" x 19")
- Extra-wide wheel base (28" diameter)
- Adjustable wide foot ring (20" diameter)
- 22" to 32" adjustable seat height
- **New** dual-wheel safety casters prevent chair kick out
- Full lifetime warranty on pneumatic mechanism (13 years on balance of components)

Meets The 10 Most Commonly Recommended Features For An Ergonomic Chair:

1. Adjustable back rest for lumbar support.
2. Minimal pressure on the coccyx.
3. Ischial tuberosities support 75% of the body's weight.
4. Pneumatically adjustable height.
5. Five-point base.
6. Casters.
7. Feet flat on the floor or on a foot rest.
8. Seat surface slanted slightly backward. Upholstery slightly porous and of rough texture.
9. Seat front edge padded and rounded.
10. Adjustable arm rest. A clearance of at least 30 cm between the undersurface of the operating room table and the top of the chair seat.



Polished 28" Aluminum Base



**FIVE  
COLORS  
TO  
CHOOSE  
FROM**



Claret  
AC-976  
AC-976A



Black  
AC-975  
AC-975A



Royal  
AC-978  
AC-978A



Greystone  
AC-979  
AC-979A



Forest Green  
AC-974  
AC-974A

NOTE: Color options are only available on non-conductive chairs.

# Anesthesia Chairs



AC-955

## Economy Anesthesia Chair *with Task Control*

AC-955 Economy Anesthesia Chair, Non-conductive, Naugahyde, \$510.00  
 AC-955A Economy Anesthesia Chair with Arms, Non-conductive, Naugahyde, \$600.00

These affordable, non-conductive chairs are equipped with **task control**, which makes these chairs more comfortable and easier to adjust.

### Features:

- Comfortable 3" thick seat that is 18<sup>1</sup>/<sub>2</sub>" x 16<sup>1</sup>/<sub>2</sub>" and has a seat board bumperguard to protect the seat edges
- Five-leg, cast aluminum base with 25" leg spread
- Pneumatic mechanism allows the seat height to adjust from 19" to 26"
- Upholstered backrest has a lumbar support curve and a plastic protective outer panel which wraps around for protection to the edges

- **New** dual-wheel safety casters prevent chair kick out



**NEW**

## Armstrong Anesthesia Chair *with Task Control*

AC-995 Armstrong Anesthesia Chair, Non-Conductive, Naugahyde, \$535.00

Task control makes this chair more comfortable, easier to adjust, and offers:

- Fixed seat tilt at a 3.5° rearward angle
- Easier backrest height adjustment
- Back tilt range of 11° forward and 15° rearward
- Backrest is attached with a steel bar that resists bending and provides firm support
- Backrest that adjusts up and down within a 4" range
- Pneumatic (air ride) adjustment mechanism with convenient fingertip

control provides adjustable seat height range from 21" to 31"

- Lifetime guarantee on air cylinder
- Independently adjustable foot ring
- Five legs for mobility and stability
- 2<sup>1</sup>/<sub>2</sub>" foam cushion, 15" seat diameter

### AC-990 Armstrong Anesthesia Chair, Conductive, \$680.00

Same as the AC-995 except that it is conductive.



AC-995



AC-925

## Foot-activated Pneumatic Stool

AC-925 Foot-activated Pneumatic Stool, Naugahyde, \$510.00

### *Hands-free Height Adjustment*

AC-920 Hand-activated Pneumatic Stool, Naugahyde, \$370.00

- Seat height adjusts from 22" to 29"
- Seat is 16" diameter x 2<sup>1</sup>/<sub>2</sub>" thick
- Foot control mechanism, located under the chair's base, is easy to operate with the press of a heel; helps keep hands free for maintaining sterile conditions
- High-resilient foam seat padding with vinyl upholstery creates a

comfortable and durable stool for all healthcare providers

- Non-chip, chrome plating and bumper guard seat make for a long-lasting, quality investment

### AC-924 Foot-activated Pneumatic Stool, \$500.00

Same as AC-925 except seat adjusts from 18" to 23".



# When Comfort Is Not Optional

AC-985 Anesthesia Chair, Black,  
Non-conductive, Naugahyde, \$660.00  
AC-985A Anesthesia Chair with Arms,  
Black, Non-conductive, Naugahyde, \$750.00

## Body frame support:

- ✓ Ergonomically designed
- ✓ Balances work posture and weight support
- ✓ Reduces seat slippage
- ✓ Eliminates pressure on the coccyx
- ✓ Adjustable back/lumbar support  
(15<sup>1</sup>/<sub>2</sub>"W x 14<sup>1</sup>/<sub>2</sub>"H)

## Custom fit adjustment:

- ✓ Six-way seat control with high-density, foam seat  
(18<sup>1</sup>/<sub>2</sub>"W x 17"D x 3" Thickness)
- ✓ Free-floating, articulating seat control with more adjustment options
- ✓ Height adjustable seat (22" - 32") and foot ring (20" Diameter)
- ✓ Five-leg, cast aluminum base with 25" leg spread
- ✓ **New** dual-wheel safety casters prevent chair kick out
- ✓ Optional armrests with width and height adjustment

## Long-lasting value:

- ✓ Full lifetime warranty on pneumatic mechanism (13 year on balance of components)
- ✓ Bumper guards to protect against damage to the vinyl upholstery
- ✓ Polished, chrome plating

**FIVE  
COLORS  
TO  
CHOOSE  
FROM**

Claret AC-986 AC-986A	Black AC-985 AC-985A	Royal AC-988 AC-988A
Greystone AC-989 AC-989A	Forest Green AC-984 AC-984A	

**NEW**



# Glove Boxes and Clearview Security™ Drug Box

## Glove Dispensers



**AGB-1 Cart-mount Glove Box, \$55.00**

5.75"W x 10"H x 37/8"D.

**AE-5359 Glove Dispenser, 6/case, \$105.00**

Wall-mountable glove dispenser in a light beige. Adjustable shelf fits different sized boxes of gloves. 14 1/2" x 6" x 4 1/2".



**AGB-2P Plexiglas Glove Box, \$45.00**

10"W x 10 3/4"H x 3 7/8"D.  
Mounting hardware not included

**AGB-3P Plexiglas Glove Box, \$60.00 (not shown)**

Holds 3 boxes of gloves. 16"W x 10"H x 3 3/4"D  
Mounting hardware not included



**AGB-2 Steel, Baked-on Enamel Glove Box, \$55.00**

Holds 2 boxes of gloves. 10"W x 10 3/4"H x 3 7/8"D.  
Mounting hardware not included

**AGB-3 Steel, Baked-on Enamel Glove Box, \$85.00**

16"W x 10"H x 3 7/8"D.

**AGB-3V Steel, Baked-on Enamel Glove Box, \$75.00**

Holds 3 boxes of gloves vertically. 5 1/2"W x 30"H x 4"D.  
Mounting hardware not included



AGB-3V

*Glove Dispensers Can Be Mounted To A-SMART® Carts*

## ClearView Security™ Drug Box

**ACV-1PB Push-Button Locking ClearView Security Drug Box with Key Override, \$325.00**

**CVS-1 ClearView Shelf, \$80.00**

**APK-1 ClearView Privacy Cover Kit, \$30.00**

**ACV-1T ClearView Blue Tray, \$19.00**

The ClearView provides a great way to secure your medications, small equipment, patients' personal valuables, RSI equipment, and so much more. The ClearView also allows for a clear view of your inventory and also allows for anesthesia cases to be setup in advance while still keeping instruments/drugs secure. The ClearView can be mounted to a wall, in a patient's room, or even on an A-SMART® Cart.

### Features:

- Natural HDPE cabinet with clear PVC door
- Corrosion-resistant
- Provides a clear view of your inventory

**Dimensions:** 7 1/2"W x 14"D x 7 3/8"H

**Weight:** 5 lbs





## Powerful Yet Gentle Portable Aspirator

*Airflow >30 lpm At The TIP Of The Catheter*  
*Vacuum Adjustable <50mmHg - >525mmHg*



### AC Suction with Battery Back-up

AE-6975 SSCOR® DUET™ AC/Battery-Powered Suction Unit with Built-in Automatic Battery Care and Protection, \$995.00

AE-6976 SSCOR® DUET™ AC/Battery-Powered Suction Unit with Retention Bracket (mountable on carts or walls), \$1,135.00

The Duet is perfect for hospital crash carts, patient transport, surgi-center backup, long term care facilities and more.

- An AC suction unit with battery back-up, the Duet provides instant, powerful suction whether or not AC power is available.
- Built-in charger maintains battery and will not overcharge. Keep it plugged in at all times for maximum battery health.
- 3 year battery replacement policy for properly cared for devices!
- When fully charged, the sealed lead acid battery provides 45 minutes +/- 10% of suctioning power, giving you the suction you need for emergency care procedures.
- Battery condition indicator.
- Universal canister bracket system allows use of the disposable canister system currently stocked in the hospital.
- Pump Airflow: Approximately 50 lpm.
- Clinical Airflow: 30+ lpm (achieved at tip of patient tubing).
- Size: 17"L x 9"H x 6"W; weighs 10.5 lbs.
- Battery Type: Sealed lead acid.

Flat back unit sits conveniently on a cart or shelf. Bracket back unit will attach to the side of the cart or to a wall saving valuable shelf and cart top space.



**AHP-3 Standard Panel for Cart-mounted DUET**  
*(Retention Bracket not included),*  
**\$40.00**



**PHP-1S Premier Panel for Cart-mounted DUET**  
*(Retention Bracket not included),*  
**\$50.00**

ETL listed, conforms to UL STD 2601-1, 2nd Ed. Conforms to IEC 60601-1 and EN 60601-1-2, and Certified to CAN/CSA STD C22.2 NO. 601.1-M90

# SSCOR® Portable Suction

## NEW SSCOR® DCell Suction™

AE-6800 DCell Suction, \$1,250.00

The DCell Suction is the industry's first and only suction device powered by off-the-shelf D Cell alkaline batteries. Designed for the hospital and medical disaster management sectors, this unique feature allows the DCell Suction to provide high quality, continuous medical vacuum without the need to recharge. After depletion, batteries can be changed in minutes, making the DCell Suction's run time almost indefinite. The DCell Suction is ideal to have for disaster management caches and stockpiles. With alkaline D Cell batteries that carry a long shelf life, this unit can be stored for extended periods of time and deployed to provide hours of suctioning time when needed.

### Specifications:

12"H X 9" W X 11"D

Weight: 9.5 lbs. with batteries installed

Battery: 12 D Cell Alkaline (not included)

Vacuum Pump: 12 V DC exceeds 525 mmHg for the first few hours and exceeds 300 mmHg for the stated run time\*

Clinical Airflow: Greater than or equal to 30LPM (based on tests performed with fresh Duracell® alkaline batteries

## Emergency Preparedness Suction Unit



AE-6800

Batteries not included

### Features:

- 6-9 hour run time\*
- With universal canister holder, the DCell accommodates most standard hospital suction canisters
- Can also be powered indefinitely from external DC or AC power
- 3 year warranty

\*Run time is defined as the period of time that the device generates more than 300 mmHg negative pressure. Total run time typically varies between 6-9 hours. These run time specifications are based upon tests using new Duracell® CopperTop batteries. Run time will vary depending on battery manufacturer, state of the charge, ambient temperature, and age of the batteries.



ASU-1

### Suction Unit Holder

ASU-1 Suction Unit Holder, \$120.00

Dimensions: 10"H x 19"W x 7<sup>3</sup>/<sub>4</sub>"D.

### SSCOR® Transporter

AE-6977 SSCOR® Transporter, \$75.00

The SSCOR® Transporter allows the DUET to hang on the side of virtually all hospital beds, gurneys, and cardiac walkers.



AE-6977

### SSCOR® Accessories

AE-6996 UL 1310 Listed AC Fixed Voltage Charger, \$71.00

AE-6995 Suction Tip, 10/box, \$22.00

AE-6974 Disposable Canister with Bacterial Filter and Shut-off Valve, 1200cc, \$4.60 ea.; 48/case, \$196.80

43200 Patient Connecting Tubing, 6 ft., \$2.35 ea.



# Broselow® Pediatric Resuscitation System

## NEW Broselow® Pediatric Emergency Tape “2017 Edition” ml units are back!

AE-4800 Broselow® Pediatric Emergency Tape, 5/pkg., \$120.00

*When time and accuracy are crucial during a pediatric emergency situation, turn to the Broselow® System which is designed with children's care and procedural accuracy in mind. Each color-coded system supplies you with all the medical dosage information (except for infusion drugs) and size-specific equipment you need for your young patients' emergency resuscitation requirements.*

### QUICK AND SIMPLE

#### 1. Measure

Measure the child with the Broselow® Tape. This gives you an accurate weight estimate, equipment sizes, and drug dosages.

#### 2. Treat

Reach into the bag for exactly the right equipment for I.V., intubation, O<sub>2</sub> delivery, or intraosseous access. There's no guessing.


The 2017 version of the Broselow Reference Tape will be compliant with the 2015 PALS Standards and will include the following:

- Pre-calculated doses in milligrams, as well as all doses in pre-calculated ml, making the entire process error free by removing the need for any calculations.



- The new tape will accommodate for obesity by allowing for a habitus adjustment if deemed appropriate to the tape weight estimation.
- When used without the habitus adjustment new tape provides an accurate estimate of Ideal Body Weight, over actual body weight, which is preferred for many emergency medications.
- Updated with additional essential medications for both EMS and the ER.
- Contains a separate section, Quick Access Meds, for those medications needed immediately.

NEW



## NEW Broselow® Pediatric Resuscitation Medication/Infusion Guide 2017 Edition

AE-4875 Broselow® Pediatric Resuscitation Medication/Infusion Guide, \$210.00

This Resuscitation/Infusion Guide is designed to provide you a quick reference to common resuscitation I.V. push drugs and vasoactive infusion drugs. The guide is sectioned off in colors that correspond with the Broselow® Tape. The guide is also a solution for those who are trying to address The Joint Commission mandate to use standard concentration drugs for infusions and eliminate the use of “Rule of 6” methods of calculating infusions. The front section provides precalculated

milliliter amounts of the most common I.V. push drugs for resuscitation. The back section provides infusion rates based on standard concentrations of common vasoactive medications. This Guide acts to further enhance the Broselow Tape. The Broselow® Pediatric Emergency System helps to reduce cognitive overload sometimes associated with pediatric resuscitations, as well as reduces the potential for medication mistakes due to calculation errors.

# Broselow® Procedural Modules

## Broselow® Pediatric Resuscitation System — Procedural Modules

(MINIMUM ORDER \$100.00)

### Oxygen Delivery Module

STOCK NO.	DESCRIPTION	QTY.	PKG.
7700RAW	Oxygen Delivery Module (Pink/Red)	2	\$31.00
7700PAW	Oxygen Delivery Module (Purple)	2	\$26.00
7700YAW	Oxygen Delivery Module (Yellow)	2	\$26.00
7700WAW	Oxygen Delivery Module (White)	2	\$26.00
7700BAW	Oxygen Delivery Module (Blue)	2	\$26.00
7700OAW	Oxygen Delivery Module (Orange)	2	\$26.00
7700GAW	Oxygen Delivery Module (Green)	2	\$26.00

### Intubation Module

7700RIN	Intubation Module (Pink/Red)	2	\$75.00
7700PIN	Intubation Module (Purple)	2	\$75.00
7700YIN	Intubation Module (Yellow)	2	\$75.00
7700WIN	Intubation Module (White)	2	\$75.00
7700BIN	Intubation Module (Blue)	2	\$75.00
7700OIN	Intubation Module (Orange)	2	\$75.00
7700GIN	Intubation Module (Green)	2	\$75.00

### I.V. Delivery Module

7700RIV	I.V. Delivery Module (Pink/Red)	2	\$72.00
7700PIV	I.V. Delivery Module (Purple)	2	\$72.00
7700YIV	I.V. Delivery Module (Yellow)	2	\$72.00
7700WIV	I.V. Delivery Module (White)	2	\$72.00
7700BIV	I.V. Delivery Module (Blue)	2	\$72.00
7700OIV	I.V. Delivery Module (Orange)	2	\$72.00
7700GIV	I.V. Delivery Module (Green)	2	\$72.00

### Intraosseous Module

7700RIO	Intraosseous Module (Pink/Red)	2	\$90.00
7700PIO	Intraosseous Module (Purple)	2	\$90.00
7700YIO	Intraosseous Module (Yellow)	2	\$90.00
7700WIO	Intraosseous Module (White)	2	\$90.00
7700BIO	Intraosseous Module (Blue)	2	\$90.00
7700OIO	Intraosseous Module (Orange)	2	\$90.00
7700GIO	Intraosseous Module (Green)	2	\$90.00

### Broselow Nylon Color Packs

7730RED	Pink/Red Pack with Four Modules	1	\$168.00
7730PUR	Purple Pack with Four Modules	1	\$168.00
7730YEL	Yellow Pack with Four Modules	1	\$168.00
7730WHI	White Pack with Four Modules	1	\$168.00
7730BLU	Blue Pack with Four Modules	1	\$168.00
7730ORG	Orange Pack with Four Modules	1	\$168.00
7730GRN	Green Pack with Four Modules	1	\$168.00

### Broselow Flying Carpet

Each Flying Carpet contains the seven Color Packs pre-filled with twenty-eight appropriately sized Procedural Modules and one Broselow Measuring Tape. *Laryngoscope Handles sold separately. (Must be used with Vital View II™ handles)*

AE-4701	Flying Carpet with Broselow™ Pediatric Emergency Tape, Seven Pouches, & Twenty-eight Modules	1	\$1,240.00
---------	--	---	------------

### Broselow Outer Shell/Organizer and Contents

Contents for the Outer Shell include one BABY BLUE™ and one PEDI BLUE™ disposable manual resuscitator; infant, child, and small adult CUFF-ABLE® disposable blood pressure cuffs; a blood pressure manometer; and two Vital View II™ LED fiber optic handles; and one small and one large Magill forceps.

AE-4703	Contents for Outer Shell only (listed above) without Organizer	1	\$655.00
AE-4704	Outer Shell only	1	\$290.00


### Replacement Procedural Modules

AE-4712	Twenty-eight Procedural Modules without Nylon Color Packs. <i>Laryngoscope Handles sold separately.</i>	1	\$920.00
4556	Replacement Vital View II™ LED fiber optic laryngoscope handle	1	\$195.00



All equipment shown comes with complete system.

"Flying Carpet" Velcro® cloth for easy access to procedural modules.



Tape gives correct drug dosages and equipment sizes.

Size-specific modules for I.V., intubation, O<sub>2</sub> delivery, and intraosseous access, color-coded to match the child's size as measured by the Broselow® Pediatric Emergency Tape.

**Broselow® Pediatric Resuscitation System**

**AE-4700 Broselow® Pediatric Resuscitation System, \$1,950.00**

Includes twenty-eight size-specific procedural modules in seven pouches, flying carpet, Broselow® Pediatric Emergency Tape, and all other equipment shown above. Comes contained in a heavy-duty ALS equipment organizer.

The Broselow® Pediatric Resuscitation System is approved for the PALS education program. For children 3 to 36 kg.

**AE-4712 Procedural Modules, \$920.00**

A complete set of twenty-eight refill modules for the Broselow® System. *Laryngoscope handles sold separately.*

**AE-4703 Equipment Pack, \$655.00**

Two disposable manual resuscitators, four cuff masks, three disposable BP cuffs, one gauge, two Magill forceps, and two Vital View II™ LED fiber optic laryngoscope handles. *Call for details.*

**NEW AE-4800 Broselow® Pediatric Emergency Tape, 5/pkg., \$120.00**




Broselow® ColorCode™ Cart, see page 107. ClearView Security™ Drug Box, see page 112.

**Broselow Tape NEW for 2017**

# Broselow® System and Cart

**STEP 1** Measure the length of your pediatric patient with the Broselow® Pediatric Emergency Tape. Child's length (shown as purple at right) falls within a color bar and weight zone on the Broselow® Tape (see page 41).



A reference at each color bar on the tape informs you of equipment sizes to perform emergency resuscitation on the child. A reference at each weight/color zone on the tape shows pre-calculated medication dosages in milligrams. Designated resuscitation equipment is contained in corresponding, color-coded equipment pouches or drawers.

**STEP 2** Find specific equipment swiftly and accurately with the Broselow System as your guide. Color-coded pouches contain size-specific equipment needed for pediatric emergency resuscitation.

Use A-SMART® to create a complete pediatric resuscitation system. Premier™ Broselow® ColorCode Carts provide color-corresponding space and locking options to stock medications and additional equipment. Example on the right shows:

## PBL-PC-9A Premier™ Broselow® Complete ColorCode Cart System

Includes PBL-PC-9 (basic cart, see page 107), AE-4712 (procedural modules), AE-4703 (equipment pack), and a single Broselow™ Pediatric Emergency Tape.

Cart shown at right with *optional* ASU-1 Suction Unit Holder, AOB-1 Oxygen Tank Bracket, PIV-1 IV Pole, PBR-3 (8) 3" Breakaway Locking Bar, PBR-6 (1) 6" Breakaway Locking Bar, and AE-6975 S-SCORT® DUET Suction Unit (see page 39).



## Broselow® Flying Carpet



**AE-4701 Broselow® Flying Carpet, \$1,240.00**

Each Flying Carpet contains the seven

Color Packs prefilled with twenty-eight appropriately sized Procedural Modules (see page 66) and one Broselow

Measuring Tape. The seven Color Packs attach to the flying carpet with Velcro™. Laryngoscope handles are not included.



# Restock Your Broselow® Cart or System

## Broselow® Procedural Modules

AE-4712 Broselow Procedural Modules, \$920.00

A complete set of twenty-eight refill procedural modules for the Broselow® System. *Laryngoscope handles sold separately.*

Each color zone includes one of each module: Intubation, O<sub>2</sub> Delivery, IV Delivery, and Intraosseous.



## NEW Broselow® Pediatric Emergency Tape 2017 Edition

see page 41.



## Broselow® Equipment Pack

AE-4703 Broselow Equipment Pack, \$655.00

Includes two color-coded disposable manual resuscitators, four cuff masks, three color-coded disposable BP cuffs, one gauge, two Magill forceps, and two Vital View II™ LED fiber optic laryngoscope handles. *Call for details.*

AE-4709 Disposable Child Resuscitator, 10/case, \$300.00

AE-4705 Disposable Infant BP Cuff (Pink/Red/Purple), 5/case, \$75.00

AE-4706 Disposable Child BP Cuff (Yellow/White/Blue), 5/case, \$75.00

AE-4707 Disposable Sm. Adult BP Cuff (Orange/Green), 5/case, \$75.00

AE-4713 Manometer for Cuff, 5/case, \$167.00

8503 Magill Forceps, Child, \$15.00

8504 Magill Forceps, Adult, \$16.00

4556 Laryngoscope Handle, \$195.00

### Replacement Supplies:

AE-4708S Disposable Infant Resuscitator, \$27.50

AE-4708 Disposable Infant Resuscitator, 10/case, \$275.00

AE-4709S Disposable Child Resuscitator, \$30.00

Visit Our Website At  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

## Broselow® Emergency Equipment Organizers

AE-4825 Mini Broselow® Equipment Organizer, \$245.00

Seven empty color-coded bags permanently attached to "Flying Carpet" cloth.

Main bag: 15"L x 10"H x 6 1/2"D.

Pouches: 12"L x 5"H x 2"D

Flying Carpet: 30"L x 12"W

AE-4800 Broselow® Pediatric Emergency Tape, 5/pkg., \$120.00

AE-4805 Broselow® Pediatric Emergency Equipment Organizer, \$410.00

Seven empty color-coded bags attached to Velcro® "Flying Carpet" cloth. Includes airway management equipment bag.

Main bag: 18"L x 11 1/2"H x 12 1/2"D.

Airway bag: 14"L x 9"H x 4 1/2"D empty.

Pouches: 14"L x 8"W.



# HEINE Laryngoscopes Handles and Charger



## HEINE Standard & Rechargeable Handles

### Standard 2.5V Fiber Optic Handles

8720 Pediatric Handle, \$177.90

8721 Adult Handle, \$192.80

8722 Stubby Handle, \$227.50

### 3.5V Rechargeable Handle System

- Offers 40% brighter light than standard halogen bulb
- Saves on hidden costs of disposable batteries and bulbs

8720R Pediatric Handle with Rechargeable NiMH Battery and without Charger (requires the 8726S for charging), \$267.90

8721R Adult Handle with Rechargeable NiMH Battery and without Charger, \$341.80

8721RI Adult Handle with Rechargeable Lithium Ion Battery and without Charger, \$420.20

8722RB Stubby Handle with Rechargeable Lithium Ion Battery and without Charger, \$282.20

8726 NT 300 Charger, \$446.80

8726S Set of Small Adapter Sleeves, \$17.20

### Replacement Bulbs

8725 2.5V Xenon Halogen Bulb, \$32.40

8725R 3.5V Xenon Halogen Bulb for Rechargeable Handles, \$44.60

## HEINE LED Laryngoscope Handle & HEINE NT 300 CHARGER

8721LED Adult Rechargeable Handle with LED Light Source, Lithium Ion Battery and Status Indicator Cap (without charger), \$732.40

8721LI Adult Rechargeable Handle with LED Light Source and Lithium Ion Battery (without charger), \$581.10

8721LM Adult Rechargeable Handle with LED Light Source and NiMH Battery (without charger), \$499.10

8721L Adult Handle with LED Light Source, \$376.70  
Uses disposable "C" batteries.

8721I Adult LED Handle Insert only, \$252.60

8720LM Pediatric Rechargeable Handle with LED Light Source and NiMH Battery (without charger), \$391.90

8720L Pediatric Handle with LED Light Source, \$297.90  
Uses disposable "AA" batteries.

8722LM Stubby Rechargeable Handle with LED Light Source and NiMH Battery (without charger), \$436.80

8722L Stubby Handle with LED Light Source, \$383.40  
Uses disposable "AA" batteries.

Introducing the new Heine Laryngoscope Handles with LED Light Source. These LED handles offer illumination that is 3X brighter than a 3.5V Halogen Bulb. LED illumination


also allows for significantly longer usage time. Plus, the LED Light Source has a lifespan of 50,000 hours and comes with a 5 year warranty.

*Heine Handles with LED Light Source are not approved for Sterrad® Sterilizers at this time.*

### Additional Status Indicator Cap Option

The Status Indicator Cap allows the user to see a more specific charge status anytime. When you turn your handle on or remove it from the charger, the LED Indicator on the bottom cap lights up for 10 seconds displaying how much battery time remains. The Indicator also remains lit at any time a laryngoscope blade is attached to the handle.

*The Status Indicator Cap comes with the 8721LED Handle. One can also be added to any Heine Adult Lithium Ion LED or 3.5V Lithium Ion Halogen Handle (special Lithium Ion battery required). Call Customer Service for details.*



8726 NT 300 Charger, \$446.80

The NT 300 Charger offers an easy-to-see visual indicator on your handles charging status. The entire collar of the NT 300 is brightly illuminated. When you place your handle into the well to be charged, the light ring will pulse, indicating your handle is being charged. Once your handle is fully charged, the light rings stay brightly lit for 4 hours indicating the full charge. This unit prevents over charging the battery and helps to guarantee extended battery life. The NT 300 is the only charger that gives you this immediate visual clue as to where your handles stand in the charging process.

The NT 300 fully charges a handle with a dead battery in a maximum of 4 hours and is compatible with all Heine rechargeable laryngoscope handles.







## NEW HEINE XP Disposable Handle Sleeve System

F-00-22-926 Disposable Handle Sleeve System, 25/case, \$301.90

Call Customer Service For Details.

- Disposable handle outer shell for Joint Commission compliance
- Uses the same illumination inserts as the Heine reusable handles providing the benefits of disposable without the compromise of poor light quality
- Up to 4 times as bright as disposable handles when used with LED illumination
- Simple inventory management
- Extremely cost effective
- Recyclable design is more environmentally responsible than fully disposable handles
- Simple to operate
- Use with Heine Classic Fiber Optic Blades (*below*) for the best performing, cost effective, and compliant system on the market

*Inserts sold separately. Choose from LED or Halogen Illumination and Disposable Battery or Rechargeable Battery Systems.*

## HEINE Classic Fiber Optic Blades

- Large, bright, 4.3mm fiber optic bundle (Macintosh blades only)
- Choice between halogen or LED light sources
- Low maintenance; easy to clean
- Compatible with a variety of high-level disinfection and sterilizing methods - call for IFUs
- Five-year performance guarantee on fiber optic bundle; lifetime warranty against manufacturer's defect

### Macintosh Blades

- 8700 Neonatal (0), \$276.20
- 8701 Infant (1), \$276.20
- 8702 Child (2), \$276.20
- 8703 Adult Med. (3), \$276.20
- 8703M Adult Med. (3; strong curve), \$276.20
- 8704 Adult Lg. (4), \$276.20
- 8705 Adult XL (5), \$301.90

### Straight Blades (HEINE Paed)

- 8730 Neonatal (0), \$276.20
- 8731 Infant (1), \$276.20

### Wisconsin Blades

- 8711 Infant (1), \$301.90
- 8712 Child (2), \$301.90
- 8713 Adult (3), \$301.90

### Miller Blades

- 8714 Premature (00), \$276.20
- 8715 Neonatal (0), \$276.20
- 8716 Infant (1), \$276.20
- 8717 Child (2), \$276.20
- 8718 Adult (3), \$276.20
- 8719 Adult Lg. (4), \$301.90



### Macintosh Flexible Tip Blades

- Tip adjusts through 70° with lever
- Lifts epiglottis in smooth, controlled motion
- Adjustable blade tip facilitates difficult intubation
- Reduced risk of trauma to teeth and tissue

8703F Flexible Tip Mac 3 Blade with Case, \$976.30

8704F Flexible Tip Mac 4 Blade with Case, \$1,059.80

8734F Flexible Tips Mac 3 and Mac 4 Blades with Case, \$2,003.90

*\*Heine Laryngoscope handles, including all components (excluding LED Light Sources): outer sheath, inner sheath, bulb, alkaline battery/rechargeable battery, and all green ISO marking elements are approved for Sterrad® Sterilizers (NX, 50, 100S, 100NX, 200).*

*\*Heine Laryngoscope handles manufactured from May 2008 on are compatible with STERRAD NX and 100NX (Date code FM and on). Handles with no date code were manufactured before May 2008 can be upgraded for compatibility with STERRAD NX and 100X sterilizers.*

# Laryngoscopes

## Blades, Handles, and Lamps

All Economy Laryngoscope Handles are the American standard hook-on style. Blades are interchangeable with each other and are compatible with those of other manufacturers.

Sizes	Macintosh (Economy)	Miller (Economy)
	\$38.00	\$38.00
Premature (0)	8600X	8615X
Infant (1)	8601X	8616X
Child (2)	8602X	8617X
Adult Med (3)	8603X	8618X
Adult Lg (4)	8604X	8619X

### Economy Laryngoscope Handles

Chrome-plated brass handles with knurled finish and American standard hook-on receptacle housing. Batteries are not included. Handles are not pictured.

8620X Pediatric Handle, \$21.00

8621X Medium Handle, \$21.00

### Laryngoscope Lamps

Bright lamps with serrated metal cases and frosted lenses.

8622 Large Lamp, \$5.00

Fits 8603X, 8604X, 8617X, 8618X, and 8619X.

8623 Small Lamp, \$5.00

Fits 8600X, 8601X, 8602X, 8615X, and 8616X.



Miller

Macintosh

### Laryngoscope Sets

8630X Macintosh Set, \$229.00

Includes medium handle, pediatric handle, and four blade sizes: 1, 2, 3, and 4. Also includes AE-1245 Laryngoscope Soft Carrying Case (shown on page 68).

8635X Miller Set, \$267.00

Includes medium handle, pediatric handle, and five Miller blade sizes: 0, 1, 2, 3, and 4. Also includes AE-1245 Laryngoscope Soft Carrying Case (shown on page 68).

## HEINE Disposable Laryngoscope Blades

(for use with Green Systems)

### Disposable Macintosh Blades

8701D Infant (1), \$5.04

8702D Child (2), \$5.04

8703D Adult Medium (3), \$5.04

8704D Adult Large (4), \$5.04

### Disposable Miller Blades

8715D Premature (0), \$5.04

8716D Infant (1), \$5.04

8717D Child (2), \$5.04

8718D Adult Medium (3), \$5.04

8719D Adult Large (4), \$5.04

- Extreme torsional stiffness, practically unbreakable
- No danger of cross contamination
- Blade shape optimized to ease the passage of the tube



*More Heine Laryngoscopes,  
pages 46-47*



## Desk Charger With LED Ophthalmoscope, Otoscope & 2 Handles

AD-445 Desk Charger With LED Ophthalmoscope, Otoscope & 2 Handles, \$895.00

### 3.5V Otoscope:

- Fiber optic light transmission for an unobstructed view
- Optical glass reduces reflection and provides 3X magnification
- Bi-directional swivel lens allows instrumentation
- Insufflator fitting for pneumatic otoscopy

### 3.5V Ophthalmoscope:

- High-performance aspherical condenser lens reduces corneal reflection
- Synthetic rubber brow rest
- Ergonomic design and integral thumb rest improves stability and allows one finger control of all functions

### 3.5V Rechargeable Handle:

- 30 hours continuous use between charges
- LED lighting technology
- Electronic rheostatic light intensity control provides infinite continuous control and auto off feature after 2 minutes

### Desk Charger:

- Colored LEDs power displays charging state
- Storage compartment accommodates a small selection of specula or spare lamps
- Includes wall bracket for permanently mounting to wall or cabinet



Includes an anti-theft locking feature. 2 year warranty on instrument heads and handle, lifetime warranty on LED lamp.

## Portable LED Diagnostic Set

AD-440 Portable LED Diagnostic Set, \$595.00



### 3.5V Otoscope:

- Fiber optic light transmission for an unobstructed view
- Optical glass reduces reflection and provides 3X magnification
- Bi-directional swivel lens allows instrumentation
- Insufflator fitting for pneumatic otoscopy

### 3.5V Ophthalmoscope:

- High-performance aspherical condenser lens reduces corneal reflection
- Synthetic rubber brow rest

- Ergonomic design and integral thumb rest improves stability and allows one finger control of all functions

### 3.5V Rechargeable Handle:

- 30 hours continuous use between charges
- LED lighting technology
- Electronic rheostatic light intensity control provides infinite continuous control and auto off feature after 2 minutes

Includes an anti-theft locking feature. 2 year warranty on instrument heads and handle, lifetime warranty on LED lamp.

## Armstrong Ophthalmoscope & Otoscope Set

AD-425 Ophthalmoscope and Otoscope Set, \$260.00

AD-426 Otoscope Only, \$139.00


AD-427 Ophthalmoscope Only, \$196.00

Ophthalmoscope includes instrument head with halogen bulb, five apertures, and seventeen lenses. Otoscope is fiber optic with removable viewing lens, disposable specula (2.75 mm

and 4.25 mm), and insufflator fitting for pneumatic otoscopy. "AA" batteries included in each handle. 2.5V halogen bulb.



# Diagnostic Tools



## HEINE Wall-mounted Ophthalmoscope/Otoscope

AD-400 Heine Wall-mounted Ophthalmoscope/Otoscope, \$1,140.10

- Newly designed magnetic On/Off mechanism turns instruments on when removed from cradle and off when returned
- Connecting system is compatible with most existing instruments
- Exclusive "Soft-Start" transformer avoids current surge at switch-on; protects bulbs in the critical power-on phase to prolong bulb life
- Rheostats located directly on handles; brightness control at your fingertips
- Locking collar on handles prevents unauthorized removal of instrument heads
- High-memory spiral cords with 10-ft. range return to shape after each use; won't sag or droop
- Includes K 180 Ophthalmoscope and K 180 Otoscope

AD-406 Optional Wall Extension Unit with Handle (shown at left with AD-414 Slit Lamp Head; Slit Lamp not included), \$320.20


AD-402 4.0mm Specula, 1000/box, \$53.50

AD-403 2.5mm Specula, 1000/box, \$53.50

## NEW HEINE EN 200 Wall Transformer with Ophthalmoscope/Otoscope

AD-500 HEINE EN 200 Wall Transformer with Ophthalmoscope/Otoscope, \$1,334.00

- Lower power consumption saves costs and is environmentally-friendly
- The integrated galvanic separation protects the handles from unwanted electrical interaction
- Space-saving housing without any openings allows optimal disinfection
- Handle switches on and off automatically when taken off or returned to the unit
- Set your brightness level once and it stays there every time you use your instrument
- When not in use and when not returned to the unit, the illumination switches off after 15 minutes
- Connectivity ready with built-in USB port
- Integrated anti-theft system prevents unauthorized instrument head removal
- Lifetime warranty against manufacturing defects



## Wall and Mobile Aneroid Sphyg

AD-310 Wall Aneroid Sphyg, \$105.00

Large, luminescent, 6" dial with markings. Cuff stores conveniently behind the case. Tube stretches to 8 feet.


AD-315 Mobile Aneroid Sphyg, \$220.00

Large, luminescent, 6" dial with markings. Features five-leg caster stand with a range of 40"-54"H. Manometer angles back 15°.





*All Latex Free*



AD-200

## Armstrong Sphyg

AD-200 Adult, \$43.00

AD-200C Child, \$43.00

AD-200I Infant, \$43.00

AD-200-O Extra Large, \$55.00

Nylon blue cuff with printed index, range markings, and artery label. 300 mmHg, no-pin-stop manometer. Ten-year calibration warranty.

808 Blood Pressure Gauge, \$16.00



AD-215

## Premium Sphyg

AD-215 Premium Armstrong Sphyg, \$55.00

AD-215I Infant, \$55.00

AD-215C Child, \$55.00

AD-215-O Large, \$60.00

AD-215-B Bariatric (44.66cm), Burgundy, \$79.00

- Lifetime calibration warranty
- Heavily chrome plated
- Navy cuff includes gauge holder
- Luminescent dial
- Includes zippered case



AD-140

AD-220

## Quad Cuff

AD-140 Quad Cuff Multiple BP Kit, \$125.00

Includes child, small adult, adult, and large adult pre-gauged cuffs with artery label and inflation system. 300 mmHg, no-pin-stop manometer with bulb. Includes nylon case. Ten-year calibration warranty.

AD-130 Tri-Cuff Kit, \$110.00

Same features as the Quad Cuff, but without child size.

AD-220 Armstrong Economy Hand-held Adult Aneroid Sphyg, \$90.00



AD-135T

## Stethoscope and Aneroid Sphyg Set

AD-135 Stethoscope and Sphyg Combo Set, Black, \$49.00

AD-135T Teal, \$49.00

AD-135R Red, \$49.00

Color-coordinated combo sets include an aneroid sphyg with 300 mmHg no-pin-stop manometer, bladder, and bulb; and a Sprague-type stethoscope in a nylon carrying case. Twenty-year calibration warranty. Five-year warranty on Sprague.



AD-265

## E-Sphyg 2 Mobile Digital Aneroid Sphyg

AD-265 E-Sphyg 2 Mobile Digital Aneroid Sphyg, \$495.00

The E-Sphyg 2 replaces the vulnerable mercury column with a backlit digital simulation that is visually familiar to every practitioner. The device features dual mode operation, manual and automatic. In manual mode the practitioner measures using the traditional auscultatory method. In auto mode, the instrument will measure systolic and diastolic pressure and pulse automatically using child, adult, and large adult cuffs (included).

# Hand-Held Pulse Oximeters

## Hand-Held Pulse Oximeter System

AD-1500 SPECTRO<sub>2</sub><sup>™</sup> | 10 Hand-Held Pulse Oximeter System, \$415.00

- PI (Pulse amplitude Index) Indicator helps identify and maintain a good test site, and enables confidence in results
- One of the largest displays on the hand-held market
- On-board sensor storage cradle stores and protects sensor (reduces damage and minimizes replacement costs)
- Compatible with all BCI<sup>®</sup> sensors
- New BCI<sup>®</sup> 20" Spot-Check Sensor (keeps frustrating cable tangles to a minimum for quick, easy patient assessments)

### Specifications:

**SpO<sub>2</sub> (Functional)**  
Range: 0-99% (1% increments)

**Pulse Rate**  
Range: 30-254 bpm  
Averaging: 8 seconds

**Pulse Signal Strength**  
Range: 0-9 indicator, indicates logarithmic strength of patient pulse

**Display:** 9-segment red LED bargraph

**Power Requirements**  
**Battery:** 4 "AA" alkaline or the optional AD-1502, 2-cell lithium ion rechargeable battery pack (with optional AD-1503, docking station/charger)

**External (from AC Charger):** 9V, 3A output; input of 100-240 VAC, 50Hz, 60Hz

### Battery Life

**Alkaline Cell Life:** Approx. 32 hours of use

**Rechargeable Battery Life:** 54 hours of continuous use (new)

**Memory Storage:** data for up to 99 patients

**Dimensions:** 3.3"W x 6.1"H x 1.7"D  
**Weight:** .75 lbs. (with "AA" batteries)

**Standards:** EN60601-1, EN60601-1-2, ISO 9919

Comes with a 3 year warranty.

### Accessories for AD 1500

#### Disposable Probes for AD-1500

AD-923 Adult, 10/box, \$157.50  
AD-924 Pediatric, 10/box, \$157.50  
AD-925 Infant, 10/box, \$157.50

#### Reusable Probes and Supplies for AD-1500

AD-922 Infant Wrap Probe, \$157.50  
AD-1011 Ear Probe, \$168.00  
AD-1013S Finger Probe, \$168.00

#### Protective Boot and Case

AD-1501 Boot for AD-1500, \$33.00



## Digit<sup>™</sup> Finger Oximeter

AD-6000 Digit Finger Pulse Oximeter, \$185.00



The Digit Finger Oximeter is a compact, portable unit that combines both monitor and sensor into one unit for reliable SpO<sub>2</sub>, pulse rate, and pulse strength measurements. It is reliable for patients from pediatric to adult. The Digit Finger Oximeter is perfect for all environments where use is needed.

The Digit Finger Oximeter features a portable, lightweight design, a one button keypad, a large LED display, a pulse strength display bar that indicates relative perfusion at the measurement site, an auto power shutdown after 8 seconds when not in use, and a low battery indicator that flashes when 30 minutes of battery time remains.

Includes two "AAA" batteries. Comes with a two-year warranty.

### Specifications:

#### SpO<sub>2</sub>

Range: 0-99%  
Accuracy: ±2 at 70-99%  
Averaging: 8 beat avg.

#### Pulse Rate

Range: 30-254 bpm  
Accuracy: ±2% or 2 bpm  
Averaging: 8 second avg.

#### Dimensions

1.7"H x 2.25"W x 1.5"D

#### Weight

3 oz.



# Dopplers and Wright™ Peak Flow Meter

## Elite Doppler System

- AD-879 Elite 200 Doppler with 2 MHz. Probe, \$770.00
- AD-880 Elite 200 Doppler with 8 MHz. Probe, \$770.00
- AD-881 Elite 200 Doppler with 5 MHz. Probe, \$770.00
- AD-882 Elite 200 Doppler with 3 MHz. Probe, \$770.00
- AD-883 Elite 200 Doppler with 5 & 8 MHz. Probes, \$970.00
- AD-884 Elite 200 Doppler with 3 & 5 MHz. Probes, \$970.00
- AD-885 Elite 200 Doppler with 3 & 8 MHz. Probes, \$970.00

The Elite 200 Doppler lets you determine fetal viability as early as eight weeks into a pregnancy and assists in monitoring arterial and venous blood flow. This high-powered doppler is also able to locate both superficial and deep vessels faster and provide strong signals with a crystal-clear, static-free sound.

The standard heart rate display on the

Elite can provide for both OB and vascular applications.

### Optional Probes:

- AD-891 2 MHz. OB Probe, Waterproof, \$375.00 Later fetal heartbeat.
- AD-892 2 MHz. OB Probe, \$320.00 Later fetal heartbeat.
- AD-893 3 MHz. OB Probe, \$320.00 Early fetal heartbeat.



Uses a 9V battery.

- AD-895 5 MHz. Vascular Probe, \$320.00 Deeper veins and arteries.
- AD-898 8 MHz. Vascular Probe, \$320.00 Peripheral vessels.

## Mini-Wright™ Peak Flow Meters

PF-339X Mini-Wright Peak Flow Meter (60-880 lpm) with NAEP/ATS26 Scale, 10/pkg., \$285.00



These compact, portable Peak Flow Meters allow an asthmatic patient's progress

to be inexpensively monitored on a daily basis at home, work, or school.

- Help predict asthma episodes
- Allow therapeutic intervention
- Easily cleaned – minimal maintenance

The Mini-Wright measures the maximum flow rate in a forced exhalation – known as the peak expiratory flow (P.E.F.). After determining a patient's best P.E.F. with his/her physician, the patient self-tests daily as instructed by the physician and

**PF-440X AFS Low-range Mini-Wright (30-400 lpm) with Sterilizable Mouthpiece, 10/pkg., \$265.00**

The AFS Low-range Mini-Wright incorporates the Aerodynamic Flow System (AFS) to bring accuracy and consistency to patients with lower peak flow requirements.

records the highest reading on a daily trend chart.

The peak flow meter's scale also conforms to the 1991 recommended technical standards by the National Asthma Education Program (NAEP); National Heart, Lung, and Blood Institute (NHLBI); and is based on American Thoracic Society (ATS) Waveform 26.

These meters come with a PF-375 sterilizable mouthpiece with universal tapered design for both adult and pediatric use.

*NOTE: Children six years and younger, or those with excessive lung damage, should use the PF-440X low range meter.*

This new, smaller Mini-Wright is specifically designed for younger patients. The scale is clear and easy to read.

## Mouthpieces for Wright Peak Flow Meter

PF-310 Disposable Mouthpieces and PF-320 Disposable Mouthpieces for Pediatrics, (must be used with PF-380), 500/box: 1 box, \$60.00; 3 boxes, \$150.00; 6 boxes, \$270.00; 12 boxes, \$480.00

Mouthpieces are not made of cardboard; eliminates dust.

PF-315 Disposable Mouthpieces with One-way Valve, 200/box, \$105.00

PF-380 Mouthpiece Adapter for Pediatric Use (allows PF-320 to be used with PF-286), \$5.00

PF-375 Universal Tapered Mouthpiece, \$5.00

# Stethoscopes



AD-600

## Adscope® Cardiology Stethoscope

- AD-600 Adscope Cardiology Stethoscope, Black, \$115.00
- AD-602 Adscope Cardiology Stethoscope, Burgundy, \$115.00
- AD-604 Adscope Cardiology Stethoscope, Dark Green, \$115.00
- AD-606 Adscope Cardiology Stethoscope, Navy, \$115.00

- Convertible chestpiece provides adult diaphragm/bell or adult/pediatric diaphragm configurations
- Provides enhanced acoustic response in both adult and pediatric diaphragm
- Stainless steel chestpiece and aural tube construction
- Extra-deep bell for unsurpassed low frequency response
- Cardiology headset with large bore stainless steel binaurals



AD-103

## Neonatal Stethoscope


AD-103 Neonatal Stethoscope, \$25.00

- External binaural
- Dual-head brass chrome-plated chestpiece
- One-piece, PVC, 24", grey tubing

## Dual Head Stethoscope

- Dual Head Stethoscope
- AD-100 Black, \$18.00
- AD-100LB Light Blue, \$18.00
- AD-100GLD Gold, \$18.00
- AD-100GR Green, \$18.00
- AD-100R Red, \$18.00
- AD-100RB Royal Blue, \$18.00
- AD-100NY Neon Yellow, \$18.00

Precision-anodized, aluminum chest piece and adjustable, chrome-plated, brass binaurals. 24" tubing.



AD-100

Visit Our Website At [www.armstrongmedical.com](http://www.armstrongmedical.com)



AD-105

## Child Stethoscope

AD-105 Child Stethoscope, \$25.00

- Dual head with swivel
- Brass chrome-plated binaural
- One-piece, PVC, 21", black tubing



AD-120

## Dual Training Stethoscope

AD-120 Dual Training Stethoscope, \$25.00

Permits students to hear exactly what the instructor is hearing. A simple turn of the stem permits a quick change to a Ford-type bell.

## Sprague Rappaport

- AD-125 Black, \$19.00
- AD-125RB Royal Blue, \$19.00
- AD-125R Red, \$19.00
- AD-125NY Neon Yellow, \$19.00
- AD-125NOR Neon Orange, \$19.00
- AD-125NBL Neon Blue, \$19.00
- AD-125N Navy, \$19.00
- AD-125FG Frosted Ice, \$19.00

Each 5-in-1 stethoscope with 22" tubing comes packed in a box with two extra sets of eartips, two diaphragm sizes, and three precision bell sizes. Patented valve design and gasket system virtually eliminate sound leakage.



AD-125



## Privacy Screens

- AE-170 Three-Section Privacy Screen, Vinyl, White, \$190.00
  - AE-170-2 Fabric, Norway, \$255.00
  - AE-170-5 Fabric, Sage, \$275.00
  - AE-170-6 Fabric, Autumn, \$275.00
- For colors, see below.*

Fully expandable aluminum base allows screen to extend to full width. All panels are anti-bacterial and meet California flame resistant requirements. Each section measures 68"H x 19"W. Panel height is 52". Aluminum <sup>3</sup>/<sub>4</sub>" tubing with stainless steel hinges.

- AE-172 Four-Section Privacy Screen, Vinyl, White, \$210.00
  - AE-172-2 Fabric, Norway, \$295.00
  - AE-172-5 Fabric, Sage, \$325.00
  - AE-172-6 Fabric, Autumn, \$325.00
- For colors, see below.*

Features a straight aluminum base. All panels are anti-bacterial and meet California flame resistant requirements. Each section measures 68"H x 19"W. Panel height is 52".



## Privacy Screens

- AE-174 Three-Section Privacy Screen with Casters, Vinyl, White, \$227.00
  - AE-174-2 Fabric, Norway, \$295.00
  - AE-174-5 Fabric, Sage, \$310.00
  - AE-174-6 Fabric, Autumn, \$310.00
- For colors, see below.*

These Privacy Screens feature 2" diameter casters for transporting convenience. They are excellent for temporary privacy in patient/resident rooms when replacing soiled cubicle curtains. All panels are anti-bacterial and meet California flame resistant requirements. Each section measures 68"H x 19"W (includes casters).



## Privacy Screens

- AE-171 Three-Section Economy Privacy Screen, Vinyl, White, \$160.00
  - AE-171-2 Fabric, Norway, \$230.00
  - AE-171-5 Fabric, Sage, \$250.00
  - AE-171-6 Fabric, Autumn, \$250.00
- For colors, see above.*

These Privacy Screens are ideal for use in patients' rooms and physicians' offices. The accordion-hinged design allows panels to fold for convenient storage. The frame is anodized aluminum. All panels are anti-bacterial and meet California flame resistant requirements. Each section measures 68"H x 19"W.



**Scales**  
see page 56

# Scales




AC-611

## Physician Digital Scale

AC-611 Physician Digital Scale, \$486.00

- LB/KG Lock-Out ends prescription errors caused by LB/KG confusion during weighing
- Fully digital platform, from weighing sensors through readout, reduces errors caused by analog parts
- Easy-to-read 1" high LCD display
- Large, low-profile platform provides safer access and stability
- EMR connectivity via USB
- User friendly interface with full compliment of digital functions including BMI
- Capacity: 500 lbs.
- Full two-year warranty



AC-615

## Remote Display Scale

AC-615 Remote Display Scale, \$291.00

- Remote LCD display on 9' cord
- Weight reading in lbs. or kg.
- 2 oz./0.05 kg. graduation
- 9V battery and AC adapter included
- One-year limited warranty
- Capacity: 400 lbs.
- Full two-year warranty



AC-600

## Balance Beam Scale

AC-600 Physician Balance Beam Scale, \$275.00

Die cast, easy-reading pound beam, heavy duty  $\frac{3}{32}$ " steel base, and steel lever system comes with measuring rod (30" to 78") and Dura-Top™ non-peel ABS platform cover.

### Convenient Features

- Capacity: 390 lbs.
- Factory calibrated
- Simple field recalibration
- Graduation is engraved and paint filled; easy to read
- Full ten-year warranty

AC-602 With Optional Casters, \$344.00



AC-635

## Digital Pediatric Seat Scale

AC-635 Digital Pediatric Seat Scale, \$496.00

- 1" easy-to-read LCD display
- Auto-zero, auto-off, zero-out, LB/KG conversion, hold/release
- Capacity: 50 lbs.
- Seat size: 24 $\frac{1}{4}$ "W x 14 $\frac{1}{4}$ "D x 18 $\frac{1}{2}$ "H




AC-650

## Digital Wheelchair Ramp Scale

AC-650 Digital Wheelchair Ramp Scale, \$2,974.00

- EMR connectivity via USB
- Wheels and folding ramps for easy portability and storage
- Wide platform for wheelchair or bariatric use
- Calculates BMI
- Capacity: 1,000 lbs.
- LB/KG Conversion, LB/KG Lock Out, Body Mass Index, Zero Out/Tare, Hold/Release, Reweigh, Auto Zero



AC-640

## Digital Pediatric Scale

AC-640 Digital Pediatric Scale, \$406.00

High-performance digital pediatric scale with motion-sensing technology can capture the weights of even the most active babies.

- Easy-to-read LCD
- Large, stable weighing tray with raised edges for baby's comfort and safety
- Weighs in pounds and kilograms
- Waterproof buttons
- Features auto-off, auto-zero, and low-battery indicator
- Tare memory
- Full two-year warranty



# Slippery Slider and Hospital Equipment

## Armstrong Slippery Slider

AE-150 Patient Slippery Slider,  
22" x 72", \$139.95

AE-155 Patient Slippery Slider,  
18" x 72", \$139.95

AE-160 Patient Slippery Slider,  
18" x 48", \$119.95

Alleviate lifting strain, transfer large patients with minimum staff, and ensure patient comfort. Armstrong Medical offers three sizes of Slippery Sliders. Each has an ultra-smooth, polyethylene surface that makes patient transfer virtually free of friction.

- Reduces overexertion related back injuries
- Usable during x-rays; Slippery Slider is radiolucent
- Allows for easy grip with cut-out hand holes
- Strong, lightweight, and thin - all models only  $\frac{3}{16}$ " thick



AC-548

## Adjustable Laptop Computer Transport Stand

AC-548 Adjustable Laptop Computer Transport Stand, \$455.00

- Stand tray securely holds laptop in place
- Heavy-gauge coated steel construction
- Foot ring easily depresses to adjust stand height 35" to 45" for standing or seated use
- Easy rolling 3" casters provide a smooth, quiet ride
- Side access handles for "grab and go" use



AC-542

AC-540

## Mobile I.V. & Infusion Stands

AC-542 Armstrong **LIGHT DUTY** Custom I.V. Stand, \$148.00

Four-pronged, adjustable (47"-84") I.V. stand has a 20" diameter base and 2" casters. Weight capacity is 22 lbs.

AC-540 Armstrong **LIGHT DUTY** Economy Mobile I.V. Stand, \$75.00

Adjustable (47"-84") and easy to roll. Weight capacity is 22 lbs. 23" diameter base and 2" casters.



AC-549

## Security Laptop Stand

AC-549 Security Laptop Stand, 595.00

- Dual-locking work surface keeps laptop and data secure and provides a large work area
- Work area adjusts from 35"-45"
- Hand levers operate a gas spring lift to adjust height for sitting or standing use
- Heavy, five-leg aluminum base for safety and stability
- 3" rubber casters for quiet use


# Schrader & Disposables

## Schrader Connectors



Schrader Medical twist-release gas connectors are safety-keyed for air. This quick connector safely handle up to 200 psi. Design prevents interchangeability of major internal parts for different gases. They are UL certified and comply to NFPA 99 (2002 edition).

- Compressed Air Check Unit**  
 AS-4126 1/4" - 18 NPT Male, \$95.00
- Safety-keyed Adapters**  
 AS-4169 1/4" - 18 NPT. Male (A), \$65.00  
 AS-4147 1/4" Hose Barb (B), \$65.00



## Disposable Resuscitators

- AE-6500 Adult Bag Mask, 6/case, \$120.00  
 AE-6510 Child Bag Mask, 12/case, \$200.00  
 AE-6520 Infant Bag Mask, 12/case, \$200.00

These disposable self-inflating resuscitation devices feature a textured bag, double-swivel elbow, and a 40" tube-style oxygen reservoir. Bags come standard with an air-cushioned mask, an integrated pressure port for monitoring

inspiratory pressure, and an exhalation diverter to help safeguard the caregiver. Additionally, the child and infant bags include a pressure release valve as an extra layer of safety for patients with the most delicate lungs.

## Disposable Pressure Infusers

AE-1580 Infusable® Disposable Pressure Infusers, 500 ml, 12/case, \$208.00

AE-1585 Infusable® Disposable Pressure Infusers, 1000 ml, 12/case, \$228.00



## Sphygms

pages 50-51



Batteries not included

## NEW SSCOR® DCell Suction™

AE-6800 DCell Suction, \$1,250.00

### Specifications:

12"H X 9" W X 11"D

Weight: 9.5 lbs. with batteries installed

Battery: 12 D Cell Alkaline (not included)  
 Vacuum Pump: 12 V DC exceeds 525 mmHg for the first few hours and exceeds 300 mmHg for the stated run time\*

Clinical Airflow: Greater than or equal to 30LPM (based on tests performed with fresh Duracell® alkaline batteries)

For more information see page 40





## Speedsplint™

*Cost Effective For Hospitals and EMS. Disposable, No Cross-Contamination, Single-Patient Use.*

The splint can go with the patient to X-Ray or even surgery and be thrown out after. There is no need to worry about lost expensive splints.

### **AE-1605 Speedsplint™ Femoral Traction Splint, \$85.00**

- Includes splint, padding, harness, and groin pads
- Fits anyone up to 6'4" tall
- Has precise traction tension control
- Saves space

### **AE-1600 Speedsplint™ Multi-Purpose Splint, 10/pkg., \$80.00**

- Single splint fits bent or straight elbow, ankle, and wrist
- Unique fold-and-splint design reduces medic's pack weight, simplifies inventory requirements, and stores flat for easy transport
- Ideal for arm injuries, as well as foot and ankle injuries

### **AE-1604 Speedsplint™ Adjustable Leg Splint, \$65.00**

- Lightweight, temporary immobilization for all leg injuries
- Provides extra support above the knee for difficult-to-splint knee injuries
- Robust ankle support eliminates undesirable leg rotation
- Fits all leg lengths
- Three adjustable sizes: small, medium, and large

# Splints



## Sam® Splint

EM-500 Sam® Splint, \$15.00

Can splint an arm, leg, or hand. Only 4 1/4" x 36". X-ray invisible, waterproof, and reusable. *New flat design for better storage.*

EM-501 Sam Splint Jr., 4 1/4" x 18", \$9.50

EM-502 Sam Wrist Splint, 4 1/4" x 9", \$5.95

EM-505 Sam Finger Splint, 12/pkg., \$17.00

*Video upon request with purchase.*



AE-1120

## Sager Extreme® Bilateral Splint

AE-1120 Sager Extreme Bilateral Traction Splint, \$505.00

- Indicated for treatment in proximal third and mid-shaft femoral fractures
- Radiolucent to all areas of femoral fracture
- Compact, robust, and lightweight
- Fits adult or child
- Treats unilateral or bilateral fractures
- Unfolds and is ready to apply in under 15 seconds
- Rapid one-person application



## Prosplints®

AE-1800 Prosplints Complete Adult Splint Kit, \$250.00

Includes wrist/forearm splint, full arm splints (small and large), full leg splints (small and large), universal combo splint, universal rib splint, hard e-collars (small and medium), and carrying case.

- Quick application and removal
- X-ray translucent/MRI compatible
- Wipes clean for reuse
- Folds for compact storage



AE-1115

## Sager Form III Splints

AE-1110 Sager Form III Complete with Case, \$345.00

AE-1115 Sager Form III Bilateral Complete with Case, \$490.00

Like all Sager Splints, Form III Splints allow the continuous documentation of traction force applied – Quantifiable Dynamic Traction. The tough ischio-perineal cushion base, combined with the slender, anatomically-shaped cushion, mean that perineal examinations and procedures may be performed with the splint in place. Extra-wide leg cushion provides improved comfort and immobilization. For ages four and up.



AE-1875

## Vacuum Splints

AE-1875 Immobilize-Vac Deluxe Extremity Splint Set, \$455.00

Includes leg splint, arm splint, wrist/ankle splint, carrying case, pump, suction unit adapter, and extra valve.

AE-1880 Immobilize-Vac Pediatric Mattress Set, \$320.00

Mattress set includes mattress, case, pump, and adapter.

- Immobilizes injured areas without external pressure
- Eliminates the need for tape, boards, or gauze



AE-1105

## Infant Sager Splint


AE-1105 Infant Sager Traction Splint, \$410.00

The AE-1105 fits an inseam of 8" to 24". Single or bilateral application.

Quantifiable Dynamic Traction: the traction decreases as the spasm releases and prevents excessive traction. Dorsalis pedis and posterior tibial pulses are accessible. No heel stand and no risk of sciatic nerve compression. For ages six and under.



# Splints and Extrication Devices



## Hare® Traction Splint

AE-276 Adult Traction Splint, \$410.00

AE-281 Pediatric Traction Splint, \$410.00

Hare Traction Splint, the original traction splinting device, realigns a fractured femur to its proper position. Protect your patient from further hemorrhage and additional nerve, vascular, bone, and muscle damage.

AE-276 and AE-281 both include ratchet traction device and "S" ring hook, heel stand and release pull ring, collet locking device, four leg support straps, ischial strap and pad, adjustable ankle strap, carrying case, and instructions.



## Olympic™ Papoose Boards

AE-560 Papoose Board, 36" Regular (ages 2-5), \$405.00

AE-570 Papoose Board, 48" Large (ages 6-12), \$550.00

AE-580 Papoose Board, 60" X-Large (teen-adult), \$730.00

An effective way to immobilize a frightened, struggling child. Three sets of closing canvas flaps with Velcro™ closures hold the patient securely to the rigid backboard. Flaps can be individually released to expose the injured area.



AE-1050

## K.E.D. Extrication Device

AE-1050 K.E.D. Extrication Device with Case, \$165.00

K.E.D. has horizontal flexibility for ease of application in cramped areas, and vertical rigidity for spinal stability and safety. It is an invaluable product for any emergency team.

- Rapid application in tight spots
- Firm support for the victim's head, neck, and spine

Further assistance can be given to the victim while the K.E.D. remains in place:

- CPR can be easily performed
- MAST trousers can be applied without interference
- Traction splint can be applied
- Will not interfere with x-rays

## Pedi-Air-Align Complete System

AE-640 Pedi-Air-Align Complete System, \$305.00

The Pedi-Air-Align's patented dual-plane head drop system allows for a patient to be completely immobilized while their airway is managed. The board accommodates the larger occipital area of a child's head and creates neutral alignment of the airway which is critical to maintaining a patient's airway. Rotomolded from high-density polyethylene (HDPE), this one-piece, seamless product allows for complete decontamination.

- Patented dual-plane head drop system
- Lifetime guarantee
- Foam-filled, lightweight, ergonomic design
- Quick-strap, hook-and-loop strapping system

Board size: 48"L x 12"W x 1.75"H

Includes board, shoulder restraint, chest strap, four body straps, two wrist restraints, two head blocks, two head block straps, and a carrying case.



AE-640

# Cervical Immobilizers & Speedsplint™



AE-590

## SpeedBlocks™ Head Immobilizer

AE-590 SpeedBlocks™ Head Immobilizer, \$55.00

Includes one set of blocks, one universal base, one set of pads, and straps.

- Disposable head straps and block base pads allow for easy, low-cost replacement
- Four quick steps for application
- Blocks slide for multi-axis adjustment
- Blocks weigh only 5.2 oz. each
- Reliable handle-locking mechanism with quick-release locks ensures confident application
- Blocks provide contoured fit for patients two and up

- Universal base and blocks with replaceable foam pads for easy cleaning
- Base independently verified for greater than 99% decontamination efficacy
- Made of HDPE durable plastic and non-absorbing, padded blocks (radiolucent; MRI and CT scan compatible)

AE-591 Strap and Pad Replacement Set, \$21.00

Includes ten foam block pads, ten Velcro® head straps, and five foam base pads (five complete replacement sets).

## Head Immobilization System

AE-1031 Head Immobilization System, \$100.00

This lightweight, easy-to-use system can be used with most spine boards, extrication collars, or as a stand-alone immobilizer. Provides access to monitoring patient's ears. System offers a metal-free design, which allows for x-rays, CT scans, and MRIs. Made of

closed-cell foam, this immobilizer will not absorb blood or other bodily fluids. Also, the immobilizer is easily sanitized between uses.

Includes two closed-cell foam head pads, an attachment base, and two reusable head straps.



AE-1031

## Best Strap™ System

AE-543 Best Strap System, \$99.00

Introducing the Best Strap System, the total strap system for fast, firm, total body immobilization. Comes with five straps for the torso/limbs and two straps that go over the patient's shoulders. All straps are color-coded for easy application and less confusion. Straps close with Velcro®.

AE-543C Carrying Case, \$21.00



Board and immobilizer not included

AE-543



AE-1600

## Speedsplint™

AE-1600 Speedsplint™ Multi-Purpose Splint, 10/pkg., \$80.00


The splint can go with the patient to X-Ray or even surgery and be thrown out after. There is no need to worry about lost expensive splints.

- Single splint fits bent or straight elbow, ankle, and wrist
- Ideal for arm injuries, as well as foot and ankle injuries

- Unique fold-and-splint design reduces medic's pack weight, simplifies inventory requirements, and stores flat for easy transport



# X-ray Translucent Boards



AE-540

## BaXstrap™ Spineboard

AE-540 BaXstrap Spineboard, \$225.00


- One of the lightest boards on the market
- No seams; easy to clean and disinfect
- Contoured design; nests compactly, provides patient comfort and rescuer ease of use
- Extra-large, raised handholds
- Radiolucent; MRI and CT compatible
- Pins are molded into board for maximum strength
- Withstands greater than 2,500 lbs.
- Includes strap holes for pediatric patients
- Can be used for water rescue
- Limited lifetime warranty

## Ultra-Vue

AE-585P (yellow), AE-585-OP (orange), AE-585-BP (blue), Ultra-Vue Board with Pins, \$215.00

AE-585 (yellow), AE-585-O (orange), AE-585-B (blue), Ultra-Vue Board, \$200.00

- Weight rating: Approximately 1,100 lbs.
- 250 lbs. float capacity
- Size: 16"W x 72"L, 1<sup>3</sup>/<sub>4</sub>" thick; 15 lbs.
- X-ray translucent, MRI and CT scan compatible
- Flat patient surface won't compromise CPR efforts or aggravate cervical spine injuries



AE-585-O

## Backboard Straps

- AE-645 Two-piece Black 5' with Swivel Speed Clip, \$18.50
- AE-649 Two-piece Black 7' with Swivel Speed Clip, \$18.50
- AE-646 Two-piece Black 5' with Buckle, \$13.00
- AE-647 One-piece Black 7' with Buckle, \$13.00
- AE-648 One-piece Black 9' with Buckle, \$14.00

All clips and buckles are metal. Add "O" to the end of the product numbers above to designate orange.

## Ultra-Spac-Sav

AE-575 Ultra-Spac-Sav w/o Pins, \$415.00

AE-575P Ultra-Spac-Sav with Pins, \$430.00

This innovative plastic folding backboard features stainless steel mechanical pivot with internal aluminum runners. Stores in half the space of standard backboards, allowing it to fit in a backpack for search & rescue and fire apparatus.

- Weight rating: 500 lbs.
- Size: 72"L x 16"W, 1.75" thick
- Folded size: 39"L x 16"W, 4" thick




AE-575

## HeadBed™ II Cervical Immobilization Device

AE-3020 HeadBed™ II Cervical Immobilization Device, \$6.00

Disposable cervical immobilization device designed for effective motion restriction of infants, children, and adults.

- Outstanding patient access for effective assessment without compromising immobilization



AE-3020

# Extrication Collars



## Stifneck® Select™

AE-3015 Stifneck Select Adjustable Extrication Collar, \$14.00

AE-3016 Stifneck Pedi-Select™ Adjustable Extrication Collar, \$14.00

Four collars in one: The Stifneck Select adjusts to match the four adult Stifneck sizes (Regular, Short, Tall, and No-Neck™), so it's the only adult collar you need.

- Minimizes inventory and storage requirements

- Rear panel opens to allow C-spine palpation and visualization
- Large trachea hole for pulse check and airway procedures
- Radiolucent; MRI and CT compatible
- Adjustment tracks ensure selected size stays in place

## Stifneck Extrication Collars



AE-3000 Stifneck Collar, Regular, \$9.00

AE-3002 Stifneck Collar, Short, \$9.00

AE-3004 Stifneck Collar, Tall, \$9.00

AE-3008 Stifneck Collar, No-Neck™, \$9.00

AE-3006 Stifneck Collar, Pediatric, \$9.00

AE-3005 Stifneck Collar, Baby No-Neck™, \$9.00

One-piece, rigid extrication collar designed for optimal patient care.

- Low-angle chin ensures stable support without limiting airway
- Color-coded straps for easy and accurate sizing
- Fold-flat design conserves storage space
- Two pediatric sizes designed for children and infants

## Philadelphia® Tracheotomy Classic Collars

All Adult Philadelphia® Collars are \$18.00. Product numbers shown in chart below.

AE-806T Complete Starter Set, \$108.00

Includes AE-821T, AE-831T, AE-832T, AE-842T, AE-843T, and AE-853T.

AE-803T Basic Starter Set, \$54.00

Includes AE-831T, AE-842T, and AE-853T.

AE-810T Philadelphia Collar, Pediatric, \$18.00

AE-809T Philadelphia Collar, Infant, \$18.00 For ages 1-3.

- Total cervical arch support
- Two-piece design
- Suitable for extrication and long-term use
- Machine washable
- Radiolucent
- Large tracheal opening for pulse monitoring and tracheotomy

	SMALL	MEDIUM	LARGE
	(10"-13" neck)	(13"-16" neck)	(16"-19" neck)
2 1/4" HIGH	AE-821T	AE-822T	AE-823T
3 1/4" HIGH	AE-831T	AE-832T	AE-833T
4 1/4" HIGH	AE-841T	AE-842T	AE-843T
5 1/4" HIGH	AE-851T	AE-852T	AE-853T





## Philly® One-Piece™ Extrication Collars

- AE-708 "Low" Size (Purple), \$9.00
- AE-702 "Short" Size (Blue), \$9.00
- AE-700 "Regular" Size (Orange), \$9.00
- AE-704 "Tall" Size (Green), \$9.00
- AE-705 "Infant" Size (Pink), \$9.00
- AE-706 "Pediatric" Size (Aqua), \$9.00

- Compact, one-piece design – no assembly required
- Lightweight, injection-molded, high-density polyethylene construction
- Pre-molded level chin support
- Mandible support adapts to the wearer
- Color-coded straps and labels for fast, accurate size identification
- Enlarged trachea opening – permits easy access for pulse monitoring and tracheotomy procedures
- Rear ventilation slots
- Comfortable foam padding
- Essentially flat for compact storage in EMS vehicles and other confined spaces
- External sizing guide
- Radiolucent – allows x-rays to be taken with the collar in place

## *Speedblocks Head*

*Immobilization, see page 62*



## Philly® Patriot® Collars

- AE-720 Patriot Collar, Adult, \$9.50
- AE-721 Patriot Collar, Pediatric, \$9.50

With the Philly Patriot Collars, simply measure the patient, adjust and lock the selected height, and apply. Either in

supine or seated position, application is quick and easy.

- Restricts flexion, extension, and rotation to immobilize the cervical spine
- Multi-height adjustments (AE-720 has four heights; AE-721 has three heights)
- Injection-molded polyethylene plastic body for improved rigidity and durability
- Oversized trachea opening for quick access to carotid pulse monitoring and emergency tracheotomies
- Pre-molded level chin support maintains patient in neutral alignment; mandible support adapts to individual patient's jaw line
- Nasal cannula holders
- Open rear ventilation panel for fluid drainage, airflow, and cervical spine palpation
- Essentially flat for efficient storage
- Water resistant for use in adverse conditions
- X-ray, CT, and MRI lucent
- Clean with mild soap and water

# Morgan® Medi-Flow Lens

## AE-1695 Morgan Medi-Flow Lens, Sterile, 12/case, \$360.00

For immediate relief of pain and definitive treatment of:

- Ocular trauma and infection
- Foreign body sensation
- Operative and pre-operative procedures

The Morgan Lens makes difficult ocular irrigation and lavage safe, simple, and effective. Early insertion of the Morgan Lens can prevent permanent ocular damage and a months of unnecessary suffering.

The Morgan Lens eliminates the need for constant supervision during irrigation or lavage.

The Morgan Lens is currently used in over 4,000 emergency departments and is successful in treating hundreds of thousands of ocular injuries.

The lens, tubing, and adapter are sealed in individual peel-open packages with instructions for use printed on the front. Disposable units are sterile and ready for use.

The lens is made of soft, low-density polyethylene.

## AE-1697 Medi-Duct, 12/pkg., \$95.00

Optional Medi-Duct is an ocular fluid-management system designed for use with the Morgan Lens. When the lens is attached, solutions are carried away from the patient for easy collection.



*Instructional training video by Dr. Morgan is available from Armstrong on a loan basis. Please call Customer Service.*



For insertion, instill topical anesthetic. Attach syringe for lavage of I.V. set for continuous flow, using solution and rate of choice. Start the solution flowing.



Have patient look down. Insert edge of the lens under the upper lid. Have the patient look up and retract the lower lid.



Release lower lid over the lens and continue flow. Tape tube and adapter to patient's forehead to prevent accidental lens removal. Absorb outflow.



For removal, have patient look up. Retract lower lid behind interior border of the lens and hold position. Have patient look down, retract upper lid, and slide lens out.



## Morgan Lens® Delivery Set

### AE-1696 Morgan Lens Delivery Set, 6/box, \$45.00

The Delivery Set makes using the Morgan Lens even easier. It is no longer necessary to jury-rig separate I.V. set-ups to irrigate two eyes at the same time. Now, the lenses are simply connected to the Delivery Set and inserted using the

normal procedure. The Morgan Lens Delivery Set provides the speed and efficiency needed when a patient's vision is on the line. The Morgan Lens is sold separately.





## Universal Precaution Compliance Kit

AE-6398 Universal Precaution Compliance Kit, \$10.00

Designed in accordance with CDC and OSHA recommendations, this kit contains everything necessary to clean up blood, spills, and bodily fluids.

Universal Precaution Kit includes:

- Disposable protective gown
- Disposable vinyl gloves (one pair)
- Disposable mask with eye shield
- One packet of Red-Z absorbent
- Disposable scoop and scraper to pick up solidified spill

- Red biohazard disposal bag
- SaniZide® pre-moistened wipe for surfaces and equipment
- One P.A.W.S.® antimicrobial skin wipe for single use
- Detailed instructions

Glove Boxes,  
see page 38



AD-700 Fingertip Digital Pulse Oximeter, \$120.00

AD-701 Lanyard, \$2.00

AD-702 Bumper, Black, \$3.00

AD-704 Bumper, Orange, \$3.00

Reads SpO<sub>2</sub> and pulse rate. Has 6 display modes and visible and audible alarms.



AD-705 Admirals™ Pediatric Fingertip Digital Pulse Oximeter, \$90.00

- Reads SpO<sub>2</sub> and pulse rate
- Accuracy: SpO<sub>2</sub> +/-3%; HR +/-2 BPM
- Six changeable display modes
- Fits fingers with thickness of 0.2-0.63"



EM-820 Penlight, 6/pkg., \$8.95

EM-800 Ring Cutter, \$26.00

EM-810 Window Punch, \$9.95

EM-860 Rubber Tourniquet, \$2.95



AE-2000 Economy Emergency Scissors, Black, \$3.95

AE-2001 Emergency Scissors, \$7.95  
High-quality, precision, surgical steel construction. Can be autoclaved. 7<sup>1</sup>/<sub>4</sub>"L. Made in Germany.



AE-2002W White, \$5.95

AE-2002Y Yellow, \$5.95

AE-2002NP Neon Pink, \$5.95


AE-2002V Purple, \$5.95

AE-2002O Orange, \$5.95

AE-2002B Blue, \$5.95

AE-2002R Red, \$5.95

Hardened, surgical steel blades with one serrated edge. Can be autoclaved. 7<sup>1</sup>/<sub>4</sub>"L.



AE-1715 JAMSHIDI® Disposable Illinois Sternal/Iliac Aspiration Needle, 15 Gauge, \$25.00

15 gauge adjustable length: 3<sup>3</sup>/<sub>8</sub>" - 17<sup>7</sup>/<sub>8</sub>"

AE-1718 JAMSHIDI® Disposable Illinois Sternal/Iliac Aspiration Needle, 18 Gauge, \$25.00

18 gauge adjustable length: 9<sup>9</sup>/<sub>16</sub>" - 1<sup>1</sup>/<sub>2</sub>"

Visit Our Website At  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

# Armstrong Trauma Cases



AE-1230B

**AE-1230B Trauma Pack, Blue, \$215.00**

**AE-1230 Trauma Pack, Orange, \$215.00**

Cavernous, zippered main compartment and four full-size exterior pockets. Front pocket opens to reveal generous tacked elastic, full-length interior pocket, and three large vertical pockets.

Size: 19"L x 14"W x 12"H



AE-1246

**AE-1246 AED/Oxygen/Trauma Backpack Plus, \$495.00**

Heavy-duty construction and task-specific design ensures that this backpack delivers ease of transport and efficient deployment of AED, Oxygen, and Trauma supplies, all while keeping a slim profile. Features 1,000 denier nylon construction with all metal hardware, padded handles, and fluid-proof bottom.

Size: 21<sup>1</sup>/<sub>2</sub>"H x 13<sup>1</sup>/<sub>2</sub>"W x 11"D



AE-1289

**AE-1290 "D" Oxygen Bag, \$155.00**

**AE-1289 "E" Oxygen Bag, \$170.00**

Zippered bag provides immediate access to the tank, which is protected by durable padding. Contains mesh pocket to hold oxygen supplies.

"D" bag: 21"L x 8<sup>1</sup>/<sub>2</sub>"W x 8<sup>1</sup>/<sub>2</sub>"H

"E" bag: 32"L x 8<sup>1</sup>/<sub>2</sub>"W x 8<sup>1</sup>/<sub>2</sub>"H



AE-1240

**AE-1240 Breathsaver Oxygen Bag, Type "D" Cylinder, \$240.00**

**AE-1240LF Latex-Free Breathsaver Oxygen Bag, Red, Type "D" Cylinder, \$265.00**

One of the most popular airway management bag in the industry. Fold-back lid has gusseted pocket with clear window. Padded main compartment holds "D" or "Jumbo D" tank. Fluid-proof Bottom casing extends up the sides for added protection. Has two full-length side pockets and one end pocket.

Outside Dim: 26"L x 13"W x 12"H



AE-1245

**AE-1245 Laryngoscope Bag, Royal Blue, \$50.00**

Padded case spreads completely open for quick access to blades, handles, and tubes.

Size: 14"L x 1<sup>1</sup>/<sub>2</sub>"W x 10"H (closed)



AE-1242

**AE-1242 Ultra Softbox Plus Trauma Bag, Orange, \$355.00**

Features a large, semi-rigid main compartment with four removable dividers. A gusseted pocket with clear window in the lid holds loose equipment or the optional AE-1243 Ultra Meds Pack. Also includes four full-size, exterior pockets.

Size: 22"L x 14"W x 12"H

**AE-1243 Ultra Meds Pack, Royal Blue, \$65.00 (not shown)**

Compact, Cordura® nylon pack protects vials and ampules. Complete, fold-open design displays all areas at once for immediate access. Diced-foam blocks allow easy customization to suit your needs. Clear, vinyl pocket holds meds list. Fits in most soft trauma bags.

Size: 13"L x 3"W x 5<sup>1</sup>/<sub>2</sub>"H

**Please Note: Supplies shown are not included.**



# Armstrong Trauma Cases



AE-4704

**AE-4704 Empty Broselow® Bag (outer shell only), \$290.00**

Outer shell only ("Flying Carpet" not included).

18"L x 11<sup>1</sup>/<sub>2</sub>"H x 12<sup>1</sup>/<sub>2</sub>"D.



AE-4825

**AE-4825 Mini Broselow® Equipment Organizer, \$245.00**

Seven empty color-coded bags permanently attached to "Flying Carpet" cloth.

Main bag: 15"L x 10"H x 6<sup>1</sup>/<sub>2</sub>"D.  
Pouches: 12"L x 5"H x 2"D  
Flying Carpet: 30"L x 12"W



AE-4805

**Broselow®**

**AE-4805 Broselow® Pediatric Emergency Equipment Organizer, \$410.00**

Seven empty color-coded bags attached to Velcro® "Flying Carpet" cloth. Includes airway management equipment bag.

Main bag: 18"L x 11<sup>1</sup>/<sub>2</sub>"H x 12<sup>1</sup>/<sub>2</sub>"D.  
Airway bag: 14"L x 9"H x 4<sup>1</sup>/<sub>2</sub>"D empty.  
Pouches: 14"L x 8"W.

*Broselow® Tape and other supplies not included.*

**AE-4800 Broselow™ Pediatric Emergency Tape, 5/pkg., \$120.00**

*See pages 41-45 for a complete description of the innovative Broselow®/Hinkle Pediatric Resuscitation System.*



AE-1294

**AE-1294 Soft Tri-Pocket Case, Orange, \$135.00**

**AE-1294B Blue Case, \$135.00**  
The ideal disaster case. Adjustable compartments and removable interior provide rapid, organized access to equipment and supplies. Three pockets are perfect for adjuncts. Three-way carry with back pack, shoulder strap, or hand strap. 20"L x 15"W x 9"H.



AE-1291

**AE-1291 Standard Trauma Bag, Orange, \$60.00**

Cordura® bag has Snaploc™ buckles for easy access. Interior tool holders and separate compartments keep equipment in place. Two wing pouches and dual pockets hold accessories. Includes shoulder strap. 17"L x 9"W x 7"H.



AE-1282

**AE-1282 Large Plastic EMS Case, Orange and White, \$190.00**

Each of the three trays in this popular EMS case opens two ways, each swings back or pulls out as a drawer. Made of super-strength plastic, the case has an extra deep bottom.

20<sup>3</sup>/<sub>8</sub>"L x 11<sup>1</sup>/<sub>2</sub>"W x 12<sup>3</sup>/<sub>4</sub>"H.



**AE-1250V Vinyl Intubation Case, Orange, \$50.00**

*(not shown)*

Case rolls and fastens in a tight bundle. Holds twelve endo tubes, a full set of oral airways, and a laryngoscope handle with four blades; or equipment of your choice. Unrolled dim.: 25" x 14".

**AE-1250BV Vinyl Intubation Case, Blue, \$50.00**

# A-SMART® Carts Table of Contents

## Introduction

Standard Steel Cart Features. . . . .	72
Premier™ Aluminum Cart Features. . . . .	73
Standard Steel Cart Dimensions. . . . .	74
Premier™ Aluminum Cart Dimensions. . . . .	75
Create Your Own Cart. . . . .	76

## Anesthesia/Procedure

<b>Wireless Auto-Locking Carts</b>	
Aluminum Wireless Auto-Locking Carts. . . . .	77-79
<b>Auto-Locking Carts</b>	
Aluminum Auto-Locking Carts. . . . .	80-81
<b>Bedside Carts</b>	
Push-Button and Auto-Locking Bedside Carts. . . . .	82-83
<b>SUPREME Anesthesia Cart Packages</b>	
Standard Steel Anesthesia Cart Package. . . . .	84
Premier™ Aluminum Anesthesia Cart Package. . . . .	85
<b>Anesthesia Accessory Packages</b>	
Standard Anesthesia Accessory Packages. . . . .	86
Premier™ Anesthesia Accessory Packages. . . . .	87
<b>Solid Color Mobile Workstations</b>	
Standard Steel Solid Color Key Locking Carts. . . . .	88
Premier™ Aluminum Solid Color Key Locking Carts. . . . .	89
<b>Two-Tone Mobile Workstations</b>	
Standard Steel Two-Tone Key Locking Carts. . . . .	90
Premier™ Aluminum Two-Tone Key Locking Carts. . . . .	91

## Emergency

<b>PREMIUM Emergency Cart Packages</b>	
Standard Steel Emergency Cart Package. . . . .	92
Premier™ Aluminum Emergency Cart Package. . . . .	93
<b>Emergency Accessory Packages</b>	
Standard Emergency Accessory Packages. . . . .	94
Premier™ Emergency Accessory Packages. . . . .	95
<b>Solid Color Emergency Carts</b>	
Standard Steel Solid Color Emergency Carts. . . . .	96
Premier™ Aluminum Solid Color Emergency Carts. . . . .	97
<b>Two-Tone Emergency Carts</b>	
Standard Steel Two-Tone Emergency Carts. . . . .	98
Premier™ Aluminum Two-Tone Emergency Carts. . . . .	99

## Specialty

<b>Specialty Carts and Accessories</b>	
Push-Button Locking Carts. . . . .	100
Electronic Touchpad Carts. . . . .	101
Electronic Touchpad Wide Cart. . . . .	102
ScopeSafe Holders. . . . .	102
Mini and Narrow Carts. . . . .	103-105
Broselow® Premier™ Aluminum Pediatric ColorCode Cart. . . . .	106-107
Premier™ Locking Bars. . . . .	108
Isolation Carts. . . . .	109
Cart Extender. . . . .	110
Security Accessories. . . . .	111-112
Trellis Systems. . . . .	112

<b>Cart Accessories. . . . .</b>	<b>110-121</b>
<b>Stainless Steel, Wire, and Utility Carts &amp; Shelving. . . . .</b>	<b>122-124</b>



PBL-AR-30 with optional Accessories, see page 99

**NEW**



PPE-AA-30 Aluminum Electronic Touchpad Cart, see page 101



# A-SMART® CARTS + AREA SALES REP = SATISFIED CUSTOMER A-SMART ADVANTAGE

## CONSIDER THESE VALUE ADDED SERVICES

# Before You Make A Longterm Investment

*Ask you local Representative for a FREE Demo of our most popular carts*

### Armstrong Medical Offers:

- Nationwide coverage (over 40 Sales Representatives)
- Pre-sale support
- Delivery support
- Post-sale support
- **The Best Overall Value Around**

### FEATURES OF A-SMART® PREMIER™ ALUMINUM CART

- **ALUMINUM Carts are 45% lighter than Steel**
- **Full Extension Ball Bearing Drawer Slides**  
*Ensures that your drawers will open and close smoothly every time.*
- **Quality Aluminum Construction**  
*One-piece outer shell construction providing you with added strength and stability. Locking components consist of zinc plated and galvanized parts.*
- **Double Side-Wall Construction**  
*Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.*
- **Two Interchangeable Locking Systems**
- **Customize Your Locking Configuration**  
*Lock individual drawers, double lock drawers, use multiple types of locks on the same cart.*
- **Slide-Thru Side Shelf**
- **Contour-Fit Plastic Top**
- **Four Vertical Aluminum Corner Mounting Tracks**  
*For easy accessory attachment (eliminates the need for holes).*
- **Stabilizing Frame With Bumper**  
*Casters are moved out from underneath to the outside of the cart. This extra wide base provides your cart with added stability.*
- **18 Colors Available**  
*Allows you flexibility in your carts appearance. Choose colors to specify department or function, or even use multiple colors on the same cart to differentiate drawers.*
- **Highest Quality Swivel Casters (two locking, one tracking)**  
*Quiet, durable, German-made casters. Smoother, easier maneuvering.*

**A-SMART® Carts Are Made In The USA**

**Armstrong Medical** INDUSTRIES INC.

575 Knightsbridge Parkway • P.O. Box 700 Lincolnshire, IL 60069-0700

800/323-4220 • 847/913-0101 • FAX: 847/913-0138  
[www.armstrongmedical.com](http://www.armstrongmedical.com) [csr@armstrongmedical.com](mailto:csr@armstrongmedical.com)

800/323-4220 • FAX: 847/913-0138 71

# A-SMART® Standard Steel Cart Features

## A-SMART® Setting The Standard

**Breakaway Locking Cart**

Breakaway Security **ABL Series**

Soft-grip Handles

Ball Bearing Drawer Slides

Double Side-wall Construction

Drawer weight capacity (all drawer sizes): 50 lbs.

Stabilizing Frame with Bumper prevents injuries. Casters are moved from under the cart to outside its perimeter. This extra-wide base provides your cart with added stability.

(1) Swivel  
(2) Swivel w/Brake  
(1) Tracking Caster

### Key Locking Cart

Key-locking Security Series



Look for the Armstrong A-SMART® label on every cart. ***It stands for quality.*** Manufactured to ISO 9001:2008 certified standards.



APT-2 Optional Plastic Top, see page 113.

Key override comes standard on Push-Button Security **APB Series**



### A-SMART® Standard Features:


- Ball bearing drawer slides
- High-quality swivel casters (two locking, one tracking)
- **TENTE** Casters
- Double Side-wall construction (see right)
- Soft-grip handles
- Stabilizing frame with bumper
- 18 Colors
- Durable powder paint
- All carts are latex free

### Double Side-Wall Construction

Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.



## Take The Strain Off Your Back. Think *LIGHT*.



Customize your locking configuration; lock individual drawers, double-lock drawers, use multiple types of locks on the same cart.

### Key Locking Bar

Has all the same features as the breakaway locking cart. Security Locking Bar PKL Series



Look for the Armstrong A-SMART® label on every cart. **It stands for quality.** Manufactured to ISO 9001:2008 certified standards.

### See the Advantages of A-SMART® PREMIER™ Aluminum Carts

- Lightweight, aluminum construction
- Two interchangeable locking systems
- Four vertical Aluminum corner mounting tracks for easy accessory attachment
- Slide-thru side shelf
- Contour-fit plastic top
- High-quality swivel casters (two locking, one tracking)
- Stabilizing frame with bumper
- Double side-wall construction
- Soft-grip handles
- 18 colors
- Durable powder paint
- Ball bearing drawer slides
- All carts are latex free



### Double Side-Wall Construction

Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.


# A-SMART® Standard Steel Cart Features

## How to CUSTOMIZE An A-SMART® Standard Steel Cart:

### 1. CHOOSE FRAME SIZE & DRAWER CONFIGURATION

What kind of locking system do you need? How much drawer space do you want?

#### ABL Series: Breakaway Locks



**ABL-6**

6 Drws.  
= 30" drawer space

ABL-6 = 6 Drws. with 30" dr. space (shown above)  
ABL-4 = 4 Drws. with 24" dr. space  
ABL-5 = 5 Drws. with 24" dr. space


#### 6-Dr. ABL Cart Dim.:

Overall: 45.40"H x 25.40"D x 34.15"W (includes bumper/handles/casters).  
Drawer size: 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
Drawer configuration: 30" total drawer height.  
6-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (4), 6<sup>1</sup>/<sub>8</sub>"H (1), 12<sup>1</sup>/<sub>2</sub>"H (1).  
Weight: 165 lbs.

#### 4-Dr. and 5-Dr. ABL Cart Dim.:

Overall: 38.90"H x 25.40"D x 34.25"W (includes bumper/handles/casters).  
Drawer size: 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
Drawer configuration: 24" total drawer height.  
4-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (2), 6<sup>1</sup>/<sub>8</sub>"H (1), 12<sup>1</sup>/<sub>2</sub>"H (1).  
5-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (1), 9<sup>3</sup>/<sub>8</sub>"H (1).  
Weight: 145 lbs.

#### AKL Series: Key Locks



**AKL-6**

6 Drws.  
= 30" drawer space

AKL-6 = 6 Drws. with 30" dr. space (shown above)  
AKL-4 = 4 Drws. with 30" dr. space  
AKL-5 = 5 Drws. with 24" dr. space

#### 4-Dr. and 6-Dr. AKL Cart Dim.:

Overall: 45.40"H x 25.40"D x 34.15"W (includes bumper/handles/casters).  
Drawer size: 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
Drawer configuration: 30" total drawer height.  
4-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (1), 9<sup>3</sup>/<sub>8</sub>"H (3).  
6-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (2), 9<sup>3</sup>/<sub>8</sub>"H (1).  
Weight: 165 lbs.

#### 5-Dr. AKL Cart Dim.:

Overall: 38.90"H x 25.40"D x 34.25"W (includes bumper/handles/casters).  
Drawer size: 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
Drawer configuration: 24" total drawer height.  
5-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (1), 9<sup>3</sup>/<sub>8</sub>"H (1).  
Weight: 145 lbs.

#### APB Series: Push-Button Locks w/Key Override



**APB-6**

6 Drws.  
= 30" drawer space

APB-6 = 6 Drws. with 30" dr. space (shown above)  
APB-5 = 5 Drws. with 24" dr. space

**Standard Push-Button Locking Carts, page 100**

#### 6-Dr. APB Cart Dim.:

Overall: 45.40"H x 25.40"D x 34.15"W (includes bumper/handles/casters).  
Drawer size: 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
Drawer configuration: 30" total drawer height.  
6-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (2), 9<sup>3</sup>/<sub>8</sub>"H (1).  
Weight: 165 lbs.

#### 5-Dr. APB Cart Dim.:

Overall: 38.90"H x 25.40"D x 34.25"W (includes bumper/handles/casters).  
Drawer size: 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
Drawer configuration: 24" total drawer height.  
5-Drawer: 2<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (1), 9<sup>3</sup>/<sub>8</sub>"H (1).  
Weight: 145 lbs.

### 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain

### 2. PICK YOUR OPTIONS

#### Choose your drawer fronts...

Choose from any one of our eighteen designer colors on the right. If you want the ABL-4 as shown with Teal drawer fronts, your cart would be ABL-T-4. Customer Service will gladly assist you with questions.

#### Four drawer sizes available:

3"    6"    9"    12"


If standard drawer arrangements don't suit your needs, slides can be moved 3" up or down (additional charge per cart). Call Customer Service for details: 800/323-4220.

How do we ship our carts? See page 75



# A-SMART® Premier™ Aluminum Cart Features

## 18 Designer Colors



## How To CUSTOMIZE An A-SMART® Aluminum Premier™ Cart:

### 1. CHOOSE FRAME SIZE



#### 24" (5-Drw) CART DIMENSIONS:

**Overall:** 39.55"H x 25.02"D x 34.05"W (includes bumper, handles, and casters).  
**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
**Drw config.:** 24" total drawer height.  
**5-Drawer:** 27<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (1), 9<sup>3</sup>/<sub>8</sub>"H (1).  
**Weight:** 78 lbs.



#### 27" (6-Drw) CART DIMENSIONS:

**Overall:** 42.95"H x 25.02"D x 34.05"W (includes bumper, handles, and casters).  
**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
**Drw config.:** 27" total drawer height.  
**6-Drawer:** 27<sup>7</sup>/<sub>8</sub>"H (4), 6<sup>1</sup>/<sub>8</sub>"H (1), 9<sup>3</sup>/<sub>8</sub>"H (1).  
**Weight:** 85 lbs.



#### 30" (6-Drw) CART DIMENSIONS:

**Overall:** 46.05"H x 25.03"D x 34.05"W (includes bumper, handles, and casters).  
**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.  
**Drw config.:** 30" total drawer height.  
**6-Drawer:** 27<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (2), 9<sup>3</sup>/<sub>8</sub>"H (1).  
**Weight:** 92 lbs.

### 2. CHOOSE A LOCK OR HAVE BOTH


#### PKL Series: Key Locks



#### PBL Series: Breakaway Locks



#### Locks fit *either or both* corner glides



See page 108 for Premier™ Locking Bars

### 3. PICK YOUR OPTIONS

#### Choose your drawer fronts...

Choose from any one of the eighteen designer colors on the left. To get the PKL-24 as shown with **Slate Blue** drawer fronts, order the PKL-SB-24. To get the PBL-27 as shown with **Cranberry** drawer fronts, order the PBL-CR-27. Call Customer Service for assistance.

#### Four drawer sizes available:

3"    6"    9"    12"

If standard drawer arrangements don't suit your needs, slides can be moved 3" up or down (additional charge per cart). Call Customer Service for details: 800/323-4220.

#### HOW DO WE SHIP OUR CARTS?

All carts are shipped via truck ground freight. Special arrangements must be made for liftgate/inside delivery. Additional delivery time may occur on customized carts. Carts may not be returned for credit.

ACCESSORIES, Pages 110-121

Visit Our Website At  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

# Get Exactly What You Need

## CREATE YOUR OWN CART



With hundreds of accessories, eighteen available colors, flexible locking options (*Wireless Auto-Locking, Wireless Auto-Locking with optional Proximity Reader, Auto-Locking, Auto-Locking with optional Proximity Reader, Push-Button Locking, Key Locking, Breakaway Locking*), and multiple cart and drawer sizes, you can customize an A-SMART® Cart to fit any purpose or department. With so many options, you can build your own cart, so that you **Get Exactly What You Need. A-SMART.**

- |                         |                          |                       |                            |
|-------------------------|--------------------------|-----------------------|----------------------------|
| Ostomy Care Cart        | Trauma Cart              | Crash Cart            | ER Cart                    |
| Toxic Shock Cart        | Computer Cart            | Isolation Cart        | Anesthesia Cart            |
| Dressing Cart           | Delivery Cart            | Endoscopy Cart        | Phlebotomy Cart            |
| Conscious Sedation Cart | Latex-Free Cart          | Bio-Terror Cart       | PICC Line Cart             |
| Bronchoscopy Cart       | Ear, Nose & Throat Cart  | Difficult Airway Cart | Epidural Cart              |
| Eye Cart                | Sexual Assault Exam Cart | Broselow Cart         | Procedure Cart             |
| Block Cart              | 250-Yard Rule Cart       | Burn Cart             | Postpartum Hemorrhage Cart |
| Supply Cart             | ER Supply Cart           | Bedside Cart          | <b>And Many More</b>       |
| Respiratory Cart        |                          |                       |                            |


**Visit Our Website At [www.armstrongmedical.com](http://www.armstrongmedical.com)**

*Take The Virtual Tour Of Our A-SMART® Carts  
To See The Benefits Of Buying An Armstrong Cart.*



WIRELESS AUTO-LOCKING  
Not MRI Compatible

## 18 Designer Colors



# High Security Cart with Wireless Management

The Armstrong Wireless Auto-Locking Cart offers added security with several programmable auto-lock features. Additional security options include badge access, dual-locking controlled substance drawer with audible alarm, and software audit trail reports for accountability. The LCD display provides the cart status information and allows for on-screen programming. With the integrated 802.11g wireless capabilities, managing the carts and users (*up to 5,000*) can be administered remotely in only minutes.

## Wireless Auto-Locking Cart

### WEL-Color-30 6-Drawer 30" Wireless Auto-Locking Cart

#### Cart Features:

- Overall cart dimensions: 48.05"H x 25.20"D x 34.45"W
- Weight: 130 lbs.
- All the quality features of an Aluminum Premier Cart (see page 73)

#### Electronic Features include:

- LCD keypad display (shows date, time of day, current battery life, and programmable menu driven display system that prompts you to each additional command)
- Integrated 802.11g wireless capabilities
- 5,000 user codes possible\*
- Supervisor code for all programming
- User code for normal, daily use
- Manual or automatic locking
- Battery operated (6 "D" cell batteries)

#### Auto-lock features include:

- Programmable auto re-lock warning with the ability to extend the open time\* - Set a time interval to warn you that the cart is going to lock. When the warning happens, with a simple press of a button, you can choose to keep the cart open
- Time of day relock\* - Set the cart to lock at a specific time of day.
- User specific relock\* - The cart can be set to a default auto-lock time. This setting allows a specific user to set their own auto-lock time interval to differ from the default.
- Ability to program actual clock at keypad\*
- Ability to limit auto-lock programming to Supervisor Code only\* - Restrict users from manipulating the auto-lock settings.



WEL-CB-30 shown with optional Controlled Substance Drawer, WEL-PR Prox Reader, PTR-1 Trellis System, ATB-4X & ATB-5X Tilt Bins, and PSB-1 Sharps Container Brackets.

**\*Please note for the Armstrong PEL and WEL carts, we suggest using alcohol only for cleaning the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.**

**Not Sold Or Serviced Outside The Continental USA.  
MORE INFORMATION, SEE PAGES 78-79**

# Aluminum WIRELESS Auto-Locking Cart

WIRELESS AUTO-LOCKING  
Not MRI Compatible

## CONTROLLED SUBSTANCE DRAWER



Get added security for controlled substances with our CS drawer. User-friendly LED screen allows for easy navigation.



WEL-MG-30CS shown with Optional WEL-PR Proximity Reader, ACL-2 Long Catheter Holder, AWB-3 Sloping Wire Basket, PTR-1 Trellis System, and ATB-4X and ATB-5X Tilt Bins.

For optimal performance, be sure to check batteries every 6 months.

Not Sold Or Serviced Outside The Continental USA.

### Factory Installed Options

#### WEL-Color-30CS 6-Drawer 30" Wireless Auto-Locking Cart with Controlled Substance Drawer

The WEL-Color-30CS incorporates all of the WEL-Color-30 features in addition to a Controlled Substance Drawer

#### Controlled Substance Drawer:

- Top 3" drawer to be a dual credential Controlled Substance drawer.
- First credential can be either a proximity badge or a pin code, and the second credential is assigned as a pin code.
- Controlled Substance drawer can be locked manually or will auto-relock at the same time interval set for the main cart auto-lock
- Controlled Substance drawer alarm can be set to alarm (audio and visual) when the drawer is left open and/or alarm for an unauthorized entry
- Additional CS drawer key override

### Other Factory Installed Options

#### WEL-PR Proximity Reader

The WEL-PR Proximity Reader will read most low frequency HID® Prox Cards (125 kHz). *The Proximity Reader can not be retrofitted on existing carts.*

#### WEL-PRI iCLASS® Reader

Contactless smart card read-only reader. Reads identification only from iCLASS® cards. *The iCLASS® Reader can not be retrofitted on existing carts.*

#### PEL-PC HID® Prox Card

PEL-PS Proximity Card Sticker  
(makes non-compatible prox cards compatible)



Simply wave the card in front of the reader to open.

\*Please note for the Armstrong PEL and WEL carts, we suggest using alcohol only for cleaning the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.

### 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain



## Accountability & Wireless Management

### Network Software

#### WEL-SW Network Software

- Operates on 802.11g frequency (at designated check-in intervals)
- Comes with (1) client access license and (1) server license

\*Please call customer service for a complete list of requirements

#### WEL-SWC Additional Single Client Access License



#### Windows Based

The Network Software allows the administrator to communicate with the Wireless Carts through the facility's 802.11g network. As users and settings change, all of the changes can be made by the administrator from the comfort of their desk. As each Wireless Cart checks into the network at the designated check-in intervals, all information will be automatically updated at each and every cart.

Many configurable options including:

- Add, edit, view, and delete users
- Control of cart settings (ex: auto-lock defaults, alarm set-up, keypress volume level, etc.)
- View audit trails which include activity as to when the cart (and CS drawer when applicable) is open, locked, key override activity and alarms
- Email & SMS (Short Message Service) Alert allows for text message, fax, and voice notification system for alarms and low battery notification with provider
- Time interval for PC/WEL cart updating frequency
- eReport notification allows system to automatically send activity reports on a daily, weekly, or monthly basis

**\*Customer is responsible for installation.  
Before purchasing software, it is imperative to contact your IT department concerning installation.**



Lock Name	User Name	Type of Access	Status	Alarm Notes	Date of Entry	Time of Entry
Room 578	N/A	SETUPMODE	Key Override CS.		11/12/2010	8:41:10 AM
Room 578	N/A	SETUPMODE	Key Override main.		11/12/2010	8:41:05 AM
Room 578	N/A	N/A	Latch1 closed, Latch2 closed.		11/12/2010	8:40:59 AM
Room 578	N/A	N/A	Unauthorized Entry alarm cleared		11/12/2010	8:40:56 AM
Room 578	Larry	PUSHBUTTON	Access granted to CS.		11/12/2010	8:40:56 AM
Room 578	Larry	PUSHBUTTON	Access granted.		11/12/2010	8:40:52 AM
Room 578	N/A	N/A	Door closed.		11/12/2010	8:40:47 AM
Room 578	N/A	N/A	Door Alarm: Unauthorized Entry	See incident report #4576	11/12/2010	8:40:44 AM
Room 578	N/A	N/A	Door opened.		11/12/2010	8:40:43 AM
Room 578	User Not Found	PUSHBUTTON	Access Denied- No rights.		11/12/2010	8:40:40 AM
Room 578	asd	SETUPMODE	Menu access granted.		11/12/2010	8:40:21 AM

# Aluminum Auto-Locking Carts

**AUTO-LOCKING**  
Not MRI Compatible

## Added Security At An Affordable Price . . .

### Aluminum Auto-Locking Carts

PEL-Color-30 6-Drawer 30" Auto-Locking Cart

PEL-Color-24 5-Drawer 24" Auto-Locking Cart

#### Features:

- Keyless Entry
- 18 colors available
- Lightweight aluminum construction
- Master code for all programming
- Access code for normal, daily use
- Up to 250 user codes possible
- Manual or automatic locking
- Auto-Lock Timer can be set from 1 to 998 minutes, or it can be shut off so that the cart never auto-locks
- Auto-Lock Timer programming can be limited to master code only
- Battery operated; never needs to be plugged into an outlet
- "Low Battery" indicator light
- Manual key override for emergency access
- All the quality features of an Aluminum Premier™ Cart (except Slide-Thru Side Shelf)

For optional software, see page 81.

**Not Sold Or Serviced Outside The Continental USA.**

**For optimal performance, be sure to check batteries every 6 months.**

#### DIMENSIONS:

##### PEL-Color-30

**Overall:** 46.05"H x 25.20"D x 34.45"W (includes bumper, handles, and casters).

**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.

**Drawer:** 30" total drawer height.

**Drawer heights:** 3"H (3), 6"H (2), 9"H (1).

**Weight:** 102 lbs.

##### PEL-Color-24

**Overall:** 39.55"H x 25.20"D x 34.45"W (includes bumper, handles, and casters).

**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.

**Drawer config.:** 24" total drawer height.

**Drawer heights:** 3"H (3), 6"H (1), 9"H (1).

**Weight:** 88 lbs.

**\*Please note for the Armstrong PEL and WEL carts, we suggest using alcohol only for cleaning the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.**

A representative sample of this Auto-Locking Cart was tested to comply with FCC Part 15.

Manufactured to ISO 9001:2008 certified standards.



PEL-SG-30 shown above with optional PEL-PR Prox Reader, AWC-1 Waste Container, and AGB-1 Glove Box Holder.



PEL-B-24

#### 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain



## ... With Optional Proximity Reader.

**Think LIGHT.**

### Optional Factory Installed Proximity Reader

PEL-PR Proximity Reader

PEL-PC HID® Prox Card

PEL-PS Proximity Card Sticker  
(makes non-compatible prox cards compatible)

Simply wave the card in front of the reader to open. Lock the cart by pressing the lock button or set the Auto-Lock Timer for anywhere from 1 to 998 minutes. The Proximity Reader will read most low frequency HID® Prox Cards (125 kHz).

The Proximity Reader is a factory installed option on the Auto-Locking Carts only. This upgrade must be noted at the time of purchase. The Proximity Reader can not be retrofitted on existing carts. Electronic capabilities include all the features of the Auto-Locking Carts.


**NOTE: If you wish to use your current HID® cards, or purchase new cards or stickers, each card must be made available to be programmed individually into each cart.**

*For optimal performance, be sure to check batteries every 6 months.*



\*Please note for the Armstrong PEL and WEL carts, we suggest using alcohol only for the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.


### Auto-Locking Cart Software



Lock Audit Trail - Total Records Displayed: 16

Lock Name	User Name	Type of Access	Status	Date of Entry	Time of Entry
magombo	sample super	PUSHBUTTON	Access Granted	07/30/04	3:00:22 PM
magombo	marcus	PUSHBUTTON	Ok - PIN Required	07/30/04	3:00:20 PM
magombo	marcus	PUSHBUTTON	Denied-Restrictions	07/30/04	2:59:32 PM
magombo	marcus	PUSHBUTTON	Denied-Restrictions	07/30/04	2:54:04 PM
magombo	sample super	PUSHBUTTON	Access Granted	07/30/04	2:52:39 PM
magombo	marcus	PUSHBUTTON	Ok - PIN Required	07/30/04	2:52:36 PM
magombo	sample super	PUSHBUTTON	Access Granted	07/30/04	2:52:22 PM
magombo	marcus	PUSHBUTTON	Ok - PIN Required	07/30/04	2:52:18 PM
magombo	sample super	PUSHBUTTON	Access Granted	07/30/04	2:51:29 PM
magombo	marcus	PUSHBUTTON	Ok - PIN Required	07/30/04	2:51:20 PM

Access Log



User Editor

### PEL-SW Auto-Locking Cart Software\*

Software is Windows™ based and requires a USB connection. It allows you to add, edit, view, and delete users, supervisors, and specific carts (see screenshot examples to the left). The access log can display the time of entry by each assigned user. You may also create user groups, and restrict users by specified times of the day or days of the week.

**\*Customer is responsible for installation. Before purchasing software, it is imperative to contact your IT department concerning installation.**

## NURSE SERVER/BEDSIDE Carts



### Aluminum Push-Button Locking Narrow Cart

AMC-3Color 6-Drawer 30" Push-Button Locking Narrow Cart

Beige shell with eleven drawer colors available (see page 105).

#### Features:

- Aluminum construction; weighs only 85 lbs.
- Stabilizing frame w/ bumper
- Plastic top
- Four 5" swivel casters (two with brake and one tracking)
- Push-button locking with key override

Overall: 43.40"H x 25.40"D x 25.30"W (includes bumper and casters).

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (3), 6"H (2), 9"H (1).

### Steel Five-Drawer 24" Push-Button Locking Mini Cart

AMC-1Color-PB 5-Drawer 24" Push-Button Locking Mini Cart with Key Override

Beige shell with eleven drawer colors available (see page 104).

Overall: 34.50"H x 20.40"D x 20.90"W.

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (3), 6"H, 9"H.

Weight: 103 lbs.

Four 3" swivel casters (two with brake and one tracking)

### Steel Four-Drawer 21" Push-Button Locking Mini Cart

AMC-21Color-PB 4-Drawer 21" Push-Button Locking Mini Cart with Key Override

Beige shell with eleven drawer colors available (see page 104).

Overall: 31.30"H x 20.40"D x 20.90"W.

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (2), 6"H, 9"H.

Four 3" swivel casters (two with brake and one tracking)

Organize Your Carts With  
Mini/Narrow Cart Trays

Bottom of page 118



## Less Running, More Patient Time Nurse Server/Bedside Carts



*Stores Under Counters*



AMC-4CB



AMC-6AA



PEL-B-24

### Aluminum Five-Drawer 24" Auto-Locking Mini Cart

AMC-4Color 5-Drawer 24"  
Auto-Locking Mini Cart

*Shown with optional Prox Reader.*

*Beige shell with eleven drawer colors  
available (see page 103).*

*See page 80 for Auto-Locking details.*

#### Dimensions:

**Overall:** 34.60"H x 20.70"D x 21.40"W.

**Drawer size:** 17"D x 15"W.

**Drawer heights:** 3"H (3), 6"H, 9"H.

**Weight:** 70 lbs.

Four 3" swivel casters (two with brake  
and one tracking)

### Aluminum Six-Drawer 30" Auto-Locking Narrow Cart

AMC-6Color 6-Drawer 30"  
Auto-Locking Narrow Cart

*Beige shell with eleven drawer colors  
available (see page 103).*

*See page 80 for Auto-Locking details.*

#### Dimensions:

**Overall:** 43.40"H x 25.40"D x 25.30"W  
(includes bumper and casters).

**Drawer size:** 17"D x 15"W.

**Drawer heights:** 3"H (3), 6"H (2),  
9"H (1).

**Weight:** 85 lbs.

Four 5" swivel casters (two with brake  
and one tracking)

### Aluminum Five-Drawer 24" Auto-Locking Cart

PEL-Color-24 Five-Drawer 24"  
Aluminum Auto-Locking Cart

All the quality features of an Aluminum  
Premier™ Cart (except Slide-Thru  
Side Shelf)

*See page 80 for Auto-Locking details.*

*Beige shell with eighteen drawer colors  
available.*

**Overall:** 39.55"H x 25.20"D x 34.45"W  
(includes bumper, handles,  
and casters).

**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 22<sup>1</sup>/<sub>4</sub>"W.

**Drawer config.:** 24" total drawer height.

**Drawer heights:** 3"H (3), 6"H (1),  
9"H (1).

**Weight:** 88 lbs.

Four 5" swivel casters (two with brake  
and one tracking)

*Auto-Locking Carts Not Sold Or Serviced Outside The Continental USA.*

# Steel Push-Button Supreme Anesthesia Cart Package

PUSH-BUTTON LOCKING

## PUSH-BUTTON SUPREME ANESTHESIA PACKAGE

Complete Steel Cart with entire Accessory Package as shown  
APB-AA-6 with ASA-1X Anesthesia Package

ATB-4X  
Tilt Bin Organizer

ATB-5X  
Tilt Bin Organizer

ADO-1  
Divided  
Organizer  
with Lid

PSR-1  
Heavy-duty Railing

AMT-2N  
Modular Drawer  
Tray with  
Ampule Holder  
*(not pictured)*

ATW-3  
Soft-grip Handles

ACL-2  
Long Catheter  
Holder

AC-670S (4)  
3" Divisible Trays  
with Dividers  
(8 short, 4 long)

AC-672S (2)  
5" Divisible Trays  
with Dividers  
(4 short, 2 long)



Look for the  
Armstrong A-SMART®  
label on every cart.  
*It stands for quality.*

**MADE IN THE USA**



ATR-6N  
Trellis System

ATD-2  
Tape Dispenser

APT-2  
Plastic Top

ALH-12  
Label Holder,  
12/pkg., 7" x 2"  
*(not pictured)*

AWB-5  
Basket with  
AAB-1  
Accessory Bar

AUH-1  
Utility Hooks,  
2/pkg.  
*(not pictured)*

AWC-1  
Waste Container  
with AAB-1  
Accessory Bar

*All Drawers Have  
Full Extensions  
Ball Bearing Slides.*

NOTE: Carts require some assembly.




# Auto-Locking **Aluminum** Supreme Anesthesia Cart Package

AUTO-LOCKING  
Not MRI Compatible

## **AUTO-LOCKING SUPREME ANESTHESIA PACKAGE**

Complete Aluminum Cart with entire Accessory Package as shown  
PEL-AA-30 with PEL-PR Prox Reader & PASA-1X Anesthesia Package



# Standard Anesthesia Accessory Packages

KEY-LOCKING

## These Accessory Packages Mount To Any Standard Steel Cart

### Deluxe Anesthesia Accessory Package

DAB-1 Deluxe Accessory Package (*Does not include cart*)

Complete Package as shown, DAB-1 and APB-AA-6 6-Drawer 30" Push-Button Locking Cart

- ATR-6N Trellis System
- ATB-4X Four-bin Tilt Bin Organizer
- ATB-5X Five-bin Tilt Bin Organizer
- PSR-1 Tubular Railing
- ASN-3 Collapsible Side Shelf
- AWC-1 Waste Container
- AAB-1 Accessory Bar
- APT-2 Plastic Top



PIV-1 I.V. Pole with Brackets, *see page 113*



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**



APB-AA-6

*Bedside Carts, see pages 82-83*

### Economy Anesthesia Accessory Package

SP-1B Economy Accessory Package (*Does not include cart*)

Complete Package as shown, SP-1B and APB-B-6 6-Drawer 30" Push-Button Locking Cart

- APT-2 Plastic Top
- AUH-1 Utility Hooks, 2/pkg.
- PSR-1 Tubular Railing
- Two AC-672S 5"H Drug Trays with four Short Dividers and two Long Dividers (fills one drawer)
- Four AC-670S 3"H Drug Trays with eight Short Dividers and four Long Dividers (fills two drawers)



Modular Drawer Trays, *see page 118*



APB-B-6

NOTE: Carts require some assembly.



## Accessory Packages To Lighten The Load. . .


**Think** LIGHT.

### Deluxe Anesthesia Accessory Package

PDAB-1 Premier Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, PDAB-1 and PEL-SB-30 6-Drawer 30" Aluminum Auto-Locking Cart

- PNS-1 Shelf Unit
- ATB-5X Five-bin Tilt Bin Organizer
- AGB-1 Glove Box
- PSB-1 Sharps Container Bracket
- AWC-1 Waste Container



### Heine LED Laryngoscope Handle

see page 46



Look for the Armstrong A-SMART® label on every cart. **It stands for quality.**

**MADE IN THE USA**



PEL-SB-30 with PEL-PR



PEL-CW-30 with PEL-PR

### Economy Anesthesia Accessory Package

PSP-1B Premier Economy Accessory Package *(Does not include cart)*

Complete Package as shown, PSP-1B and PEL-CW-30 6-Drawer 30" Aluminum Auto-Locking Cart

- AGB-1 Glove Box
- AWC-1 Waste Container
- Four AC-670S 3" Drug Trays with eight Short Dividers and four Long Dividers
- Two AC-672S 5" Drug Trays with four Short Dividers and two Long Dividers



Premier™ Locking Bars, see page 108

See page 73 for Premier™ Cart features and page 75 for dimensions.

# Standard **Steel** Solid Color Mobile Workstations

KEY-LOCKING

## Solid Colors, Solid Features

All A-SMART® Key Locking Carts come with double side-wall construction, ball bearing drawer slides, a standard stabilizing frame/bumper, and four all-swivel casters (two locking and one tracking).



AA-6

### 6-Drawer 30" Standard Mobile Workstations

AA-6 Mobile Workstation,  
Deep Blue

Drawer sizes: 3"H (3), 6"H (2), 9"H



AN-6

### 6-Drawer 30" Standard Mobile Workstations

AN-6 Mobile Workstation,  
Beige

Drawer sizes: 3"H (3), 6"H (2), 9"H



AN-5

### 5-Drawer 24" Standard Mobile Workstations

AN-5 Mobile Workstation,  
Beige

Drawer sizes: 3"H (3), 6"H, 9"H



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**

*Want more accessories? See pages 110-121*



PIV-1 I.V. Pole with Brackets, see page 113



ATR-6N Standard Trellis Rail System, see page 112



ASN-3 Side Shelf, see page 113



APT-2 Standard Plastic Top, see page 113

Standard Anesthesia Packages, pages 84, 86  
Cart Accessories, pages 110-121

See page 72 for Standard Cart features and page 74 for dimensions.



# Premier™ Aluminum Solid Color Mobile Workstations

KEY-LOCKING  
Not MRI Compatible

## Premier™ Carts Are Easier To Push and STOP. . .

### Think *LIGHT*.



PKL-B-27

#### 6-Drawer 27" Premier Mobile Workstation

PKL-B-27 Mobile Workstation,  
Beige

Drawer sizes: 3"H (4), 6"H, 9"H



AEF-6 Plastic Divider Set  
for 6" drawer, see page 120



PKL-B-24

#### 5-Drawer 24" Premier Mobile Workstation

PKL-B-24 Mobile Workstation,  
Beige

Drawer sizes: 3"H (3), 6"H, 9"H



PKL-B-30

#### 6-Drawer 30" Premier Mobile Workstation

PKL-B-30 Mobile Workstation,  
Beige


Drawer sizes: 3"H (3), 6"H (2), 9"H



Tilt Bins, see page 117

Visit Our Website At  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

#### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Key locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Premier™ Anesthesia Packages, pages 85, 87

Cart Accessories, pages 110-121

# Standard **Steel** Two-Tone Mobile Workstations

KEY-LOCKING

## Colors To Match Your Decor

When your colors are added to the quality construction of an Armstrong A-SMART® Cart, you've got a Mobile Workstation like no other. **A-SMART®**



AKL-BG-4

### 4-Drawer 30" Standard Mobile Workstation

AKL-BG-4 Mobile Workstation, Blue Grey

Drawer sizes: 3"H, 9"H (3)



AKL-SB-6

### 6-Drawer 30" Standard Mobile Workstation

AKL-SB-6 Mobile Workstation, Slate Blue

Drawer sizes: 3"H (3), 6"H (2), 9"H



AKL-MG-5

### 5-Drawer 24" Standard Mobile Workstation

AKL-MG-5 Mobile Workstation, Moss Gray

Drawer sizes: 3"H (3), 6"H, 9"H



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**



Modular Drawer Trays, see page 118



Plastic Divider Sets, see page 120



ClearView Security™ Drug Box and Shelf, see page 112

### 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain

Standard Anesthesia Packages, pages 84, 86


Mobile Workstation Accessories, pages 110-121



# Premier™ Aluminum Two-Tone Mobile Workstations

KEY-LOCKING  
Not MRI Compatible

## 18 Designer Colors



## Maneuverability At It's Best. . .

**Think LIGHT.**



### 6-Drawer 27" Premier Mobile Workstation

PKL-CB-27 Mobile Workstation, Taupe

Drawer sizes: 3"H (4), 6"H, 9"H

### 6-Drawer 30" Premier Mobile Workstation

PKL-HG-30 Mobile Workstation, Hunter Green

Drawer sizes: 3"H (3), 6"H (2), 9"H

### 4-Drawer 30" Premier Mobile Workstation

PKL-Y-4 Mobile Workstation, Yellow

Drawer sizes: 3"H, 9"H (3)



Look for the Armstrong A-SMART® label on every cart. **It stands for quality.**

**MADE IN THE USA**


### 5-Drawer 24" Premier Mobile Workstation

PKL-SB-24 Mobile Workstation, Slate Blue

Drawer sizes: 3"H (3), 6"H, 9"H



### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Key locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Premier™ Anesthesia Packages, pages 85, 87

Cart Accessories, pages 110-121

# Standard **Steel** Premium Emergency Cart Package

BREAKAWAY LOCKING

## Standard Premium Emergency Package

AR-6 with APE-1X Emergency Package





# Premier™ Premium Aluminum Emergency Cart Package

BREAKAWAY LOCKING  
Not MRI Compatible

## Premier™ Premium Emergency Package

PBL-AR-30 with PAPE-1X Emergency Package

### Double Side-Wall Construction

PDS-1  
Shelf Unit

ATW-3  
Soft-grip Handles

ADV-10N  
Modular Dividers  
for ADT-1

ADT-1  
Drug Tray  
with Lid

AC-670S (4)  
3" Drug Trays  
with Dividers  
(8 short, 4 long)

AC-672S (2)  
5" Drug Trays  
with Dividers  
(4 short, 2 long)

ASU-1  
Suction Unit Holder



PCB-1  
Cardiac Board w/Brackets

PIV-1  
I.V. Pole

APS-5  
Plastic Seals  
(not pictured)

APT-5  
Plastic Top

Slide-thru  
Side Shelf

ALH-12  
Label Holders

AOB-1  
Oxygen Tank  
Bracket

AOS-1  
Hospital Grade  
Outlet Strip

PBR-30  
Breakaway Bar



Look for the Armstrong  
A-SMART® label on every  
cart. *It stands for quality.*

**MADE IN THE USA**

# Standard Emergency Cart Accessory Packages

## BREAKAWAY LOCKING

### Deluxe Emergency Accessory Package

EB-1 Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, EB-1 with AR-6 6-Drawer 30" Emergency Cart

- ASN-3 Collapsible Side Shelf
- PIV-1 I.V. Pole with Brackets
- APT-2 Plastic Top
- PSR-1 Heavy-duty Railing
- AOB-1 Oxygen Tank Bracket
- ACB-1 Cardiac Board with Brackets
- APS-5 Plastic Breakaway Seals, 100/pkg.
- AUH-1 Utility Hooks, 2/pkg
- Two AC-672S 5"H Trays with four AC-685S (short) Dividers & two AC-690S (long) Dividers (fills one drawer)
- Four AC-670S 3"H Trays with eight AC-675S (short) Dividers and four AC-680S (long) Dividers (fills two drawers)



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**



See page 74 for Standard Cart features and page 76 for dimensions.

NOTE: Carts require some assembly.

### Economy Emergency Accessory Package

BE-1 Economy Accessory Package *(Does not include cart)*

Complete Package as shown, BE-1 with ABL-AR-6 6-Drawer 30" Emergency Cart

- PIV-1 I.V. Pole with Bracket
- PSR-1 Heavy-duty Railing
- ACB-1 Cardiac Board with Brackets
- APT-2 Plastic Top
- APS-5 Plastic Breakaway Seals, 100/pkg.
- AOB-1 Oxygen Tank Bracket



ADS-IN, Defib One Shelf with Swivel Top, page 116



AE-6976 Cart-mounted DUET with Retention Bracket *(requires AHP-3), see page 39*



**BREAKAWAY LOCKING**  
Not MRI Compatible

## Deluxe Emergency Accessory Package

PEB-1 Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, PEB-1 with PAR-30 6-Drawer 30" All Aluminum Emergency Cart

- PIV-1 I.V. Pole with Brackets
- PSR-1 Tubular Steel Railing
- AOB-1 Oxygen Tank Bracket
- PCB-1 Cardiac Board with Brackets
- APS-5 Breakaway Seals, 100/pkg.
- Two AC-672S 5"H Trays with four AC-685S (short) Dividers and two AC-690S (long) Dividers
- Four AC-670S 3"H Trays with eight AC-675S (short) Dividers and four AC-680S (long) Dividers



PCB-1 Cardiac Board



ASU-1 Suction Unit Holder, see page 113

Mount either of these accessory packages on any A-SMART® Premier™ Aluminum Emergency Cart



PAR-30

## Economy Emergency Accessory Package

PBE-1 Economy Accessory Package *(Does not include cart)*

Complete Package as shown, PBE-1 with PBL-AB-30 6-Drawer 30" All Aluminum Emergency Cart

- PIV-1 I.V. Pole with Brackets
- PCB-1 Cardiac Board with Brackets *(shown above)*
- AOB-1 Oxygen Tank Bracket
- APS-5 Breakaway Seals, 100/pkg.



Look for the Armstrong A-SMART® label on every cart. ***It stands for quality.***

**MADE IN THE USA**

See page 73 for Premier™ Cart features and page 75 for dimensions



PBL-AB-30

# Standard **Steel** Solid Color Emergency Carts

BREAKAWAY LOCKING

## Manufactured To ISO 9001:2008 Certified Standards. *A-SMART*®.



ABL-B-6

### 6-Drawer 30" Standard Emergency Cart

ABL-B-6 Emergency Cart, Beige  
Drawer sizes: 3" (4), 6", 12"



ABL-B-5

### 5-Drawer 24" Standard Emergency Cart

ABL-B-5 Emergency Cart, Beige  
Drawer sizes: 3" (3), 6", 9"



ABL-B-4

### 4-Drawer 24" Standard Emergency Cart

ABL-B-4 Emergency Cart, Beige  
AR-4 Emergency Cart, Red  
*(not pictured)*  
Drawer sizes: 3" (2), 6", 12"



AR-6

### 6-Drawer 30" Standard Emergency Cart

AR-6 Emergency Cart, Red  
Drawer sizes: 3" (4), 6", 12"



AR-4

### 4-Drawer 24" Standard Emergency Cart

AR-4 Emergency Cart, Red  
Drawer sizes: 3" (2), 6", 12"



AR-27

### 6-Drawer 27" Standard Emergency Cart

AR-27 Emergency Cart, Red  
Overall Dim.: 42.20"H x 25.40"D x  
34.25"W (includes bumper, handles,  
and casters)  
Drawer sizes: 3" (4), 6", 9"  
Weight: 160 lbs.

See page 72 for Standard Cart features  
and page 74 for dimensions.

NOTE: Carts require some assembly.



# Premier™ Aluminum Solid Color Emergency Carts

BREAKAWAY LOCKING

Not MRI Compatible

## Drawers Always Slide Smoothly With Our Double Side-Wall Construction. . . Think **LIGHT.**



PBL-B-27

### 6-Drawer 27" Premier Emergency Cart

PBL-B-27 Emergency Cart, Beige  
Drawer sizes: 3"H (4), 6"H, 9"H



PBL-B-24

### 5-Drawer 24" Premier Emergency Cart

PBL-B-24 Emergency Cart, Beige  
Drawer sizes: 3"H (3), 6"H, 9"H



PAR-30

### 6-Drawer 30" Premier Emergency Carts

PAR-30 Emergency Cart, Red  
PBL-B-30 Emergency Cart, Beige (not pictured)  
Drawer sizes: 3"H (3), 6"H (2), 9"H



AC-670S and AC-672S, see page 120




PSR-1 Tubular Railing, see page 113



Look for the Armstrong A-SMART® label on every cart.  
**It stands for quality.**

**MADE IN THE USA**

### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Premier™ Emergency Packages, pages 93, 95

Cart Accessories, pages 110-121

# Standard **Steel** Two-Tone Emergency Carts

BREAKAWAY LOCKING

**Armstrong A-SMART® Two-Tone Carts**  
come with your choice of 18 colored drawer fronts.

*18 Designer Colors*



ABL-L-5

### 5-Drawer 24" Standard Emergency Cart

ABL-L-5 Emergency Cart, Lavender

Drawer sizes: 3" (3), 6", 9"



ABL-T-4

### 4-Drawer 24" Standard Emergency Cart

ABL-T-4 Emergency Cart, Teal

Drawer sizes: 3" (2), 6", 12"



ABL-TC-6

### 6-Drawer 30" Standard Emergency Cart

ABL-TC-6 Emergency Cart, Terra Cotta

Drawer sizes: 3" (4), 6", 12"



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**

### 3-Drawer 24" Standard Emergency Cart

ABL-B-3 Emergency Cart, Beige

Drawers: 6"H, 9"H (2)



ABL-B-3

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain



**AHP-3 Panel for Cart-mounted DUET Suction Unit**  
(Retention Bracket not included),  
see page 114



**AE-6976 Cart-mounted DUET Suction Unit with Retention Bracket**  
(requires AHP-3), see page 39



**ASN-3 Standard Side Shelf**,  
see page 113

Standard Emergency Packages, pages 92, 94

Emergency Cart Accessories, pages 110-121




# Premier™ Aluminum Two-Tone Emergency Carts

BREAKAWAY LOCKING  
Not MRI Compatible

## Carts Built To Last Longer!

Think *LIGHT*.

### 18 Designer Colors



PBL-AR-30

### 6-Drawer 30" Premier Emergency Cart

PBL-AR-30 Emergency Cart, Red  
Drawer sizes: 3" (3), 6" (2), 9"



PBL-SG-27

### 6-Drawer 27" Premier Emergency Cart

PBL-SG-27 Emergency Cart, Sage Green  
Drawer sizes: 3" (4), 6", 9"



PBL-SB-24

### 5-Drawer 24" Premier Emergency Cart

PBL-SB-24 Emergency Cart, Slate Blue  
Drawer sizes: 3" (3), 6", 9"



AMT-2N Modular Tray,  
see page 118



ASU-1 Suction Unit Holder  
(shown with AE-6975 DUET Suction Unit), see page 113



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**

### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Premier™ Emergency Packages, pages 93, 95

Cart Accessories, pages 110-121

# Standard **Steel** Push Button Cart

**PUSH-BUTTON**  
Not MRI Compatible

## Security With The Push Of A Button



APB-CB-6

### Standard **Steel** Push-Button Locking Carts

#### APB-Color-6 Standard Six-Drawer 30" Push-Button Locking Cart

Drawer sizes: 3"H (3), 6"H (2), 9"H.  
Weight: 170 lbs.

#### APB-Color-5 Standard Five-Drawer 24" Push-Button Locking Cart (*below*)

Drawer sizes: 3"H (3), 6"H, 9"H. Weight: 150 lbs.

This system comes with a key override to bypass Push-Button System if needed.

18 colors available.

**NOTE:** Push-Button control knob must be turned manually to lock and unlock.

APB-CB-6 shown at left with optional APT-2 Plastic Top and PSR-1 Tubular Railing.



PSR-1 Tubular Railing, see page 113



PIV-1 IV Pole with Brackets, see page 113



APB-B-5

### 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain




# NEW Aluminum ELECTRONIC Touchpad Carts

ELECTRONIC TOUCHPAD LOCKING  
Not MRI Compatible

## NEW Aluminum ELECTRONIC Touchpad Cart

### 18 Designer Colors



### PPE-Color-30 Six-Drawer 30" Aluminum Electronic Touchpad Cart

Drawer sizes: 3"H (3), 6"H (2), 9"H

Weight: 92 lbs.

Overall dimensions: 46.34"H x 25.24"D x 34.45" W

#### Features:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Plastic top
- Accessory panels (3)
- Stabilizing frame with bumper
- Aluminum mounting tracks
- 18 colors available

**NOTE: No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.**



*\* Please note for the Armstrong Electronic Touchpad Carts, we suggest using alcohol only for the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.*



PPE-AA-30 shown with optional PSR-1 Railing, ATD-2 Tape Dispenser, PTR-1 Trellis System, and ATB-4X & ATB-5X Tilt Bins

### PPE-Color-24 Five-Drawer 24" Aluminum Electronic Touchpad Cart

Drawer sizes: 3"H (3), 6"H, 9"H

Weight: 78 lbs.

Overall dimensions: 40.42"H x 25.24"D x 34.45" W

#### Features:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Plastic top
- Accessory panels (3)
- Stabilizing frame with bumper
- Aluminum mounting tracks
- 18 colors available

**NOTE: No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.**

# NEW Aluminum ELECTRONIC Touchpad Wide Cart & ScopeSafe™

ELECTRONIC TOUCHPAD LOCKING AND KEY-LOCKING



AKE-B-6

## NEW Aluminum ELECTRONIC Touchpad WIDE Cart

AKE-B-6 6-Drawer 30" Electronic Touchpad Wide Cart, Beige

A-SMART® Wide Carts feature:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Extra-wide drawers: 7" wider than standard cart drawers
- Seamless, one-piece plastic top
- Standard stabilizing frame with bumper
- 5" all-swivel casters (two locking, one tracking)
- Beige drawer fronts in beige frame
- Carts weigh 118 lbs.

**Dimensions:**

**Overall size:** 45.73"H x 25.24"D x 40.70"W (includes bumper, handles, and casters).

**Drawer size:** 17<sup>1</sup>/<sub>16</sub>"D x 29<sup>1</sup>/<sub>4</sub>"W.

**Actual Top Dimensions:** 32<sup>1</sup>/<sub>2</sub>"W x 17"D

**Available drawer configuration:**

**Six-Drawer:** 27<sup>7</sup>/<sub>8</sub>"H (3), 6<sup>1</sup>/<sub>8</sub>"H (2), 9<sup>3</sup>/<sub>8</sub>"H (1).

**NOTE: No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.**

A-SMART WIDE CART SYSTEMS ACCEPT ALL THE ACCESSORIES THAT FIT ON THE SIDES OF OUR STANDARD MOBILE WORKSTATIONS. SEE PAGES 110-121 FOR MORE OPTIONS.

*\* Please note for the Armstrong Electronic Touch Pad Carts, we suggest using alcohol only for cleaning the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.*



ATR-9N


ATR-9N Wide Cart Trellis Rail System, \$195.00

ATR-9X Additional Wide Cart Hanging Bar, \$42.00

ANS-9 Narrow Shelf Unit for Wide Cart, \$330.00

Dim.: 10" x 34<sup>1</sup>/<sub>2</sub>". Adds up to 21" of height to cart.

Wt. capacity: 50 lbs.



ANS-9

## ScopeSafe™

**Key-Locking Storage for Delicate Scopes**

ASH-4 ScopeSafe™ Locking Scope Holder, \$805.00

- Holds four scopes
- Has a padded interior for cushioned protection
- Contains multiple vents for ample ventilation
- Contains a protective plexiglas door front
- Adjustable hooks and brackets
- Attaches easily to wall or cart
- Built with a durable steel frame
- Overall size: 48"H x 6"D x 16"W
- Weight: 44 lbs.

ASH-2 ScopeSafe™ Locking Scope Holder, \$665.00 (not pictured)

- Holds two scopes
- Same high-quality features as ASH-4
- Compact: 48"H x 6"D x 9"W
- Weight: 31 lbs.
- Adjustable hooks and brackets



ASH-4

*Takes two Accessory Panels to mount to a Premier Cart. Two ASH-2 ScopeSafes can be mounted to one side of a cart. **Not compatible with Mini or Narrow Carts.** Please note: ASH-2 and ASH-4 will not hold colonoscopes.*




# Aluminum Mini and Narrow Auto-Locking Carts

**AUTO-LOCKING**  
Not MRI Compatible


## Spacesavers

### That Give You The Added Benefit Of Controlled Security

#### 11 Designer Colors



SG - Sage Green



MG - Moss Grey



TC - Terra Cotta



AA - Deep Blue



B - Beige



BG - Blue Grey



HG - Hunter Green



T - Teal



Y - Yellow



CB - Taupe



CW - Cherry Woodgrain

### Aluminum 24" Auto-Locking Mini Cart

(with *Optional Factory Installed Proximity Reader*)

AMC-4Color 5-Drawer 24" Auto-Locking Mini Cart

PEL-PR Optional Proximity Reader

Overall dimensions: 34.60"H x 20.70"D x 21.40"W.

#### Features:


- Keyless Entry
- 11 colors available
- Lightweight aluminum construction
- Master code for all programming
- Access code for normal, daily use
- Up to 250 user codes possible
- Manual or automatic locking
- Auto-Lock Timer can be set from 1 to 998 minutes, or it can be shut off so that the cart never auto-locks

- Battery operated; never needs to be plugged into an outlet
- "Low Battery" indicator light
- Manual key override for emergency access

Simply wave the card in front of the reader to open. Lock the cart by pressing the lock button or set the Auto-Lock Timer for anywhere from 1 to 998 minutes. The Proximity Reader will read most low frequency HID® Prox Cards (125 kHz).

**\*Please note for the Armstrong PEL and WEL carts, we suggest using alcohol only for cleaning the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.**

AMC-4SG  
with *optional*  
PEL-PR  
Prox Reader



AMC-6CB with *optional*  
PEL-PR Prox Reader,  
AGB-1 Glove Box Holder,  
and AWC-1 Waste  
Container and ATW-3  
Push Handles



### Aluminum 30" Auto-Locking Narrow Cart

(with *Optional Factory Installed Proximity Reader*)

AMC-6Color Six-Drawer 30" Auto-Locking Narrow Cart

PEL-PR Optional Proximity Reader

Overall: 43.40"H x 25.40"D x 25.30"W (includes bumper and casters).

### Optional Factory Installed Proximity Reader

PEL-PR Proximity Reader

PEL-PC HID® Prox Card

PEL-PS Proximity Card Sticker (*makes non-compatible prox cards compatible*)

The Proximity Reader is a factory installed option on the Auto-Locking Carts only. This upgrade must be noted at the time of purchase. The Proximity Reader can not be retrofitted on existing carts. Electronic capabilities include all the features of the Auto-Locking Carts.

**NOTE: If you wish to use your current HID® cards, or purchase new cards or stickers, each card must be made available to be programmed individually into each cart.**

*Auto-Locking Cart Software, see page 81.*

*Auto-Locking Carts Not Sold Or Serviced Outside The Continental USA.*

# STEEL Mini Carts

PUSH-BUTTON LOCKING AND KEY LOCKING



## Spacesaver

### Steel Two-Tone Mini Carts

Ideal For Treatments, Procedures, & the Bedside. *Small enough to put under counters.*

Compact, five-drawer beige shell Mini Carts include a plastic top, ball bearing drawer slides, corner bumper guards, and 3" swivel casters\* (two locking, one tracking). Mini Carts also features colored drawer fronts in a beige frame.


Overall: 34.50"H x 20.40"D x 20.90"W.

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (3), 6"H, 9"H.

Weight: 103 lbs.

#### 11 Designer Colors



AMC-1SG-PB



AMC-1MG



AMC-21CW-PB

### Steel Push-Button 5-Drawer Locking Mini Cart

AMC-1Color-PB 5-Drawer 24" Push-Button Locking Mini Cart

### Steel Key Locking 5-Drawer Mini Cart

AMC-1Color 5-Drawer 24" Key Locking Mini Cart

### Steel 4-Drawer Push-Button 21" Locking Mini Cart

AMC-21Color-PB 4-Drawer 21" Push-Button Locking Mini Cart

Overall: 31.30"H x 20.40"D x 20.90"W.

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (2), 6"H, 9"H.

Weight: 85 lbs.

*\*Optional 5" Casters with Stabilizing Frame, call customer service.*



# Aluminum Narrow Carts

PUSH-BUTTON LOCKING AND KEY LOCKING  
Not MRI Compatible


## Great For All Those Tight Spaces.

Think *LIGHT*.

### 11 Designer Colors



SG - Sage Green



MG - Moss Grey



TC - Terra Cotta



AA - Deep Blue



B - Beige



BG - Blue Grey



HG - Hunter Green



T - Teal



Y - Yellow



CB - Taupe



CW - Cherry Woodgrain

### Features:

- Lightweight aluminum construction; weighs only 85 lbs.
- Stabilizing frame with bumper
- Four 5" swivel casters (two with brake and one tracking)
- Plastic top
- Eleven colors available
- Ball bearing drawer slides

### Dimensions:

**Overall:** 43.40"H x 25.40"D x 25.30"W (includes bumper and casters).

**Drawer configuration:** 30" of total drawer height.

**Drawer size:** 17"D x 15"W.

**Drawer heights:** 3"H (3), 6"H (2), 9"H (1).

**Weight:** 85 lbs.

**Auto-Locking  
Narrow Cart,**  
see page 103



**Aluminum 30" Key  
Locking Narrow Cart**  
AMC-2Color Six-Drawer 30"  
Key Locking Narrow Cart



**Aluminum 30"  
Push-Button Locking  
Narrow Cart**

AMC-3Color Six-Drawer  
30" Push-Button Locking  
Narrow Cart

### Optional **MINI** & **NARROW** Cart Accessories

APH-1 Push Handle, \$60.00

MDM-1B Keyboard Drawer, \$120.00  
Weight capacity: 30 lbs.

ATR-7N Three-rail Trellis System,  
\$185.00

ATR-8N Three-rail Trellis System with  
Four Hooks, \$215.00

AMT-4N Full Drawer Tray with Nine  
Ampule Dividers, Ampule Tray, Rail,  
and Four Adjustable Dividers, \$90.00  
(Fits one per drawer), see page 118.

AMT-5N Full Drawer Tray with  
Two Rails and Four Adjustable  
Dividers, \$70.00 (Fits one per drawer),  
see page 118.

AEF-9M Plastic Divider Set for 9"  
Mini Cart Drawer, \$74.00  
Two 16<sup>3</sup>/<sub>8</sub>"L rails and four 14<sup>3</sup>/<sub>16</sub>"L  
x 8<sup>1</sup>/<sub>4</sub>"H dividers.

AEF-6M Plastic Divider Set for 6"  
Mini Cart Drawer, \$70.00  
Two 16<sup>3</sup>/<sub>8</sub>"L rails and four 14<sup>3</sup>/<sub>16</sub>"L  
x 5<sup>1</sup>/<sub>2</sub>"H dividers.

AEF-3M Plastic Divider Set for 3"  
Mini Cart Drawer, \$54.00  
Two 16<sup>3</sup>/<sub>8</sub>"L rails and four 14<sup>3</sup>/<sub>16</sub>"L  
x 2<sup>1</sup>/<sub>2</sub>"H dividers.

See pages 110-121 for more  
accessories.

# Broselow® Pediatric Resuscitation System

## Broselow® Pediatric Resuscitation System



### Broselow® Pediatric Resuscitation System

AE-4700 Broselow® Pediatric Resuscitation System, \$1,950.00

Includes twenty-eight size-specific procedural modules in seven pouches, flying carpet, Broselow® Pediatric Emergency Tape, and all other equipment shown at left. All contained in a heavy-duty ALS equipment organizer.

The Broselow® Pediatric Resuscitation System is approved for the PALS education program. For children 3 to 36 kg.

### NEW Broselow® Pediatric Emergency Tape “2017 Edition” mLs units are back!

AE-4800 Broselow® Pediatric Emergency Tape, 5/pkg., \$120.00

The 2017 version of the Broselow Reference Tape will be compliant with the 2015 PALS Standards

- Pre-calculated doses in milligrams, as well as all doses in pre-calculated ml, making the entire process error free by removing the need for any calculations.
- The new tape will accommodate for obesity by allowing for a habitus adjustment if deemed appropriate to the tape weight estimation.

#### 1. Measure

Measure the child with the Broselow® Tape. This gives you an accurate weight estimate, equipment sizes, and drug dosages.

#### 2. Treat


Reach into the bag for exactly the right equipment for I.V., intubation, O<sub>2</sub> delivery, or intraosseous access. There's no guessing.

See page 41 for more features



### NEW Broselow® Pediatric Resuscitation Medication/Infusion Guide “2017 Edition”

AE-4875 Broselow® Pediatric Resuscitation Medication/Infusion Guide, \$210.00



This Resuscitation/Infusion Guide is designed to provide you a quick reference to common resuscitation I.V. push drugs and vasoactive infusion drugs. The guide is sectioned off in colors that correspond with the Broselow® Tape. The guide is also a solution for those who are trying to address The Joint Commission mandate to use standard concentration drugs for infusions and eliminate the use of “Rule of 6” methods of calculating infusions. The front section provides precalculated

milliliter amounts of the most common I.V. push drugs for resuscitation. The back section provides infusion rates based on standard concentrations of common vasoactive medications. This Guide acts to further enhance the Broselow Tape. The Broselow® Pediatric Emergency System helps to reduce cognitive overload sometimes associated with pediatric resuscitations, as well as reduces the potential for medication mistakes due to calculation errors.

### STOCK YOUR CART

Fill the drawers with the Broselow® Pediatric Resuscitation System, size-specific pediatric equipment for a complete resuscitation system.

AE-4712 Procedural Modules, \$920.00

A complete set of twenty-eight modules for the seven Broselow® Color Packs or the Broselow® ColorCode™ Cart.

Modules for Intraosseous Access, O<sub>2</sub> Delivery, Intubation, and I.V. Delivery.

AE-4703 Equipment Pack, \$655.00

Contains all the equipment that comes in the original Broselow® bag: two disposable manual resuscitators, four cuff masks, three disposable blood pressure cuffs, one gauge, two Magill forceps, and two Vital View II™ LED fiber optic laryngoscope handles.



# Premier™ Aluminum Broselow® ColorCode™ Cart

BREAKAWAY LOCKING  
Not MRI Compatible

*Use in conjunction with your  
Broslow® Pediatric Resuscitation System.*

## Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks



PBR-3 Three-inch Breakaway Locking Bar, \$45.00

PBR-6 Six-inch Breakaway Locking Bar, \$45.00

Alternative locking options help simplify inventory or limit access to sensitive materials. See page 108 for more locking options.



PBL-PC-9



Use OPTIONAL 3" drawer locks so that only the accessed drawers need restocking.

PBL-PC-9 shown with OPTIONAL Accessory Package PC-9E:

Includes AOB-1 Oxygen Tank Bracket, ASU-1 Suction Unit Holder, PIV-1 IV Pole, PCB-1 Cardiac Board (not shown), PBR-3 (8) 3" Breakaway Locking Bar, PBR-6 (1) 6" Breakaway Locking Bar, and APS-5 Plastic Seals. Also shown, AE-4800 Broslow® Pediatric Emergency Tape (sold separately).

PBL-PC-9 BASIC A-SMART® Premier™ Broslow® ColorCode™ Cart (pictured left)

Includes Nine-Drawer A-SMART Premier Breakaway Lock Cart (30" total drawer space), slide-thru side shelf, 30" locking bar, stabilizing frame/bumper system, seamless plastic top, soft-grip push handles, three accessory attachment panels, and all-swivel casters (two locking and one tracking).

PBL-PC-9A COMPLETE A-SMART Premier™ Broslow® ColorCode™ Cart System (not shown)

System includes PBL-PC-9 Basic ColorCode Cart, as well as AE-4712 Broslow® Procedural Modules, AE-4703 Broslow® Equipment Pack, and a Broslow® Pediatric Emergency Tape (see pages 41-45).

Premier™ Emergency Packages, pages 85, 87

Cart Accessories, pages 110-121

## Restock Only The Accessed Drawers. Save Money. Save Time.

Don't waste money checking crash cart drawers that haven't even been opened. With A-SMART® Premier™ Carts, individual drawer seals can be arranged so that only the accessed drawers need to be checked for inventory, a small investment that pays off with every use.

- Double lock medication drawers
- Improve equipment security
- Ready for changing needs

When carts are shared by departments or personnel, A-SMART Premier's locking options allow each to be locked with a

separate key. Never wonder about who else has access to your equipment.

Maybe you don't need this kind of flexibility now, but you might later. With A-SMART Premier, you can update your locking system at any time. Need to change a mobile workstation into a crash cart? No problem. Need to add a key lock to an existing cart? Done. Only A-SMART Premier offers you that kind of flexibility.

Locks may be placed on either corner glide for flexible security. **Double locking** drawers can help control inventory and access.



### Breakaway Locking Bars

Mix lock sizes and types for precise inventory control.

Size	Prod#	Price
3"	PBR-3	\$45.00
6"	PBR-6	\$45.00
9"	PBR-9	\$45.00
15"	PBR-15	\$45.00
21"	PBR-21	\$45.00
24"	PBR-24	\$45.00
27"	PBR-27	\$45.00
30"	PBR-30	\$45.00

### Key Locking Bars

Different sized locking bars help control access.

Size	Prod#	Price
3"	PKE-3	\$65.00
6"	PKE-6	\$65.00
9"	PKE-9	\$65.00
24"	PKE-24	\$75.00
27"	PKE-27	\$75.00
30"	PKE-30	\$75.00

### Drawer Specific

Drawer-specific locks simplify inventory, saving money in restocking time.



Look for the Armstrong A-SMART® label on every cart. **It stands for quality.**

**MADE IN THE USA**



# Two-Tone Aluminum & Steel Isolation Carts

KEY LOCKING & BREAKAWAY LOCKING  
Not MRI Compatible

## Isolation Carts



AKL-Y-400IP

### AKL-Y-400IP STANDARD Steel Key Locking 4-Drawer 30" Isolation Cart Package

Includes AKL-Y-4 Four-Drawer Cart, PSR-1 Tubular Railing, APT-2 Plastic Top, AUH-1 Utility Hooks, and ACR-1 Coat Rack.

For Standard Cart features, see page 72.  
For Standard Cart dimensions, see page 74.



PKL-Y-4

### PKL-Y-4 Basic PREMIER™ Aluminum Key Locking 4-Drawer 30" Isolation Cart

### PBL-Y-4 Basic PREMIER™ Aluminum Breakaway Locking 4-Drawer 30" Isolation Cart

Drawers: 3"H, 9"H (3)

For Premier™ Cart features, see page 73.  
For Premier™ Cart dimensions, see page 75.



AMC-1Y

### AMC-1Y 5-Drawer 24" Key Locking Steel Mini Cart

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (3), 6"H, 9"H.

Weight: 103 lbs.



AMC-2Y

### AMC-2Y 6-Drawer 30" Key Locking Aluminum Narrow Cart

Drawer size: 17"D x 15"W.

Drawer heights: 3"H (3), 6"H (2), 9"H (1).

Weight: 85 lbs.



AKL-Y-4

### AKL-Y-4 Basic STANDARD Steel Key Locking 4-Drawer 30" Isolation Cart

Drawers: 3"H, 9"H (3)

For Standard Cart features, see page 72.  
For Standard Cart dimensions, see page 74.



ABL-Y-3

### ABL-Y-3 Basic STANDARD Steel Breakaway Locking 3-Drawer 24" Isolation Cart

Drawers: 6"H, 9"H (2)

For Standard Cart features, see page 72.  
For Standard Cart dimensions, see page 74.

# Cart Extender and Glove Boxes



## Cart Extender

ACE-1 Armstrong Cart Extender\*, Front-open, \$560.00

- Adds 4.8 cu. ft. to cart
- Key lock to protect contents
- Includes three adjustable shelves
- Each shelf fits AC-670S or AC-672S drawer trays (one per shelf)
- Safety foot adds stability
- Shipping Weight: 52 lbs.
- Dimensions: 34<sup>3</sup>/<sub>8</sub>"H x 18"D x 13<sup>1</sup>/<sub>2</sub>"W


\*Cart Extenders only fit A-Smart Carts with 30" of drawer space. Does not fit on a Narrow, Mini, or Wireless Cart.

## Looking For Extra Wide Drawers?

See page 102 for  
**NEW** Aluminum ELECTRONIC Touchpad WIDE Cart



## Glove Dispensers



**AGB-1 Cart-mount Glove Box, \$55.00**  
5.75"W x 10"H x 3<sup>7</sup>/<sub>8</sub>"D.

**AGB-3 Steel, Baked-on Enamel Glove Box, \$85.00**  
16"W x 10"H x 3<sup>7</sup>/<sub>8</sub>"D.

**AE-5359 Glove Dispenser, 6/case, \$105.00**  
Wall-mountable glove dispenser in a light beige. Adjustable shelf fits different sized boxes of gloves.  
14<sup>1</sup>/<sub>2</sub>" x 6" x 4<sup>1</sup>/<sub>2</sub>".

**AGB-2P Plexiglas Glove Box, \$45.00**  
10"W x 10<sup>3</sup>/<sub>4</sub>"H x 3<sup>7</sup>/<sub>8</sub>"D.  
*Mounting hardware not included*



**AGB-2 Steel, Baked-on Enamel Glove Box, \$55.00**  
Holds 2 boxes of gloves. 10"W x 10<sup>3</sup>/<sub>4</sub>"H x 3<sup>7</sup>/<sub>8</sub>"D.  
*Mounting hardware not included*

**AGB-3P Plexiglas Glove Box, \$60.00 (not shown)**  
Holds 3 boxes of gloves. 16"W x 10"H x 3<sup>3</sup>/<sub>4</sub>"D  
*Mounting hardware not included*

**AGB-3V Steel, Baked-on Enamel Glove Box, \$75.00**  
Holds 3 boxes of gloves vertically.  
5<sup>1</sup>/<sub>2</sub>"W x 30"H x 4"D.  
*Mounting hardware not included*





## Auto-Locking Wall Cabinet with Prox Reader

### PEL-IC Auto-Locking Wall Cabinet with Prox Reader

The Auto-Locking Wall Cabinet with Prox Reader is perfect for extra security of items inside or outside of patients' rooms or nurses stations. With this cabinet, you get all the great locking features of our Auto-Locking Carts and Prox Reader (see pages 80-81).

- Aluminum construction
- Two adjustable shelves
- Key override
- Includes hardware for wall or base mounting
- Battery powered (6 "D" batteries)
- Optional AC adapter

Dimensions: 16<sup>1</sup>/<sub>2</sub>"W x 27"H x 12"D  
Weight: 20 lbs.

PEL-AC Optional AC Adapter



PEL-IC

*Security of items inside OR outside of patients' rooms*



PEL-CT

## Auto-Locking Countertop Cabinet with Prox Reader

### PEL-CT Auto-Locking Countertop Cabinet with Prox Reader

The Auto-Locking Countertop Cabinet with Prox Reader is perfect for extra security of items inside or outside of patients' rooms or nurses stations. With this cabinet, you get all the great locking features of our Auto-Locking Carts and Prox Reader (see pages 80-81).

- Aluminum construction
- 3" drawer

- Includes mounting hardware
- Key override
- Battery powered (6 "D" batteries)
- Optional AC adapter available

Dimensions: 18<sup>1</sup>/<sub>2</sub>"W x 7"H x 18"D  
Drawer size: 17"D x 15"W  
Weight: 18 lbs.

PEL-AC Optional AC Adapter

## Auto-Locking Countertop Cabinet Trays



AMT-4N

### AMT-4N Auto-Locking Countertop Cabinet Tray, \$90.00

Full drawer tray with nine ampule dividers, ampule tray, rail, and four adjustable dividers. Fits one per cabinet. 14<sup>3</sup>/<sub>16</sub>"W x 16"D x 2<sup>1</sup>/<sub>2</sub>"H

### AMT-5N Auto-Locking Countertop Cabinet Tray, \$70.00

Full drawer tray with two rails and four adjustable dividers. Fits one per cabinet. 14<sup>3</sup>/<sub>16</sub>"W x 16"D x 2<sup>1</sup>/<sub>2</sub>"H



AMT-5N

# Security Box and Trellis Systems


## ClearView Security™ Drug Box and Shelf

ACV-1PB Push-Button Locking ClearView Security Drug Box with Key Override, \$325.00

CVS-1 ClearView Shelf, \$80.00

APK-1 ClearView Privacy Cover Kit, \$30.00

ACV-1T ClearView Blue Tray, \$19.00



ACV-1PB

CVS-1

The ClearView provides a great way to secure your medications, small equipment, patients' personal valuables, RSI equipment, and so much more. The ClearView also allows for a clear view of your inventory and also allows for anesthesia cases to be setup in advance while still keeping instruments/drugs secure. The ClearView can be mounted to a wall, in a patient's room, or even on an A-SMART® Cart.

### Features:

- Natural HDPE cabinet with clear PVC door
- Corrosion-resistant
- Provides a clear view of your inventory

**Dimensions:** 7<sup>1</sup>/<sub>2</sub>"W x 14"D x 7<sup>3</sup>/<sub>8</sub>"H

**Weight:** 5 lbs



ACV-1PB with APK-1

*Want to keep prying eyes away? Just install the APK-1 ClearView Privacy Cover Kit.*



*Attach the ClearView Security Drug Box to an A-SMART® Cart with the CVS-1 ClearView Side Shelf.*

## Trellis Systems



PTR-1

**PTR-1 PREMIER™ Trellis System, \$130.00**

Adds up to 21" to cart height. Bars are height-adjustable. Includes two uprights and three cross bars.

**PTR-5 PREMIER™ Trellis System with Hooks, \$170.00 (not pictured)**

Adds up to 21" to cart height. Bars are height-adjustable. Includes two uprights and three cross bars. Includes one bar with five hooks in the top rail.



ATR-5N

**ATR-5N STANDARD Trellis Rail System with Five Hooks, \$225.00**

Same features as the ATR-6N (see right), but with five hooks on the top rail for hanging supplies.

**PHB-1 Additional Hanging Bar for Standard or Premier Trellis Systems, \$42.00**



ATR-6N

**ATR-6N STANDARD Trellis Rail System, \$185.00**

The top bar is 21<sup>1</sup>/<sub>2</sub>" above top of cart.

**ATR-9N Wide Cart Trellis Rail System, \$195.00**

Fits Armstrong's Wide AKE Series Carts (see page 102).





PSR-1


**PSR-1 Tubular Railing, \$135.00**

Made of electro-polished stainless steel. Extends 3<sup>1</sup>/<sub>4</sub>" above the top of the cart.

**APT-2 STANDARD Plastic Top, \$65.00**

Seamless ABS plastic protects cart metal and provides easy-to-clean work surface. Rounded corners mean infection has nowhere to hide.

Dimensions: 18<sup>1</sup>/<sub>2</sub>" x 27".



APT-2

**APT-5 PREMIER™ Plastic Top, \$65.00**

(not pictured)

Dimensions: 18<sup>1</sup>/<sub>2</sub>" x 27".

**PIV-1 I.V. Pole with Brackets, \$105.00**

Adjustable: extends 9<sup>1</sup>/<sub>2</sub>" to 39<sup>1</sup>/<sub>2</sub>" above cart top.



PIV-1

**ASN-3 STANDARD Side Shelf, \$160.00**

**Simple One-Hand Slow Shelf Retraction**

Adds 1.8 square feet of temporary work surface when extended and folds neatly away when not in use.

Dimensions: 17<sup>3</sup>/<sub>4</sub>"D x 14<sup>7</sup>/<sub>8</sub>"W.

Weight Limit: 25 lbs.



ASN-3

**PSN-3 PREMIER™ Side Shelf, \$160.00**

**Simple One-Hand Slow Shelf Retraction**

Same as ASN-3 (see above).



ACB-1

**ACB-1 STANDARD Plastic Cardiac Board with Brackets, \$82.00**

Heavy-duty plastic board with handle that mounts to back of cart for easy access.

Dimensions: 15.5" x 21.5" x .45".



PCB-1

**PCB-1 PREMIER™ Cardiac Board with Brackets, \$82.00**

Same as ACB-1 (see above).

**ASU-1 Suction Unit Holder, \$120.00**

Dimensions: 10"H x 19"W x 7<sup>3</sup>/<sub>4</sub>"D.



PSN-3



ASU-1

**APC-1 Armstrong Wire Punch Card Tray, White, \$40.00**

Will hold up to thirty punch cards per tray.

Dimensions: 7"H x 14"D x 7"W.

Great for fluid transfer.



APC-1

**ADO-1 Divided Organizer with Lid, \$110.00**

Dimensions: 2<sup>3</sup>/<sub>4</sub>"H x 7<sup>1</sup>/<sub>2</sub>"W x 7<sup>1</sup>/<sub>2</sub>"D.



ADO-1

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 114


# Cart Accessories



ASB-2

**ASB-2 Armstrong Security Box, \$155.00**

Lockable box for portable use or inside a 3"H drawer. Designed to fit beside an AC-670S Drug Tray. The security box measures 11"W x 16"D x 2<sup>5</sup>/<sub>8</sub>"H.



ATD-2

**ATD-2 Tape Dispenser for Trellis Rail, \$75.00**

includes two trellis brackets

**ATD-3 Tape Dispenser for Cart Top, \$55.00**

5 tape brackets, each bracket width is 2<sup>1</sup>/<sub>4</sub>" (not pictured)

**AWB-1 Wire Basket\*, \$45.00**

Plastic-coated wire basket. 15"L x 12"W. 4"H in back, 3"H in front. Sloping top, flat bottom.

**AWB-3 Wire Chart Holder\*, \$35.00**

Plastic-coated wire chart holder. 16"L x 4"D x 10"H.

**AWB-5 Sloping Wire Basket\*, \$35.00**

Plastic-coated wire basket measures 12" square. 8"H in back, 4"H in front. Sloping top, flat bottom.



AWB-3

**AWC-1 Waste Container\*, includes 2 ATB-1 brackets, \$33.00**

13<sup>5</sup>/<sub>8</sub> quarts. 11<sup>3</sup>/<sub>8</sub>" x 8<sup>1</sup>/<sub>4</sub>" x 12<sup>1</sup>/<sub>8</sub>".

**AWC-2 Fire Retardant Waste Container\*, includes 2 ATB-1 brackets, \$65.00**

14 quarts. 12<sup>1</sup>/<sub>4</sub>" x 8<sup>1</sup>/<sub>4</sub>" x 11<sup>1</sup>/<sub>2</sub>". (not pictured)

**AWC-3 Large Waste Container\*, includes 2 ATB-1 brackets, \$35.00**

28 quarts. 14<sup>3</sup>/<sub>8</sub>" x 10<sup>1</sup>/<sub>4</sub>" x 15". (not pictured)

**ACV-1PB Push-Button Locking ClearView Security Drug Box with Key Override, \$325.00**

**CVS-1 ClearView Shelf, \$80.00**

**PHP-1 PREMIER™ Accessory Panels, 2/pkg., \$40.00**

16" x 4".

**PHP-5 PREMIER™ Accessory Panels for Back of Cart, 2/pkg., \$50.00**

23" x 4".

**AAB-1 STANDARD Accessory Bar Bracket, \$21.00**

Required for hanging some accessories on this page. Attaches to two locations on either side of cart. Size: 18" x 1".

(For A-SMART Premier Carts, substitute **PHP-1** Accessory Panel).

**PHP-1S PREMIER™ Panel for Cart-mounted DUET, \$50.00**

(not pictured)

**AHP-3 STANDARD Panel for Cart-mounted DUET, \$40.00**

(not pictured)

\*Requires AAB-1 or PHP-1 to hang (see above)



ACV-1PB with CVS-1



PHP-1

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.



AND-1



AWB-1



AWB-5



AWC-1



AAB-1

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see this page.



# Cart Accessories



ACS-2

**\*ACS-2 Short Catheter Holder, \$103.00**

- Has metal construction with neutral beige finish
- Requires no screws to mount
- 19"H x 8"W x 4"D
- Divided into two compartments with non-removable divider

**\*ACS-3 Short Catheter Holder with Clear Front, \$98.00**

- Clear front enables you to see inventory
- Front panel removes for easy cleaning
- Requires no screws to mount
- 19"H x 4"W x 4"D



ACS-3



AOB-1

**\*ACL-2 Long Catheter Holder, \$145.00**

- Fits A-SMART® Carts with 30" drawer space
- Has metal construction with neutral beige finish
- Requires no screws to mount
- 28"H x 8"W x 4"D
- Divided into two compartments with non-removable divider

**\*ACL-3 Long Catheter Holder with Clear Front, \$102.00**

- Fits A-SMART Carts with 30" drawer space
- Clear front enables you to see inventory
- Front panel removes for easy cleaning
- Requires no screws to mount
- 28"H x 4"W x 4"D



ACL-2

**\*AOB-1 Oxygen Tank Bracket, \$70.00**

Fits "D" and "E" sized tanks.



AOS-1

**AFT-1 STANDARD Fluid Tray System for Cart with 24" Drawer Space, 3 trays/pkg., \$165.00**

**AFT-1P PREMIER™ Fluid Tray System for Cart with 24", 27", or 30" Drawer Space, 3 trays/pkg., \$190.00 (not pictured)**

**AFT-2 STANDARD Fluid Tray System for Cart with 30" Drawer Space, 3 trays/pkg., \$165.00**

Three height-adjustable trays on rail system. Inner tray dimensions: 17<sup>1</sup>/<sub>8</sub>"L x 4<sup>1</sup>/<sub>2</sub>"D x 4"H.



AFT-1

**AOS-1 Hospital-grade Outlet Strip, \$135.00**

Six outlets and 6-ft cord.

**AOS-15 Hospital-grade Outlet Strip, \$145.00**

Six outlets and 15-ft cord.



PSB-1

**\*PSB-1 PREMIER™ Sharps Container Brackets, \$32.00**

Fits Sharps containers that are 2 gal. or 3 gal.

**ASB-1N STANDARD Sharps Container Bracket, \$72.00**

Fits Sharps containers that are 2 gal. or 3 gal.




ASB-1N

*\*May require PHP-1 panels to mount to a Premier™ Cart.*

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 114

# Shelving Units



ACD-1N PCD-1

**ACD-1N STANDARD Defib Corner Shelf Unit, \$235.00** (not compatible with Narrow or Mini Carts)  
Same features as PCD-1 (see below).

**PCD-1 PREMIER™ Defib Corner Shelf Unit, \$260.00**

Adds 9<sup>1</sup>/<sub>2</sub>" to height of cart. Swivel top rotates 360°, locks every 45°. Adjustable width: 13" to 18". Depth: 12". Capacity: 25 lbs.



PDS-1

**PNT-2 PREMIER™ Two-shelf Unit, Narrow, \$375.00** (not pictured)

Dim.: 10<sup>1</sup>/<sub>4</sub>" x 27<sup>1</sup>/<sub>2</sub>". Adds up to 21"H. Units include two uprights and two shelves.



ABS-1N

**PNS-1 PREMIER™ One-shelf Unit, Narrow, \$265.00**

Dim.: 10<sup>1</sup>/<sub>4</sub>" x 27<sup>1</sup>/<sub>2</sub>". Adds up to 21"H to cart. Capacity: 55 lbs.

**PDS-1 PREMIER™ One-shelf Unit, Deep, \$285.00**

Dim.: 14<sup>5</sup>/<sub>8</sub>" x 27<sup>7</sup>/<sub>16</sub>". Adds up to 21"H to cart. Capacity: 55 lbs.



PNS-1

**PDT-2 PREMIER™ Two-shelf Unit, Deep, \$385.00** (not pictured)

Dim.: 14<sup>5</sup>/<sub>8</sub>" x 27<sup>7</sup>/<sub>16</sub>". Adds up to 21"H to cart. Units include two uprights and two shelves. Capacity: 55 lbs.



ATN-1N

**ABS-1N STANDARD One-shelf Unit, Deep, \$285.00**

Adds 20" to height of cart. Outer dimensions: 14<sup>5</sup>/<sub>8</sub>"D x 27<sup>7</sup>/<sub>16</sub>"W. Capacity: 75 lbs.

**ABS-2N STANDARD One-shelf Unit, Narrow, \$275.00**

Dim.: 10<sup>1</sup>/<sub>4</sub>" x 27<sup>1</sup>/<sub>2</sub>". Adds up to 20"H to cart. Capacity: 55 lbs.



ABS-2N

**ATN-1N STANDARD Two-shelf Unit, Deep, \$420.00**

Outer dimensions: 14<sup>5</sup>/<sub>8</sub>"D x 27<sup>7</sup>/<sub>16</sub>"W. Adds 23" to height of cart. Capacity: 75 lbs.

**ADS-1N STANDARD Armstrong Defib One-shelf Unit with Swivel Top, \$560.00**

Outer dimensions: 14<sup>1</sup>/<sub>16</sub>" x 26<sup>7</sup>/<sub>8</sub>". Adds 12" to height of cart. Swivel top can rotate 360°. Capacity: 50 lbs.



ADS-1N PDS-2

**PDS-2 PREMIER™ Armstrong Defib One-shelf Unit with Swivel Top, \$575.00**

Same as ADS-1N (see above).

## Cart Trellis Systems

See pages 112

All products are compatible with both Standard & Premier Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 114.



## Tilt Bin Organizers Maximize Storage

### Tilt Bins Organizers

Designed for modular assembly, Tilt Bins organize supplies to match your needs. Clear Tilt Bins leave stock visible for quick identification, yet are securely arranged.

Each bin has a slot for labeling and swings open for easy access. Bins are removable for cleaning or refilling.

**NOTE:** A Trellis Rail System or Shelving Unit (except ACD-1N, PCD-1, PLS-1, and ALS-1) is required to mount Tilt Bins to an A-SMART Cart. See page 112 for Trellis Rail Systems and page 116 for Shelving Units.

#### ATB-6X Six-Bin Unit, \$27.00

Each bin measures: 4<sup>1</sup>/<sub>2</sub>"H x 5<sup>1</sup>/<sub>4</sub>"D

#### ATB-5X Five-Bin Unit, \$32.00

Each bin measures: 6<sup>1</sup>/<sub>2</sub>"H x 5<sup>1</sup>/<sub>4</sub>"D

#### ATB-4X Four-Bin Unit, \$45.00

Each bin measures: 8<sup>1</sup>/<sub>8</sub>"H x 6<sup>5</sup>/<sub>8</sub>"D

#### ATB-3X Three-Bin Unit, \$60.00

Each bin measures: 9<sup>1</sup>/<sub>2</sub>"H x 7<sup>3</sup>/<sub>4</sub>"D

Each Tilt Bin Organizer has a total width of 22<sup>1</sup>/<sub>2</sub>".

### Tilt Bin Dividers


Divide each bin into two sections.

ATB-4D Dividers for ATB-4X, 4/pkg., \$15.00

ATB-3D Dividers for ATB-3X, 3/pkg., \$18.00



All products are compatible with both Standard and Premier™ Carts unless otherwise stated.



A combination of bin sizes can be used to accommodate multiple-sized supplies or equipment.

**Need More Organization  
For All Your Supplies?  
See pages 118-120.**



### Lockable Tilt Bin Organizers



**NOTE:** A Trellis Rail System or Shelving Unit (except ACD-1N, PCD-1, PLS-2, and ALS-2) is required to mount Tilt Bins to an A-SMART Cart. See page 112 for Trellis Rail Systems and page 116 for Shelving Units. Accessories (including Lockable Tilt Bins) can't be mounted directly above Lockable Tilt Bins.

*Tilt Bin dimensions shown above.*

#### Key Lock Tilt Bins

ATB-5XK Five-Bin Unit, \$115.00

ATB-4XK Four-Bin Unit, \$135.00

ATB-3XK Three-Bin Unit, \$145.00

#### Breakaway Lock Tilt Bins

ATB-5XL Five-Bin Unit, \$75.00

ATB-4XL Four-Bin Unit, \$90.00

ATB-3XL Three-Bin Unit, \$100.00

# Modular Drawer Trays

## A Place For Everything...



AMT-1N

**AMT-1N Modular Tray, \$75.00**  
Full drawer tray includes two vertical rails and eight adjustable dividers with card holder. Dimensions: 21<sup>1</sup>/<sub>2</sub>"W x 15<sup>1</sup>/<sub>4</sub>"D x 2<sup>5</sup>/<sub>16</sub>"H.



AMT-6N

**AMT-2N Modular Tray, \$100.00**  
Same as the AMT-1 (see above), except with an ampule organizer, eleven adjustable dividers, and attached label holders. Dimensions: 21<sup>1</sup>/<sub>2</sub>"W x 15<sup>1</sup>/<sub>4</sub>"D x 2<sup>5</sup>/<sub>16</sub>"H.




AMT-2N

**AMT-3N Full Drawer Tray with Three Ampule Holders, \$140.00**  
Thirty-three dividers, three ampule trays, and four securing clips. Dimensions: 21<sup>1</sup>/<sub>2</sub>"W x 15<sup>1</sup>/<sub>4</sub>"D x 2<sup>5</sup>/<sub>16</sub>"H.



AMT-7N

**AMT-6N Modular Drawer Tray, \$105.00**  
Durable, plastic modular tray fits in a 6"D drawer. Includes two vertical rails and six dividers. Dimensions: 21<sup>1</sup>/<sub>2</sub>"W x 15<sup>1</sup>/<sub>4</sub>"D x 5"H.



AMT-3N

**AMT-7N Modular Tray, \$90.00**  
Full drawer tray, eighteen ampule dividers, and seven long dividers. Dimensions: 21<sup>1</sup>/<sub>2</sub>"W x 15<sup>1</sup>/<sub>4</sub>"D x 2<sup>5</sup>/<sub>16</sub>"H.




AMT-8N

**AMT-8N Modular Tray, \$90.00**  
Full drawer tray, twenty-nine ampule dividers, and four long dividers. Dimensions: 21<sup>1</sup>/<sub>2</sub>"W x 15<sup>1</sup>/<sub>4</sub>"D x 2<sup>5</sup>/<sub>16</sub>"H.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

## Mini/Narrow Cart Trays



AMT-4N

**AMT-4N Full Drawer Tray for Mini Cart, \$90.00**  
Full drawer tray with nine ampule dividers, ampule tray, rail, and four adjustable dividers. Fits one per drawer. 14<sup>3</sup>/<sub>16</sub>"W x 16"D x 2<sup>1</sup>/<sub>2</sub>"H



AMT-5N

**AMT-5N Full Drawer Tray for Mini Cart, \$70.00**  
Full drawer tray with two rails and four adjustable dividers. Fits one per drawer. 14<sup>3</sup>/<sub>16</sub>"W x 16"D x 2<sup>1</sup>/<sub>2</sub>"H



## ...Everything In Its Place.

### ADV-1N Modular Drawer Divider Kit, \$50.00

Subdivides a 3"H drawer to better organize your supplies. Includes two 22" rails, five 15 1/2" rails, ten 6 1/2" dividers, and three "L" clips.



**IMMEDIATE**, accurate access

**FASTER** emergency response

**EASIER**, better inventory control



### ADT-1 Crash Cart Drawer Tray with Lid, \$50.00

#### Control Your Drawer Inventory

Clear acrylic lid with thumb hole can be secured with locking seal.

Inside dim.: 21"W x 15"D x 2 1/2"H.

Outside dim.: 22"W x 16"D x 2 3/4"H.



### ADV-10N Modular Divider Kit for ADT-1 Drawer Tray, \$35.00

Fits inside ADT-1 Drawer Tray. Includes two 20 5/8"L grooved plastic rails, four 14"L rails, six 3 15/16"L dividers, and six 5 1/2"L dividers.



### ADV-6N Modular Drawer Divider Kit for 6", 9", and 12" Drawers, \$75.00

5" high modular kit fits 6", 9", or 12"H drawer. Includes two 22" rails, four 15 1/2" rails, and six 6 5/8" dividers.



### APS-1 Pull Tite™ Seal, 100/pkg., \$21.00

#### Easy-to-use Drawer Inventory Control

9 7/8" polypropylene seals fit the ADT-1 Drawer Tray. Consecutively numbered.



### ADB-1 Security Bag, 20/pkg., \$19.00

26 1/2" x 18 1/2". Fits over 3" trays (ADT-1, AMT-1, AMT-2, and AMT-3 Trays). Perforated center shows you when it's been used.

### ADB-2 Security Bag, 20/pkg., \$21.50

32 1/2" x 21 1/2". Fits over 5" tray.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

# Drawer Trays for Stock Exchange



AC-670S



ACC-655S



AC-665

AC-672S

## Drug Trays

**AC-670S Divisible Drug Tray, 3"D, \$16.00**

Sand color. Inside dimensions: 14<sup>7</sup>/<sub>8</sub>"L x 9<sup>1</sup>/<sub>4</sub>"W. Holds up to eleven short dividers and seven long dividers (*not included*). Drawers hold a maximum of two trays, except for the Mini and Narrow Carts.

**AC-675S Short Divider for AC-670S, \$3.00** Max. slots: 11

**AC-680S Long Divider for AC-670S, \$3.40** Max. slots: 7

**AC-672S Divisible Drug Tray, 5"D, \$19.50**

Same as the AC-670S (*see left*), but made for deeper drawers. Drawers hold a maximum of two trays.

**AC-685S Short Divider for AC-672S, \$3.40** Max. slots: 11

**AC-690S Long Divider for AC-672S, \$4.10** Max. slots: 7

**ACC-655S Clear Acrylic Tray Cover, 2/pkg., \$20.00** (*see above*)

See-through cover fits AC-670S and AC-672S Divisible Drug Trays. Thumb notch ensures easy access.

**AC-665 Drug Tray, Small, 2<sup>1</sup>/<sub>2</sub>"D, \$13.50** (*see above*)

Blue color. Inside dimensions: 9<sup>1</sup>/<sub>4</sub>"L x 6<sup>5</sup>/<sub>8</sub>"W. Drawers hold a maximum of four trays.

**AC-675 Long Divider for AC-665, \$3.30** Max. slots: 5

**AC-677 Short Divider for AC-665, \$2.80** Max. slots: 7



ASB-2

## Armstrong Security Box

**ASB-2 Armstrong Security Box, \$155.00**

### Double-Lock Security

This Lockable Security Drug Box can be portable or fastened to the inside of a 3"D drawer. The Security Box is designed

to fit two per drawer, or next to our AC-670S 3"D Drug Tray.

Dimensions: 11"W x 16"D x 2<sup>5</sup>/<sub>8</sub>"H.

## Plastic Divider Sets

**AEF-3 Plastic Divider Set for 3" Drawer, \$55.00**

Two 22"L rails and six 16<sup>1</sup>/<sub>4</sub>"L x 2<sup>5</sup>/<sub>8</sub>"H dividers.



AEF-6

**AEF-3L Long Plastic Divider Set for 3" Drawer, \$55.00**

Four 22"L rails and two 16<sup>1</sup>/<sub>4</sub>"L x 2<sup>5</sup>/<sub>8</sub>"H dividers.

**AEF-6 Plastic Divider Set for 6" Drawer, \$70.00**

Two 22"L rails and six 16<sup>1</sup>/<sub>4</sub>"L x 5<sup>5</sup>/<sub>8</sub>"H dividers.

**AEF-6L Long Plastic Divider Set for 6" Drawer, \$74.00**

Four 22"L rails and two 16<sup>1</sup>/<sub>4</sub>"L x 5<sup>5</sup>/<sub>8</sub>"H dividers.

**AEF-9 Plastic Divider Set for 9" Drawer, \$73.00**

Two 22"L rails and six 16<sup>1</sup>/<sub>4</sub>"L x 8<sup>3</sup>/<sub>4</sub>"H dividers.



AEF-6L

**AEF-9L Long Plastic Divider Set for 9" Drawer, \$76.00**

Four 22"L rails and two 16<sup>1</sup>/<sub>4</sub>"L x 8<sup>3</sup>/<sub>4</sub>"H dividers.


All products are compatible with both Standard and Premier™ Carts unless otherwise stated.




# Cart Accessories



**ACG-1 Caster Guard,**  
4/set, \$210.00  
For 5" dia. wheels.  
Durable with living hinge,  
locking tab, closure  
positioning/support pins.



**ATD-2 Tape Dispenser for  
Trellis Rail,** \$75.00  
includes two trellis brackets  
**ATD-3 Tape Dispenser for  
Cart Top,** \$55.00  
(not pictured)



**APH-1 STANDARD Push  
Handle,** \$60.00  
Does not fit A-SMART  
Premier or AMC-7 Carts.



**ATW-3 Soft-grip Twin  
Push Handles,** \$55.00



**ACC-24 Navy Cart Cover** for  
24" A-SMART® Cart, \$150.00  
**ACC-27 Navy Cart Cover** for  
27" A-SMART Cart, \$150.00  
**ACC-30 Navy Cart Cover** for  
30" A-SMART Cart, \$150.00  
Nylon covers with Velcro tabs.



**NEW DESIGN**  
**AQC-1 Silent Knight  
Tablet Crusher,** \$115.00  
Includes fifty pouches.  
**AQC-2 Silent Knight  
Pouch,** 1000/box, \$55.00



**APS-5 Plastic Seals,**  
100/pkg., \$32.00  
Consecutively numbered.



**ADL-1 Drawer Liner,**  
\$22.00  
Black vinyl and latex free.  
17" x 22". Fits standard  
size drawers.



**ALH-1 Label Holders and  
Cards (blank),** 100/pkg.,  
\$48.00 3" x 1"  
**ALH-2 Label Holders,**  
100/pkg., \$85.00 7" x 2"  
**ALH-12 Label Holders,**  
12/pkg., \$10.00 7" x 2"



**APS-1 Pull-tite Seals,**  
100/pkg., \$21.00  
97/8" polypropylene seals  
open without using scissors.  
Consecutively numbered.



**AUH-1P PREMIER Utility  
Hooks,** 2/pkg., \$22.00  
Mount in corner tracks.



**AUH-1 STANDARD Utility  
Hooks,** 2/pkg., \$21.00  
Mount to sides or back  
of cart.

*All products are compatible with both Standard and Premier™ Carts unless otherwise stated.*

*To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 114.*

# Utility Carts



AUC-400



AUC-412



AUC-200

## AUC-400 Stainless Steel Utility Cart, \$230.00

Cart comes with three shelves, handle with "easy-grip" edging, and smooth-rolling 3 1/2" casters.

- Weight capacity: 300 lbs.
- Overall size: 16"D x 24"W x 32"H
- Heavy-duty, 22 gauge shelves
- Welded construction
- 16 gauge legs

## AUC-412 Stainless Steel Cart with Lock, \$385.00

- Three shelves
- 4" drawer with lock
- Weight capacity: 300 lbs.
- Overall size: 16"D x 21"W x 30"H
- Heavy-duty, 22 gauge shelves
- Welded construction
- 16 gauge legs

## AUC-200 High-Density Poly Cart, \$125.00

- Three sturdy chemical resistant polyethylene plastic shelves
- American made with lifetime warranty
- Shelves have 1/4" lip to prevent objects from falling off
- 4" swivel casters, two with locking brake
- Dimensions: 24"W x 18"D x 34"H
- Weight capacity: 300 lbs.



AST-5



AST-10



AUC-410

## AST-5 Stainless Steel Utility Cart with Single Cabinet, \$298.00

- Sturdy, welded, stainless steel construction
- Fully-enclosed cabinet with hinged door and magnetic catch
- Weight capacity: 300 lbs.
- Four 3 1/2" swivel casters
- Cart dimensions: 16"D x 21"W x 30"H
- One key lock

## AST-10 Stainless Steel Utility Cart with Dual Cabinet, \$397.00

- Sturdy, welded, stainless steel construction
- Fully-enclosed cabinets with hinged door and magnetic catch
- Weight capacity: 300 lbs.
- Four 3 1/2" swivel casters
- Cart dimensions: 16"D x 21"W x 30"H
- Two key locks

## AUC-410 Stainless Steel Utility Cart, \$185.00

- Three 22 gauge shelves
- Weight capacity: 300 lbs.
- Polished stainless steel
- 16 gauge legs
- Welded construction
- Cart dimensions: 16"D x 21"W x 30"H



# Utility Carts



AUC-300

## AUC-300 Poly Cart with Tub-Style Shelves, \$150.00

- Three sturdy tub-style chemical resistant polyethylene plastic shelves
- American made with lifetime warranty
- Tub depth is 2<sup>3</sup>/<sub>4</sub>" with a 12" shelf clearance
- 4" swivel casters, two with locking brake
- Dimensions: 24"W x 18"D x 34"H
- Weight capacity: 300 lbs.



AUC-100

## AUC-100 Poly Cart, \$120.00

- 16" x 20" molded shelves with raised lip on all sides
- Aluminum legs with rounded corners and clear protective finish
- 3" swivel casters
- Weight capacity: 150 lbs.
- 29" working height
- Easy to assemble



AUC-150

## AUC-150 Z Cart, \$349.00

This specialty cart is perfect for use in hospitals or a doctor's office.

- Welded steel construction
- Durable white powder-coat finish
- 16" x 22" shelf and drawer
- Shipped fully assembled
- Top shelf: 22" x 14" x 8"
- 200-lb. weight capacity
- Dimensions: 16"D x 22"W x 40"H



ACW-5

## ACW-5 Mail Cart, \$315.00

This cart comes with two removable baskets, adjustable hanging file folder guide rails for use with legal or letter sized folders, and non-marking casters. Contents are visible from all sides.

### Top basket:

17"W x 25"L x 11"H

### Bottom basket:

17"W x 25"L x 7"H

### Cart:

23"W x 30"L x 39"H



ACW-1

## ACW-1 General Purpose Wire Push Cart, \$400.00

These versatile carts feature three 18" x 36" shelves that can be adjusted or removed to suit your needs. Rubber wheels provide smooth, quiet operation.

- Three adjustable shelves
- Cart size: 18"W x 36"L x 38<sup>1</sup>/<sub>2</sub>"H
- Four 5" rubber casters
- Two push handles
- Zinc-plated finish



AUC-550S

## AUC-550S Wire Linen Cart with Cover, \$705.00

This sturdy, wire linen cart with a solid steel bottom shelf features a Velcro-fastened nylon cart cover. Dimensions: 18"W x 24"L x 60"H.

- Four swivel casters
- Three spacious shelves
- Nylon cover which opens easily and fastens out of the way
- Zinc-plated finish
- Weight capacity: 400 lbs.
- Shipping weight.: 55 lbs.
- Some assembly required

# Wire Shelving

## Mobile Wire Shelving

### Mobile Wire Shelving Units

These heavy-duty, zinc-plated steel units provide convenient, portable storage for linens, equipment, and other items. Units come with four shelves (two adjustable), poles, and four swivel casters (two with brake). Shipped unassembled.

Catalog No.	Shelf Size	Height	Price
ACW-1836	18" x 36"	68"	\$495.00
ACW-1836C	Wire Cart Cover		\$240.00
ACW-1848	18" x 48"	68"	\$550.00
ACW-1848C	Wire Cart Cover		\$250.00

Catalog No.	Shelf Size	Height	Price
ACW-1860	18" x 60"	68"	\$610.00
ACW-1860C	Wire Cart Cover		\$270.00
ACW-2436	24" x 36"	68"	\$560.00
ACW-2436C	Wire Cart Cover		\$255.00
ACW-2448	24" x 48"	68"	\$620.00
ACW-2448C	Wire Cart Cover		\$300.00
ACW-2460	24" x 60"	68"	\$710.00
ACW-2460C	Wire Cart Cover		\$320.00



ACW-1836



ACT-1836

### 3-Sided Mobile Cart

Post-mounted panels eliminate the need for a top shelf. 3" x 3" mesh back and side enclosures keep materials from falling during movement. Weight capacity is 1,200 lbs. Height of cart is 70".

ACT-1836 3-Sided Mobile Cart, 18" x 36" Shelves, \$595.00
ACT-1836C 3-Sided Mobile Cart Cover, \$230.00

ACT-1848 3-Sided Mobile Cart, 18" x 48" Shelves, \$650.00
ACT-1848C 3-Sided Mobile Cart Cover, \$200.00
ACT-2448 3-Sided Mobile Cart, 24" x 48" Shelves, \$795.00
ACT-2448C 3-Sided Mobile Cart Cover, \$200.00

*Wire carts & shelf units come in a wide range of additional sizes. Please call for ordering details.*

## Wire Shelving

### AW-4618 Wire Shelving Unit, \$390.00

Convenient storage for supplies, linens, and equipment. Starter unit includes four posts and four 18" x 48" adjustable shelves. Zinc-plated finish.

Overall size: 18"W x 48"D x 74"H.

### AW-4618A Wire Add-On Shelving Unit, \$370.00

Attaches to an existing AW-4618 Wire Shelving Unit. Includes two posts and four 18" x 48" adjustable shelves. Zinc-plated finish. **NOTE:** Requires an AW-4618 for installation (listed left). Other sizes available.

Overall size: 18"D x 48"W x 74"H.



AW-4618



- A**
- A-SMART Auto-Locking Carts . . . . . 80-81, 83, 103
  - A-SMART Bedside Carts . . . . . 82-83
  - A-SMART Mini Cart . . . . . 82-83, 103-104
  - A-SMART Narrow Cart . . . . . 82-83, 103, 105
  - A-SMART Premier Carts . . . . . 73-109
    - Accessories . . . . . 108, 110-121
    - Anesthesia . . . . . 85-91
    - Electronic Touchpad . . . . . 101
    - Emergency . . . . . 93-99, 106-107
    - Isolation . . . . . 109
    - Locking Accessories . . . . . 120
    - Mobile Workstations . . . . . 97-103
    - Pediatric Emergency . . . . . 118-119
  - A-SMART Standard Carts . . . . . 84-133
    - Accessories . . . . . 122-133
    - Anesthesia . . . . . 96-102
    - Emergency . . . . . 104-110
    - Isolation . . . . . 109
    - Mobile Workstations . . . . . 84-90
    - Push-Button . . . . . 100
  - A-SMART Wireless Auto-Locking Carts . . . . . 77-79
  - Actar 911 CPR Manikins . . . . . 6-7
    - Cadet Pac (child) . . . . . 6
    - Infantry (infant) . . . . . 7
    - Patrol (adult) . . . . . 6
    - Squadron (adult) . . . . . 6
  - Actar AED Trainer . . . . . 5
  - Actar D-fib CPR Manikins . . . . . 4
  - AED Trainers . . . . . 5, 15
  - Airway Management Trainers . . . . . 15, 16-17
  - ALS Manikins . . . . . 13-16
  - Anesthesia Carts . . . . . 84-91
  - Anesthesia Chairs . . . . . 35-37
  - Arm, Arterial Puncture . . . . . 21
  - Arm, I.V. Training . . . . . 22, 30
  - Arrhythmia Training . . . . . 12-13
  - Arterial Puncture Arm . . . . . 21
  - Aspiration Needles . . . . . 23, 67
  - Auscultation Manikin . . . . . 19
  - Auto-Locking Cart . . . . . 77-81, 83, 103
- B**
- Backboards, Emergency . . . . . 63
  - Bags, EMS & Trauma . . . . . 68-69
  - Bedside Carts . . . . . 82-83
  - Blood Pressure Equipment . . . . . 51
  - Broselow ColorCode Cart . . . . . 106-107
  - Broselow Pediatric Emergency Tape . . . . . 41-45, 106
  - Broselow Pediatric Equip. Bag . . . . . 45, 69
  - Broselow® Pediatric Resuscitation System . . . . . 52-56, 106-107
- C**
- Cart Accessories . . . . . 108, 110-121
    - Accessory Panels . . . . . 114
    - Cart Extender . . . . . 110
    - Catheter Holders . . . . . 115
    - Drawer Trays & Dividers . . . . . 111, 118-120
    - Fluid Tray System . . . . . 115
    - Glove Box . . . . . 110
    - I.V. Pole . . . . . 113
    - Locking . . . . . 108
    - Outlet Strip . . . . . 115
    - Plastic Top . . . . . 113
    - Railings . . . . . 113
    - Scope Holders . . . . . 102
    - Sharps Container Bracket . . . . . 115
    - Shelving Units . . . . . 116
    - Suction Unit Holders . . . . . 113
    - Tilt-Bin Organizers . . . . . 117
    - Trellis System . . . . . 112
    - Wire Baskets . . . . . 114
  - Carts** . . . . . 70-114
    - Accessories . . . . . 108, 110-121
    - Anesthesia . . . . . 84-91
    - Auto-Locking . . . . . 77-81, 83, 103
    - Bedside . . . . . 82-83
    - Electronic Touchpad . . . . . 101-102
    - Emergency . . . . . 92-99, 106-107
    - Isolation . . . . . 109
    - Keyless Locking . . . . . 67-71, 101-104
    - Mini . . . . . 82-83, 103-104
    - Mobile Wire Shelving Unit . . . . . 124
    - Narrow . . . . . 82-83, 103, 105
    - Pediatric Emergency . . . . . 106-107
    - Stainless Steel . . . . . 122-123
    - Wide . . . . . 102
    - Wireless Auto-Locking . . . . . 77-79
  - Casualty Simulation Kits . . . . . 24
  - Catheterization Simulators . . . . . 32
  - Catheter Holders, Cart-Mounted . . . . . 115
  - Central Lines Injection Trainer . . . . . 21
  - Central Lines Insertion Trainer . . . . . 21
  - Cervical Collars . . . . . 64-65
  - Chairs, Anesthesia . . . . . 35-37
  - Check Units . . . . . 58
  - Child Tape, Broselow . . . . . 41-45, 106
  - Code Cart, Pediatric . . . . . 106-107
  - Code Carts . . . . . 92-99, 106-107
  - Collars, Cervical . . . . . 64-65
  - ColorCode Ped. Emerg. Cart . . . . . 106-107
  - CPR Adjuncts . . . . . 11
  - CPR Backboard . . . . . 11
  - CPR Manikins . . . . . 4, 6-10
    - Adult . . . . . 4, 6, 8, 10
    - Actar . . . . . 4, 6-7
    - Child . . . . . 7-10
    - Group Training . . . . . 4, 6-7
    - Infant . . . . . 7, 9-10
    - Laerdal . . . . . 9-10, 13
  - CPR Masks . . . . . 33-34
  - Crash Cart, Pediatric . . . . . 106-107
  - Crash Carts . . . . . 92-99, 106-107
  - Crisis Manikins . . . . . 12-13
- D**
- DCell Suction . . . . . 40
  - Defib Manikin . . . . . 4, 12-13
  - Defib Pad Adapters . . . . . 12
  - Diagnostic Tools . . . . . 41-45, 47-50
  - Disposable Laryngoscope Blades . . . . . 48
  - Dopplers . . . . . 53
  - Drawer Cover . . . . . 114
  - Drawer Trays . . . . . 111, 118-120
  - Drug Box . . . . . 38, 112
  - Duet Suction Unit . . . . . 39
- E**
- Electronic Touchpad Carts . . . . . 101-102
  - Emergency Airway Systems . . . . . 46-48
  - Emergency Bags . . . . . 68-69
  - Emergency Cart Packages . . . . . 92-97, 107
  - Emergency Carts . . . . . 92-99, 106-107
  - Emergency Supplies . . . . . 67
  - Emergency System, Pediatric . . . . . 41-45, 106-107
  - EMS Bags . . . . . 68-69
  - Extrication Collars . . . . . 64-65
  - Eye Trauma System . . . . . 66
- F**
- Face Shields, CPR . . . . . 33-34
  - Face Shields, Manikin . . . . . 11
  - Fluid Tray System for Carts . . . . . 115

# Index

## G

Gas Connectors. . . . . 58  
Glove Dispensers. . . . . 38, 110

## H

Head Immobilizer. . . . . 62-63  
Heine Laryngoscopes. . . . . 46-47

## I

I.V. Trainers. . . . . 22, 30  
Immobilization. . . . . 59-65  
Injectable Training Arm. . . . . 20, 22  
Injection/Infusion Simulator, Pediatric. . . . . 23  
Intradermal Injection Simulator. . . . . 20  
Intramuscular Injection Trainer. . . . . 24  
Intraosseous Infusion Simulators. . . . . 23  
Intraosseous Needles. . . . . 23, 67  
Intubation Trainers. . . . . 15, 17-20  
Isolation Carts. . . . . 109  
I.V. Pole, Cart-Mounted. . . . . 113  
I.V. Stand, Mobile. . . . . 57  
I.V. Training Simulators. . . . . 22, 30

## L

Laerdal AED Trainer. . . . . 15  
Laerdal ALS Manikins. . . . . 13  
Laerdal CPR Manikins. . . . . 9-10  
Laerdal Intubation Trainers. . . . . 17, 19  
Laryngoscopes. . . . . 46-48  
Lockable Drawer Cover. . . . . 114  
Locking Accessories, Cart. . . . . 108

## M

Manikin  
    AED Training. . . . . 5, 12-16  
    ALS. . . . . 12-16  
    Crisis. . . . . 12-13  
    CPR. . . . . 4, 6-10  
    Defib. . . . . 4, 12-16  
    Geriatric. . . . . 27-29  
    Hospital Training. . . . . 27-29  
    Intubation. . . . . 15, 17-20  
    Nursing. . . . . 27-29  
    Patient Care. . . . . 27-29  
    Rescue. . . . . 25  
    Trauma. . . . . 25  
    X-ray Positioning. . . . . 26  
Manikin Storage Cases. . . . . 11  
Microkey CPR Barrier Device. . . . . 33  
Micromask CPR Mask. . . . . 33  
Microshield CPR Barrier Device. . . . . 33

Mini Carts. . . . . 82-83, 103-104  
Mobile Shelving. . . . . 124  
Mobile Workstations. . . . . 84-91  
Modular Drawer Trays. . . . . 111, 118-120  
Morgan Lens. . . . . 66

## N

Narcotics Drawer Cover. . . . . 114  
Narrow Cart. . . . . 82-83, 103-105  
Newborn Simulator. . . . . 31  
Nursing Manikins. . . . . 27-29

## O

O.B. Manikins. . . . . 30-31  
Ophthalmoscope. . . . . 49-50  
Otoscope. . . . . 49-50  
Oxygen Bags. . . . . 68-69

## P

Pacing/Defib Simulator. . . . . 12  
Patient Care Manikins. . . . . 27-29  
Patient Transfer Device. . . . . 57  
Peak Flow Meters. . . . . 53  
Pediatric Emergency Carts. . . . . 106-107  
Pediatric Emergency System  
    . . . . . 41-45, 106-107  
Pediatric Measuring Tape. . . . . 41-45  
Pelvic Exam Model. . . . . 32  
Personal Protection Equipment. . . . . 33-34  
Phantom, X-ray. . . . . 26  
Pill Crusher. . . . . 121  
Pneumothorax Simulator. . . . . 22  
Premier Carts (see A-SMART). . . . . 73-121

    Accessories. . . . . 108, 110-121  
    Anesthesia Packages. . . . . 85, 87  
    Emergency Packages. . . . . 93, 95, 107  
Premier ColorCode Cart. . . . . 106-107  
Pressure Infusers. . . . . 58  
Privacy Screens. . . . . 55  
Proximity Reader. . . . . 81  
Pulse Oximeters, Hand-Held. . . . . 52, 67  
Punch Card Tray. . . . . 113

## R

Radiography Manikin. . . . . 26  
Rescue Manikins. . . . . 25  
Resusci Patient Face Shields. . . . . 34  
Resuscitation Masks. . . . . 33-34  
Resuscitators. . . . . 58

## S

S-SCORT Suction Units. . . . . 39-40  
Scales. . . . . 56

Schrader Connectors, Gas Outlets. . . . . 58  
Screens, Privacy. . . . . 55  
Security Bag. . . . . 119  
Security Box, Cart. . . . . 112, 114  
Shelving, Mobile. . . . . 124  
Shelving Units for Carts. . . . . 116  
Speedsplint Disposable Splints. . . . . 59, 62  
Sphygmomanometers. . . . . 50-51  
Splints. . . . . 59-62  
Stainless Steel Carts. . . . . 122-123  
Standard Carts (see A-SMART). . . . . 72-121  
    Accessories. . . . . 110-121  
    Anesthesia Packages. . . . . 84, 86  
    Emergency Packages. . . . . 92, 94  
Stethoscopes. . . . . 54  
Stool, Anesthesia. . . . . 36  
Suction Units. . . . . 39-40

## T

Tape, Broselow Ped. Emergency. . . . . 41-45  
Tilt Bins. . . . . 117  
Trauma Cases. . . . . 68-69  
Trauma Manikin. . . . . 25  
Trellis Systems, Cart. . . . . 112

## U

Universal Precaution Kit. . . . . 67  
Utility Carts. . . . . 122-123

## W

Wide Carts. . . . . 102  
Wire Baskets. . . . . 114  
Wire Shelving. . . . . 124  
Wireless Auto-Locking Carts. . . . . 77-79  
Wound Care Simulator. . . . . 24

## X

X-ray Positioning Manikin. . . . . 26



# Armstrong Medical Order Form

575 Knightsbridge Parkway • P.O. Box 700 • Lincolnshire, IL 60069-0700 • 847/913-0101 • Toll Free: 800/323-4220 • FAX: 847/913-0138

**A service charge of \$20.00 will be added to orders under \$100.00.**

**Payment terms are net 30 DAYS upon receipt of valid purchase order from Government facilities, hospitals, and other institutions with credit approval.**

**Visa, MasterCard, American Express and Discover orders gladly accepted.**

<b>BILL TO:</b> <hr/> Name _____ Title _____ Organization _____ Street _____ City _____ State _____ ZIP Code _____ Phone _____ Fax _____ Email _____ Is the above address new? <input type="checkbox"/> Yes <input type="checkbox"/> No _____ I have ordered from Armstrong before: <input type="checkbox"/> Yes <input type="checkbox"/> No	<b>SHIP TO:</b> <hr/> Name _____ Title _____ Organization _____ Street _____ City _____ State _____ ZIP Code _____ Phone _____ Fax _____ Email _____ Special Delivery Instructions _____ _____
---	---

**Purchase Order #** \_\_\_\_\_ **Customer Code #** \_\_\_\_\_

Quantity	Catalog Number	Description	Price Each	Total Price

When paying in advance, please call 800/323-4220 and ask our Customer Service Staff for the correct tax and shipping charges.

*Orders paid by credit card may be subject to a 3% processing fee.*

Visa  MasterCard  American Express  Discover

Credit Card # \_\_\_\_\_

CID# \_\_\_\_\_ Exp. Date \_\_\_\_\_

Name as it appears on card \_\_\_\_\_

Signature \_\_\_\_\_

Print Name \_\_\_\_\_

Title \_\_\_\_\_ Date \_\_\_\_\_

MERCHANDISE TOTAL	
AMOUNT FOR TAX*	
SHIPPING CHARGES	
3% CREDIT CARD PROCESSING FEE	
<b>TOTAL</b>	

**\*IF TAX EXEMPT, PLEASE ATTACH A COPY OF YOUR CERTIFICATE.**

# Ordering and Shipping Information

## CORPORATE OFFICE

Armstrong Medical Industries, Inc. • 575 Knightsbridge Parkway • P.O. Box 700 • Lincolnshire, IL 60069-0700  
847/913-0101 • 800/323-4220 • FAX: 847/913-0138  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

Your opinion counts, let us hear from you at: [csr@armstrongmedical.com](mailto:csr@armstrongmedical.com)

Hours: 8:00 - 4:30 CST\*

## PLACING YOUR ORDER

- When placing your order with us, please refer to your customer code number, which appears on all of our invoices and packing lists.
- Please write your customer code number on the outside of your catalog for immediate reference when placing phone orders.
- All federal government orders, other than credit card, must be in writing. Please fax your order to us or send in by mail. Thank you.
- Although they may not be listed individually, Customer Service can help you order special parts; please just ask.

## TERMS OF SALE

- Net 30 days for regular customers with established credit.
- Visa, MasterCard, American Express, and Discover accepted.
- Orders paid by credit card may be subject to a 3% processing fee
- A service charge of \$20.00 will be added to orders under \$100.00.
- Service or handling charges are not refundable.
- Terms: payment in advance or net 30 days based upon credit approval.
- C.O.D. shipments are not available.
- Delinquent accounts assessed a finance charge equal to an annual rate of 18%.

## CATALOG PRICING

- Prices do not include assembly
- Products are priced "each" unless otherwise noted.
- Our prices are as current as possible, but many of our suppliers issue price changes without notice. Therefore, all prices are subject to change without prior notice.
- When calling our Customer Service Department for any reason, please ask if your price listings are current. Since all catalogs will be numbered or dated, please refer to the date.
- This catalog is available free of charge to any qualified individual, organization, or institution.

## SHIPPING TERMS

- Our standard shipping terms are prepaid and added to the invoice.
- In stock, small parcels are delivered within 7-10 days.
- Cart delivery time confirmed at the time of order.
- Standard shipping for small packages is FedEx® Ground. Expedited shipping is available for most items at an additional cost.

## NATIONAL SALES FORCE

- Armstrong Medical Industries, Inc. supports our customers with a National Sales Force to assist with product demonstrations or any technical questions about our products. Please ask our Customer Service Department for the name of the representative in your area.

## RETURNED GOODS

- All returned goods must have an authorization number assigned by our Customer Service Department. Returns must be within 30 days from date of original shipment.
- Please refer to your invoice number and customer code number when phoning in your request for returning merchandise.
- All returned goods must be sent prepaid. Any return that is not prepaid will be refused.
- Diagnostic, resuscitation equipment, sterilized items, non-stock items, customized orders, A-SMART® Carts, and supply carts are non-returnable. Our Customer Service Department will be glad to answer any questions you may have regarding returns.

## RESTOCKING CHARGE

- A 15% restocking charge may be charged to your account for returned merchandise.
- An exception to this policy would be if we shipped incorrect merchandise. If this is the case, please call our Customer Service Department and report the incorrect shipment within 10 days. The same policy regarding an authorization number applies.

## DAMAGED GOODS

- If you receive a small package that you believe is damaged, we ask that you call your delivery service for an inspection of the damaged articles.
- Please save all boxes and packing material to prove that the items were packed properly.
- **All small parcel damage claims must be reported within 10 days after receipt of merchandise. All truck order damage claims must be reported within 7 days. All moving van damage claims must be noted at the time of delivery.**
- Photos required for all damage claims of both packaging and product damage.
- If you receive a shipment via truck, inspect all cartons and note visible damage on Bill of Lading **BEFORE SIGNING YOUR NAME.**

*\*When ordering, please remember that the lunch hours are always our busiest times. Thank you for your patience.*



# Products



*You Can*



*Depend  
On*



Prices and products are subject to change without notice.  
©Copyright 2017 Armstrong Medical Industries, Inc. All rights reserved.