

# *Light* Aluminum Premier™ Carts



***Armstrong Medical***

A - S M A R T® C A R T S

More Standard Features, More Flexibility,

# MORE CART FOR YOUR MONEY

## The Benefits Of Buying An Armstrong **ALUMINUM** A-SMART® Premier™ Cart ... **Think LIGHT.**

### The Benefits Of Buying **ALUMINUM** Instead of Steel

#### Steel

76cm Steel Cart weighs: **74.8kg**

69cm Steel Cart weighs: **72.6kg**

61cm Steel Cart weighs: **65.8kg**

Liquids can cause rust

Less maneuverability

Risk of accidental injury increases  
with additional weight

#### Aluminum

76cm Aluminum Cart weighs: **41.7kg**

69cm Aluminum Cart weighs: **38.5kg**

61cm Aluminum Cart weighs: **35.4kg**

Reduces risk of rusting

Easier maneuverability

Easier to push, turn, and **STOP**; avoid injuries  
to patients, visitors, nurses, and doctors

### FEATURES OF A-SMART® PREMIER™ **ALUMINUM** CART

- **ALUMINUM** Carts are 45% lighter than Steel
- **Full Extension Ball Bearing Drawer Slides**  
*Ensures that your drawers will open and close smoothly every time.*
- **Quality Aluminum Construction**  
*One-piece outer shell construction providing you with added strength and stability. Locking components consist of zinc plated and galvanized parts.*
- **Double Side-Wall Construction**  
*Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.*
- **Two Interchangeable Locking Systems**
- **Customize Your Locking Configuration**  
*Lock individual drawers, double lock drawers, use multiple types of locks on the same cart.*
- **Slide-Thru Side Shelf**
- **Contour-Fit Plastic Top**
- **Four Vertical Aluminum Corner Mounting Tracks**  
*For easy accessory attachment (eliminates the need for holes).*
- **Stabilizing Frame With Bumper**  
*Casters are moved out from underneath to the outside of the cart. This extra wide base provides your cart with added stability.*
- **18 Colors Available**  
*Allows you flexibility in your carts appearance. Choose colors to specify department or function, or even use multiple colors on the same cart to differentiate drawers.*
- **Highest Quality Swivel Casters (two locking, one tracking)**  
*Quiet, durable, German-made casters. Smoother, easier maneuvering.*
- **Durable Powder Paint**  
*Fully powder coated paint throughout. Prevents chips, scratches, fading, and wearing.*



Look for the Armstrong A-SMART® label on every cart. ***It stands for quality.***

Manufactured to ISO 9001:2015 certified standards.



# **Armstrong A-SMART®**

## **Electronic Touchpad Carts**

Security At The Push Of A Button



# Get Exactly What You Need

# CREATE YOUR OWN CART



With hundreds of accessories, eighteen available colors, flexible locking options (*Electronic Touchpad, Key Locking, Breakaway Locking*), and multiple cart and drawer sizes, you can customize an A-SMART® Cart to fit any purpose or department. With so many options, you can build your own cart so that you **Get Exactly What You Need**.

#### **A-SMART.**

Ostomy Care Cart  
Toxic Shock Cart  
Dressing Cart  
Conscious Sedation Cart  
Bronchoscopy Cart  
Eye Cart  
Block Cart  
Supply Cart  
Respiratory Cart

Trauma Cart  
Computer Cart  
Delivery Cart  
Latex-Free Cart  
Ear, Nose &  
Throat Cart  
Sexual Assault Exam Cart  
ER Supply Cart  
Crash Cart

Isolation Cart  
Endoscopy Cart  
Bio-Terror Cart  
Difficult Airway Cart  
Broselow Cart  
Burn Cart  
Bedside Cart  
ER Cart  
Anesthesia Cart

Phlebotomy Cart  
PICC Line Cart  
Epidural Cart  
Procedure Cart  
**And Many More**

*See The Benefits of Buying An Armstrong Cart*

**[www.armstrongmedical.com](http://www.armstrongmedical.com)**



# Table of Contents



## **Armstrong Medical Industries, Inc.**

Corporate Headquarters  
575 Knightsbridge Parkway • P.O. Box 700  
Lincolnshire, IL 60069-0700  
847/913-0101 • FAX: 847/913-0138  
[www.armstrongmedical.com](http://www.armstrongmedical.com)

<b>Introduction</b>	
Standard Steel Cart Features.....	4
Premier™ All Aluminum Cart Features.....	5
Standard Steel Cart Dimensions.....	6
Premier™ All Aluminum Cart Dimensions.....	7
<b>PREMIUM Emergency Cart Packages</b>	
Standard Steel Emergency Cart Package.....	8
Premier™ All Aluminum Emergency Cart Package.....	9
<b>Emergency Accessory Packages</b>	
Standard Emergency Accessory Packages.....	10
Premier™ Emergency Accessory Packages.....	11
<b>Solid Color Emergency Carts</b>	
Standard Steel Solid Color Emergency Carts.....	12
Premier™ All Aluminum Solid Color Emergency Carts.....	13
<b>Two-Tone Emergency Carts</b>	
Standard Steel Two-Tone Emergency Carts.....	14
Premier™ All Aluminum Two-Tone Emergency Carts.....	15
<b>Broselow® Premier™ Aluminum Pediatric ColorCode Cart.....</b>	<b>16</b>
<b>Electronic Touchpad Bedside Carts.....</b>	<b>17</b>
<b>SUPREME Anesthesia Cart Packages</b>	
Standard Steel Anesthesia Cart Package.....	18
Premier™ All Aluminum Anesthesia Cart Package.....	19
<b>Anesthesia Accessory Packages</b>	
Standard Anesthesia Accessory Packages.....	20
Premier™ Anesthesia Accessory Packages.....	21
<b>Solid Color Mobile Workstations</b>	
Standard Steel Solid Color Key Locking Carts.....	22
Premier™ All Aluminum Solid Color Key Locking Carts.....	23
<b>Two-Tone Mobile Workstations</b>	
Standard Steel Two-Tone Key Locking Carts.....	24
Premier™ All Aluminum Two-Tone Key Locking Carts.....	25
<b>Specialty Carts and Accessories</b>	
Premier™ Locking Bars.....	26
Trellis Systems.....	26
Electronic Touchpad Locking Carts.....	27
Mini and Narrow Carts.....	28-29
ScopeSafe™ Holders & Wide Cart.....	30
Isolation Carts.....	31
<b>Cart Accessories.....</b>	<b>32-Inside Back Cover</b>



**Look for the Armstrong A-SMART® label  
on every cart. It stands for quality.**

## A-SMART® Setting The Standard

**Breakaway Locking Cart**  
Breakaway Security ABL Series

Soft-grip Handles

Double Side-wall Construction

Stabilizing Frame with Bumper prevents injuries. Casters are moved from under the cart to outside its perimeter. This extra-wide base provides your cart with added stability.

Ball Bearing Drawer Slides

(1) Swivel  
(2) Swivel w/Brake  
(1) Tracking Caster

**Electronic Touchpad Locking Cart**

**NEW** Key override comes standard on Electronic Touchpad Security ATE Series


**Key Locking Cart**

APT-2 Optional Plastic Top, see page 35.



Look for the Armstrong A-SMART® label on every cart.  
**It stands for quality.**  
**MADE IN THE USA**

### A-SMART® Standard Features:

- Ball bearing drawer slides
- High-quality swivel casters (two locking, one tracking)
-  Casters
- Double Side-wall construction (see right)
- Soft-grip handles
- Stabilizing frame with bumper
- 18 Colors
- Durable powder paint
- All carts are latex free

### Double Side-Wall Construction

Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.



# Take The Strain Off Your Back. **Think LIGHT.**



## See the Advantages of A-SMART® PREMIER™ **Aluminum** Carts

- Lightweight, aluminum construction
- Two interchangeable locking systems
- Four vertical **Aluminum** corner mounting tracks for easy accessory attachment
- Slide-thru side shelf
- Contour-fit plastic top
- High-quality swivel casters (two locking, one tracking)
- Stabilizing frame with bumper
- Double side-wall construction
- Soft-grip handles
- 18 colors
- Durable powder paint
- Ball bearing drawer slides
- All carts are latex free



## Double Side-Wall Construction

Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.

# Standard Steel Cart Features

## How to Customize An A-SMART® Standard Steel Cart:

### 1. CHOOSE FRAME SIZE

*What kind of locking system do you need? How much drawer space do you want?*

#### ABL Series: Breakaway Locks




#### 6-Dr. ABL Cart Dim.:

**Overall:** 117.5cm High x 63.5cm Deep x 85.1cm Wide  
(includes bumper/handles/casters).  
**Drawer size:** 43.3cm Deep x 56.5cm Wide  
**6-Drawer:** 7.6cm (4), 15.6cm (1), 31.8cm (1).  
**Weight:** 74.8kg

#### 4-Dr. and 5-Dr. ABL Cart Dim.

**Overall:** 101cm High x 65.3cm Deep x 85.1cm Wide  
(includes bumper/handles/casters).  
**Drawer size:** 43.3cm Deep x 56.5cm Wide  
**4-Drawer:** 7.3cm (2), 15.6cm (1), 32.2cm (1).  
**5-Drawer:** 7.3cm (3), 15.6cm (1), 23.8cm (1).  
**Weight:** 65.8kg

#### AKL Series: Key Locks



#### 4-Dr. and 6-Dr. AKL Cart Dim.:


**Overall:** 117.5cm High x 63.5cm Deep x 85.1cm Wide  
(includes bumper/handles/casters).  
**Drawer size:** 43.3cm Deep x 56.5cm Wide  
**4-Drawer:** 7.3cm (1), 23.8cm (3).  
**6-Drawer:** 7.3cm (3), 15.6cm (2), 23.8cm (1).  
**Weight:** 74.8kg

#### 5-Dr. AKL Cart Dim.:

**Overall:** 101cm High x 65.3cm Deep x 85.1cm Wide  
(includes bumper/handles/casters).  
**Drawer size:** 43.3cm Deep x 56.5cm Wide  
**5-Drawer:** 7.3cm (3), 15.6cm (1), 23.8cm (1).  
**Weight:** 65.8kg

#### ATE Series:

Electronic Touchpad Locks w/Key Override



#### 6-Dr. ATE Cart Dim.:

**Overall:** 116.4cm High x 64.1cm Deep x 86.6cm Wide  
(includes bumper/handles/casters).  
**Drawer size:** 43.3cm Deep x 56.5cm Wide  
**6-Drawer:** 7.3cm (3), 15.6cm (2), 23.8cm (1).  
**Weight:** 77.1kg

#### 5-Dr. ATE Cart Dim.:

**Overall:** 99.8cm High x 64.1cm Deep x 86.6cm Wide  
(includes bumper/handles/casters).  
**Drawer size:** 43.3cm Deep x 56.5cm Wide  
**5-Drawer:** 7.3cm (3), 15.6cm (1), 23.8cm (1).  
**Weight:** 68kg

### 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain

### 2. Pick Your Options

#### Choose your drawer fronts...

Choose from any one of our eighteen designer colors on the right. If you want the ABL-4 as shown with Teal drawer fronts, your cart would be ABL-T-4.


#### Four drawer sizes available:

**7.3cm 15.6cm 23.8cm 31.8cm**

If standard drawer arrangements don't suit your needs, slides can be moved 7.3cm up or down (additional charge per cart).



## 18 Designer Colors



## How To Customize An A-SMART® All Aluminum Premier™ Cart:

### 1. Choose Frame Size



#### 61cm (5-Drw) CART DIMENSIONS:

**Overall:** 99.1cm High x 63.5cm Deep x 83.2cm Wide (includes bumper, handles, and casters).

**Drawer size:** 43.3cm Deep x 56.5cm Wide.

**5-Drawer:** 7.3cm (3), 15.6cm (1), 23.8cm (1).

**Weight:** 35.4kg



#### 68.6cm (6-Drw) CART DIMENSIONS:

**Overall:** 107.3cm High x 63.5cm Deep x 83.2cm Wide (includes bumper, handles, and casters).

**Drawer size:** 43.3cm Deep x 56.5cm Wide.

**6-Drawer:** 7.3cm (4), 15.6cm (1), 23.8cm (1).

**Weight:** 38.6kg



#### 76.2cm (6-Drw) CART DIMENSIONS:

**Overall:** 114.9cm High x 63.5cm Deep x 83.2cm Wide (includes bumper, handles, and casters).

**Drawer size:** 43.3cm Deep x 56.5cm Wide.

**6-Drawer:** 7.3cm (3), 15.6cm (2), 23.8cm (1).

**Weight:** 41.7kg

### 2. Choose a Lock or Have Both

#### PKL Series: Key Locks



#### PBL Series: Breakaway Locks



Locks fit either or both corner glides



### 3. Pick Your Options

#### Choose your drawer fronts...

Choose from any one of the eighteen designer colors on the left. To get the PKL-24 as shown with **Hunter Green** drawer fronts, order the PKL-HG-24. To get the PBL-27 as shown with **Red** drawer fronts, order the PBL-AR-27.

#### Four drawer sizes available:

**7.3cm 15.6cm 23.8cm 31.8cm**

If standard drawer arrangements don't suit your needs, slides can be moved 7.3cm up or down (additional charge per cart).



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

**MADE IN THE USA**

**ACCESSORIES,  
Pages 32-Inside Back Cover**

# Standard Premium Emergency Cart Package

BREAKAWAY LOCKING

## Standard Premium Emergency Package


AR-6 with APE-1X Emergency Package





## Premier™ Premium Emergency Package PBL-AR-30 with PAPE-1X Emergency Package

### Double Side-Wall Construction



# Standard Emergency Cart Accessory Packages

## BREAKAWAY LOCKING

### Deluxe Emergency Accessory Package

EB-1 Deluxe Accessory Package (*Does not include cart*)

Complete Package as shown, EB-1 with AR-6 Six-Drawer 76cm Emergency Cart

- ASN-3 Collapsible Side Shelf
- PIV-1 I.V. Pole with Brackets
- APT-2 Plastic Top
- PSR-1 Heavy-duty Railing
- AOB-1 Oxygen Tank Bracket
- ACB-1 Cardiac Board with Brackets
- APS-5 Plastic Breakaway Seals, 100/pkg.
- AUH-1 Utility Hooks, 2/pkg
- Two AC-672S 12.7cm High Trays with four AC-685S (short) Dividers & two AC-690S (long) Dividers (fills one drawer)
- Four AC-670S 7.3cm High Trays with eight AC-675S (short) Dividers and four AC-680S (long) Dividers (fills two drawers)



Modular Drawer Dividers  
*see page 40*



Look for the  
Armstrong A-SMART® label  
on every cart.  
*It stands for quality.*  
**MADE IN THE USA**




### Economy Emergency Accessory Package

BE-1 Economy Accessory Package (*Does not include cart*)

Complete Package as shown, BE-1 with ABL-AR-6 Six-Drawer 76cm Emergency Cart

- PIV-1 I.V. Pole with Bracket
- PSR-1 Heavy-duty Railing
- ACB-1 Cardiac Board with Brackets
- APT-2 Plastic Top
- APS-5 Plastic Breakaway Seals, 100/pkg.
- AOB-1 Oxygen Tank Bracket



Tilt Bins, *see page 37*

NOTE: Carts require some assembly.

*See page 4 for Standard Cart features  
and page 6 for dimensions.*




## Deluxe Emergency Accessory Package

PEB-1 Deluxe Accessory Package (*Does not include cart*)

Complete Package as shown, PEB-1 with PAR-30 Six-Drawer 76cm All Aluminum Emergency Cart

- PIV-1 I.V. Pole with Brackets
- PSR-1 Tubular Steel Railing
- AOB-1 Oxygen Tank Bracket
- PCB-1 Cardiac Board with Brackets
- APS-5 Breakaway Seals, 100/pkg.
- Two AC-672S 12.7cm High Trays with four AC-685S (short) Dividers and two AC-690S (long) Dividers
- Four AC-670S 7.3cm High Trays with eight AC-675S (short) Dividers and four AC-680S (long) Dividers



ADV-1 Modular Drawer Divider Kit, see page 39



PAR-30

Mount either of these accessory packages on any A-SMART® Premier™ All Aluminum Emergency Cart



Look for the Armstrong A-SMART® label on every cart.  
*It stands for quality.*  
**MADE IN THE USA**



PBL-AB-30

## Economy Emergency Accessory Package

PBE-1 Economy Accessory Package (*Does not include cart*)

Complete Package as shown, PBE-1 with PBL-AB-30 Six-Drawer 76cm All Aluminum Emergency Cart

- PIV-1 I.V. Pole with Brackets
- PCB-1 Cardiac Board with Brackets (*shown above*)
- AOB-1 Oxygen Tank Bracket
- APS-5 Breakaway Seals, 100/pkg.

See page 5 for Premier™ Cart features and page 7 for dimensions

# Standard Steel Solid Color Emergency Carts

BREAKAWAY LOCKING

## Manufactured To ISO 9001:2015 Certified Standards. **A-SMART®.**



ABL-B-6

### Six-Drawer 76cm Standard Emergency Cart

**ABL-B-6 Emergency Cart, Beige**

Drawer sizes: 7.3cm (4), 15.6cm, 31.8cm

*Double-wall construction and ball bearing  
drawer slides add to the quality of all  
Armstrong Carts.*



ABL-B-5

### Five-Drawer 61cm Standard Emergency Cart

**ABL-B-5 Emergency Cart, Beige**

Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm

*Built-in, stabilizing frame/bumper gives  
Armstrong Carts the stability for  
emergency situations.*



ABL-B-4

### Four-Drawer 61cm Standard Emergency Cart

**ABL-B-4 Emergency Cart, Beige  
AR-4 Emergency Cart, Red**

*(not pictured)*

Drawer sizes: 7.3cm (2), 15.6cm, 31.8cm



AR-6

### Six-Drawer 76cm Standard Emergency Cart

**AR-6 Emergency Cart, Red**

Drawer sizes: 7.3cm (4), 15.6cm, 31.8cm

*Double-wall construction and ball bearing  
drawer slides add to the quality of all  
Armstrong Carts.*



AR-4

### Four-Drawer 61cm Standard Emergency Cart

**AR-4 Emergency Cart, Red**

Drawer sizes: 7.3cm (2), 15.6cm, 31.8cm



AR-27

### Six-Drawer 69cm Standard Emergency Cart

**AR-27 Emergency Cart, Red**

Overall Dim.: 109.2cm High x 63.5cm  
Deep x 85.1cm Wide

*(includes bumper, handles, and casters)*

Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm  
Weight: 72.6kg

*See page 4 for Standard Cart features  
and page 6 for dimensions.*

*NOTE: Carts require some assembly.*

# Drawers Always Slide Smoothly With Our Double Side-Wall Construction. . . **Think** *LIGHT*.



PBL-B-27

## Six-Drawer 69cm Premier Emergency Cart

PBL-B-27 Emergency Cart, Beige  
Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm



PBL-B-24

## Five-Drawer 61cm Premier Emergency Cart

PBL-B-24 Emergency Cart, Beige  
Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm



PAR-30

## Six-Drawer 76cm Premier Emergency Carts

PAR-30 Emergency Cart, Red  
PBL-B-30 Emergency Cart, Beige  
(not pictured)  
Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm



AC-670S and AC-672S,  
see page 40



PSR-1 Tubular Railing,  
see page 35




Look for the  
Armstrong A-SMART® label  
on every cart.

*It stands for quality.*

**MADE IN THE USA**

### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks



# Standard Steel Two-Tone Emergency Carts

BREAKAWAY LOCKING

**Armstrong A-SMART® Two-Tone Carts**  
come with your choice of 18 colored drawer fronts.



ABL-L-5

**Five-Drawer 61cm  
Standard Emergency Cart**

ABL-L-5 Emergency Cart,  
Lavender

Drawer sizes: 7.3cm (3), 15.6cm,  
23.8cm



ABL-T-4

**Four-Drawer 61cm  
Standard Emergency Cart**

ABL-T-4 Emergency Cart,  
Teal

Drawer sizes: 7.3cm (2), 15.6cm,  
31.8cm



ABL-TC-6

**Six-Drawer 76cm Standard  
Emergency Cart**

ABL-TC-6 Emergency Cart,  
Terra Cotta

Drawer sizes: 7.3cm (4), 15.6cm,  
31.8cm



Look for the Armstrong  
A-SMART® label on every  
cart. *It stands for quality.*

**MADE IN THE USA**

**Three-Drawer 61cm  
Standard Emergency Cart**

ABL-B-3 Emergency Cart  
Beige

Drawers: 15.6cm, 23.8cm (2)



ABL-B-3

## 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain



**AHP-3 Panel for Cart-mounted  
DUET Suction Unit**  
(Retention Bracket not included),  
see page 34




**PSR-1 Tubular Railing,**  
see page 35



**ASN-3 Standard Side Shelf,**  
see page 35

18 Designer  
Colors



# Carts Built To Last Longer!

**Think** LIGHT.



PBL-AR-30

## Six-Drawer 76cm Premier Emergency Cart

**PBL-AR-30 Emergency Cart, Red**  
Drawer sizes: 7.3cm (3), 15.6cm (2),  
23.8cm



PBL-SG-27

## Six-Drawer 69cm Premier Emergency Cart

**PBL-SG-27 Emergency Cart,  
Sage Green**  
Drawer sizes: 7.3cm (4), 15.6cm,  
23.8cm



PBL-SB-24

## Five-Drawer 61cm Premier Emergency Cart

**PBL-SB-24 Emergency Cart,  
Slate Blue**  
Drawer sizes: 7.3cm (3), 15.6cm,  
23.8cm



**AMT-2N Modular Tray,**  
see page 38



**PDS-1 One shelf Unit, Deep,**  
see page 36



Look for the  
Armstrong A-SMART® label  
on every cart.  
*It stands for quality.*  
**MADE IN THE USA**

## Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster **TENTE Casters**
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

# Premier™ Aluminum ColorCode™ Cart

BREAKAWAY LOCKING

Not MRI Compatible

## Use The ColorCode™ Cart In Pediatric Resuscitation Emergencies.




PBL-PC-9

PBL-PC-9 **BASIC** A-SMART®  
**Premier™** Broselow® ColorCode™  
Cart (pictured above)

Includes Nine-Drawer A-SMART Premier Breakaway Lock Cart (76cm total drawer space), slide-thru side shelf, 76cm locking bar, stabilizing frame/bumper system, seamless plastic top, soft-grip push handles, three accessory attachment panels, and all-swivel casters (two locking and one tracking).



Use OPTIONAL 7.3cm drawer locks so that only the accessed drawers need restocking.



**TENTE** Casters

PBL-PC-9 shown above with  
**OPTIONAL** Accessory Package  
PC-9E: Includes AOB-1 Oxygen Tank  
Bracket, ASU-1 Suction Unit Holder, PIV-1

IV Pole, PCB-1 Cardiac Board (not shown),  
PBR-3 (8) 7.3cm Breakaway Locking Bar,  
PBR-6 (1) 15.6cm Breakaway Locking Bar,  
and APS-5 Plastic Seals.


### Greater Continuity, Accuracy, and Efficiency

PBR-3 7.3cm Breakaway  
Locking Bar  
PBR-6 15.6cm Breakaway  
Locking Bar

Alternative locking options help  
simplify inventory or limit access to  
sensitive materials. See page 26 for  
more locking options.



## 18 Designer Colors



# NURSE SERVER/BEDSIDE Carts



AMC-3E-BG



AMC-1E-SG



PPE-AA-5



AMC-21E-CW



**Stores Under  
Counters**

## Aluminum Six-Drawer 76cm Electronic Touchpad Locking Narrow Cart

**AMC-3E-Color Six-Drawer 76cm Electronic Touchpad Locking Narrow Cart**

*Beige shell with eleven drawer colors available (see page 29).*

### Features:

- Aluminum construction; weighs only 38.5kg
- Stabilizing frame w/ bumper
- Plastic top
- Four 12.7cm swivel casters (two with brake and one tracking)
- Electronic Touchpad locking with key override

### Dimensions:

**Overall:** 111.3cm High x 64.1cm Deep x 62.7cm Wide (includes bumper and casters).

**Drawer configuration:** 76cm of total drawer height.

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (3), 15.6cm (2), 23.8cm (1).

## Aluminum Five-Drawer 61cm Electronic Touchpad Locking Cart

**PPE-Color-24 Five-Drawer 61cm Electronic Touchpad Locking Cart** *Beige shell with eighteen drawer colors available.*

**Drawer sizes:** 7.3cm (3), 15.6cm (1), 23.8cm (1).

**Weight:** 35.4kg

Four 12.7cm swivel casters (two with brake and one tracking)

*See page 27 for more details.*

## Steel Five-Drawer 61cm Electronic Touchpad Locking Mini Cart

**AMC-1E-Color Five-Drawer 61cm Electronic Touchpad Locking Mini Cart with Key Override**

*Beige shell with eleven drawer colors available (see page 28).*

**Overall:** 89.1cm High x 52.4cm Deep x 54.2cm Wide.

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (3), 15.6cm (1), 23.8cm (1).

**Weight:** 46.7kg

Four 7.6cm swivel casters (two with brake and one tracking)

## Steel Four-Drawer 53cm Electronic Touchpad Locking Mini Cart

**AMC-21E-Color Four-Drawer 53cm Electronic Touchpad Locking Mini Cart with Key Override**

*Beige shell with eleven drawer colors available (see page 28).*

**Overall:** 81.6cm High x 52.3cm Deep x 54.2cm Wide.

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (2), 15.6cm (1), 23.8cm (1).

**Weight:** 38.6kg

Four 7.3cm swivel casters (two with brake and one tracking)


# Standard Steel Supreme Anesthesia Cart Package

ELECTRONIC TOUCHPAD

## STANDARD SUPREME ANESTHESIA PACKAGE

Complete Steel Cart with entire Accessory Package as shown

ATE-AA-6 with ASA-1X Package




NOTE: Carts require some assembly.

# Premier™ Aluminum Supreme Anesthesia Cart Package

ELECTRONIC TOUCHPAD  
Not MRI Compatible

## PREMIER™ SUPREME ANESTHESIA PACKAGE

Complete Aluminum Cart with entire Accessory Package as shown  
PKL-AA-30 with PASA-1X Package



Look for the  
Armstrong A-SMART® label  
on every cart.  
*It stands for quality.*  
**MADE IN THE USA**

NOTE: Carts require some assembly.



## These Accessory Packages Mount To Any Standard Steel Cart

### Deluxe Anesthesia Accessory Package

DAB-1 Deluxe Accessory Package (*Does not include cart*)

Complete Package as shown, DAB-1 and ATE-AA-6 Six-Drawer 76cm Electronic Touchpad Locking Cart

- ATR-6N Trellis System
- ATB-4X Four-bin Tilt Bin Organizer
- ATB-5X Five-bin Tilt Bin Organizer
- PSR-1 Tubular Railing
- ASN-3 Collapsible Side Shelf
- AWC-1 Waste Container
- AAB-1 Accessory Bar
- APT-2 Plastic Top



PIV-1 I.V. Pole with  
Brackets, *see page 35*



Look for the Armstrong  
A-SMART® label on every  
cart. *It stands for quality.*  
**MADE IN THE USA**



ATE-AA-6

### *Electronic Touchpad Carts, see page 27*



ATE-B-6

### Economy Anesthesia Accessory Package

SP-1B Economy Accessory Package (*Does not include cart*)

Complete Package as shown,  
SP-1B and ATE-B-6 Six-Drawer  
76cm Electronic Touchpad Locking Cart

- APT-2 Plastic Top
- AUH-1 Utility Hooks, 2/pkg.
- PSR-1 Tubular Railing
- Two AC-672S 12.7cm Drug Trays with four Short Dividers and two Long Dividers (fills one drawer)
- Four AC-670S 7.3cm Drug Trays with eight Short Dividers and four Long Dividers (fills two drawers)



Modular Drawer Trays,  
*see page 38*

*NOTE: Carts require some assembly.*

# Accessory Packages To Lighten The Load. . .

**Think** *LIGHT.*

## Deluxe Anesthesia Accessory Package

PDAB-1 Premier Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, PDAB-1 and PKL-SB-30 Six-Drawer 76cm All Aluminum Mobile Workstation

- PNS-1 Shelf Unit
- ATB-5X Five-bin Tilt Bin Organizer
- AGB-1 Glove Box
- PSB-1 Sharps Container Bracket
- AWC-1 Waste Container



**ScopeSafe™**  
*Protect Your  
Delicate  
Scopes*

ScopeSafe™  
Locking Scope  
Holder,  
*see page 30*



PKL-SB-30

## Economy Anesthesia Accessory Package

PSP-1B Premier Economy Accessory Package  
*(Does not include cart)*

Complete Package as shown, PSP-1B and PKL-CW-30 Six-Drawer 76cm All Aluminum Mobile Workstation

- AGB-1 Glove Box
- AWC-1 Waste Container
- Four AC-670S 7.3 cm Drug Trays with eight Short Dividers and four Long Dividers
- Two AC-672S 12.7 cm Drug Trays with four Short Dividers and two Long Dividers



PKL-CW-30



Premier™  
Locking  
Bars,  
*see page 26*

*See page 5 for Premier™ Cart features and page 7 for dimensions.*

# Standard Steel Solid Color Mobile Workstations

KEY-LOCKING

## Solid Colors, Solid Features

All A-SMART® Key Locking Carts come with double side-wall construction, ball bearing drawer slides, a standard stabilizing frame/bumper, and four all-swivel casters (two locking and one tracking).



### Six-Drawer 76cm Standard Mobile Workstations

AA-6 Mobile Workstation,  
Deep Blue

Drawer sizes: 7.3cm (3), 15.6cm (2),  
23.8cm



### Six-Drawer 76cm Standard Mobile Workstations

AN-6 Mobile Workstation,  
Beige

Drawer sizes: 7.3cm (3), 15.6cm (2),  
23.8cm



### Five-Drawer 61cm Standard Mobile Workstations

AN-5 Mobile Workstation,  
Beige

Drawer sizes: 7.3cm (3), 15.6cm,  
23.8cm



Look for the  
Armstrong A-SMART® label  
on every cart.  
*It stands for quality.*

**Want more accessories?** See pages 32-Inside Back Cover



PIV-1 I.V. Pole with  
Brackets, see page 35



ATR-6N Standard Trellis  
Rail System, see page 26



ASN-3 Side Shelf,  
see page 35



APT-2 Standard Plastic  
Top, see page 35

Standard Anesthesia Packages, pages 18, 20  
Cart Accessories, pages 32-Inside Back Cover

See page 4 for Standard Cart features and  
page 6 for dimensions.



## Premier™ Carts Are Easier To Push and STOP. . .

**Think** *LIGHT.*



**Six-Drawer 69cm Premier Mobile Workstation**

PKL-B-27 Mobile Workstation,  
Beige

Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm



**Five-Drawer 61cm Premier Mobile Workstation**

PKL-B-24 Mobile Workstation,  
Beige

Drawer sizes: 7.3cm (3), 15.6cm,  
23.8cm



**Six-Drawer 76cm Premier Mobile Workstation**

PKL-B-30 Mobile Workstation,  
Beige

Drawer sizes: 7.3cm (3), 15.6cm (2),  
23.8cm



AEF-6 Plastic Divider Set  
for 6" drawer, see page 38



PDS-1 One-shelf Unit, Deep,  
see page 36


Visit Us At:

[www.armstrongmedical.11](http://www.armstrongmedical.11)



Look for the  
Armstrong A-SMART® label  
on every cart.  
***It stands for quality.***

### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

# Standard Steel Two-Tone Mobile Workstations

KEY-LOCKING

## Colors To Match Your Decor

*When your colors are added to the quality construction of an Armstrong A-SMART® Cart, you've got a Mobile Workstation like no other. **A-SMART®***



### Four-Drawer 76cm Standard Mobile Workstation

AKL-BG-4 Mobile Workstation,  
Blue Grey  
Drawer sizes: 7.3cm, 23.8cm (3)



### Six-Drawer 76cm Standard Mobile Workstation

AKL-SB-6 Mobile Workstation,  
Slate Blue  
Drawer sizes: 7.3cm (3),  
15.6cm (2), 23.8cm



### Five-Drawer 61cm Standard Mobile Workstation

AKL-MG-5 Mobile Workstation,  
Moss Gray  
Drawer sizes: 7.3cm (3), 15.6cm,  
23.8cm



Look for the Armstrong A-SMART® label  
on every cart. ***It stands for quality.***  
**MADE IN THE USA**



Standard Side Shelf,  
see page 35



Modular Drawer Trays,  
see page 38



Wire Baskets, see page 34

## 18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain


Standard Anesthesia Packages, pages 18, 20

Accessories, pages 32-Inside Back Cover

# Maneuverability At It's Best. . .

**Think** *LIGHT.*

## 18 Designer Colors



PKL-CB-27



PKL-HG-30



PKL-Y-4

### Six-Drawer 69cm Premier Mobile Workstation

PKL-CB-27 Mobile Workstation, Taupe

Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm

### Six-Drawer 76cm Premier Mobile Workstation

PKL-HG-30 Mobile Workstation, Hunter Green

Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm

### Four-Drawer 76cm Premier Mobile Workstation

PKL-Y-4 Mobile Workstation, Yellow

Drawer sizes: 7.3cm, 23.8cm (3)



Look for the  
Armstrong A-SMART® label  
on every cart.

*It stands for quality.*

**MADE IN THE USA**

[www.armstrongmedical.com](http://www.armstrongmedical.com)

### Five-Drawer 61cm Premier Mobile Workstation

PKL-SB-24 Mobile Workstation, Slate Blue

Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm



PKL-SB-24

### Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster
- Tracking caster
- Key locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks





## Restock Only The Accessed Drawers. Save Money. Save Time.

Don't waste money checking crash cart drawers that haven't even been opened. With A-SMART® Premier™ Carts, individual drawer seals can be arranged so that only the accessed drawers need to be checked for inventory, a small investment that pays off with every use.

- Double lock medication drawers
- Improve equipment security
- Ready for changing needs

When carts are shared by departments or personnel, A-SMART Premier's locking options allow each to be locked with a separate key. Never wonder about who else has access to your equipment.

Maybe you don't need this kind of flexibility now, but you might later. With A-SMART Premier, you can update your locking system at any time. Need to change a mobile workstation into a crash cart? No problem. Need to add a key lock to an existing cart? Done. Only A-SMART Premier offers you that kind of flexibility.



*Locks may be placed on either corner glide for flexible security. Double locking drawers can help control inventory and access.*



### Key Locking Bars

*Different sized locking bars help control access.*


Size	Prod#
7.3cm	PKE-3
15.6cm	PKE-6
23.8cm	PKE-9
61cm	PKE-24
69cm	PKE-27
76cm	PKE-30



### Breakaway Locking Bars

*Mix lock sizes and types for precise inventory control.*

Size	Prod#
7.3cm	PBR-3
15.6cm	PBR-6
23.8cm	PBR-9
38.1cm	PBR-15
53.3cm	PBR-21
61cm	PBR-24
69cm	PBR-27
76cm	PBR-30



#### PTR-1 PREMIER™ Trellis System

Adds up to 53.3cm to cart height. Bars are height-adjustable. Includes two uprights and three cross bars.

#### PTR-5 PREMIER™ Trellis System with Hooks (not pictured)

Adds up to 53.3cm to cart height. Bars are height-adjustable. Includes two uprights and three cross bars. Includes one bar with five hooks in the top rail.



#### ATR-5N STANDARD Trellis Rail System with Five Hooks

Same features as the ATR-6N (see right), but with five hooks on the top rail for hanging supplies.

#### PHB-1 Additional Hanging Bar for Standard or Premier Trellis Systems



#### ATR-6N STANDARD Trellis Rail System

The top bar is 54.6cm above top of cart.


#### ATR-9N Wide Cart Trellis Rail System

Fits Armstrong's Wide AKW Series Carts (see page 30).

# NEW Aluminum Electronic Touchpad Touchpad Carts

ALUMINUM AND STEEL ELECTRONIC TOUCHPAD  
Not MRI Compatible

## 18 Designer Colors



## NEW Aluminum ELECTRONIC Touchpad Cart

### PPE-Color-30 76cm Six-Drawer Aluminum Electronic Touchpad Cart

Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm  
Weight: 41.7 kg.

### PPE-Color-24 61cm Five-Drawer Aluminum Electronic Touchpad Cart

Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm.  
Weight: 34.4 kg.

#### Features:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Plastic top
- Accessory panels (3)
- Stabilizing frame with bumper
- Aluminum mounting tracks
- 18 colors available

**NOTE:** No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.

*\* Please note for the Armstrong Electronic Touchpad Carts, we suggest using alcohol only for the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.*



PPE-AA-30

PPE-AA-30 shown with optional PSR-1 Railing, ATD-2 Tape Dispenser, PTR-1 Trellis System, and ATB-4X & ATB-5X Tilt Bins



ATE-CB-6

## Standard Steel Electronic Touchpad Locking Carts

### ATE-Color-6 76cm Standard Six-Drawer Electronic Touchpad Locking Cart

Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm.  
Weight: 77kg

### ATE-Color-5 61cm Standard Five-Drawer Electronic Touchpad Locking Cart

Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm.  
Weight: 68kg

This system comes with a key override to bypass Electronic Touchpad System if needed.

18 colors available.

ATE-CB-6 shown at left with optional APT-2 Plastic Top and PSR-1 Tubular Railing.



AMC-1SG-E

**\*Optional 12.7cm Casters with Stabilizing Frame**



AMC-1B

# Spacesaver

## Steel Two-Tone Mini Carts

Ideal For Treatments, Procedures, & the Bedside. *Small enough to put under counters.*

### Steel Electronic Touchpad Locking Mini Cart

AMC-1E-Color Five-Drawer 61cm Electronic Touchpad Locking Mini Cart

Compact, five-drawer beige shell Mini Carts include a plastic top, ball bearing drawer slides, corner bumper guards, and 3" swivel casters\* (two locking, one tracking). Mini Carts also features colored drawer fronts in a beige frame.

**Overall:** 89.1cm High x 52.4cm Deep x 54.2cm Wide.

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (3), 15.6cm, 23.8cm

**Weight:** 46.7kg

#### Features:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Plastic top
- 11 colors available

**NOTE: No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.**

### Steel Key Locking Mini Cart

AMC-1Color 61cm Five-Drawer Key Locking Mini Cart

**Overall:** 88.9cm High x 52.1cm Deep x 53.3cm Wide.

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (3), 15.6cm, 23.8cm

**Weight:** 46.7kg

### Steel 53cm Electronic Touchpad Locking Mini Cart

AMC-21E-Color Four-Drawer 53cm Electronic Touchpad Locking Mini Cart

**Overall:** 81.6cm High x 52.4cm Deep x 54.2cm Wide.



AMC-21E-CW

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (2), 15.6cm, 23.8cm.

**Weight:** 38.5kg

#### 11 Designer Colors



SG - Sage Green



MG - Moss Grey



TC - Terra Cotta



AA - Deep Blue



B - Beige



BG - Blue Grey



HG - Hunter Green



T - Teal



Y - Yellow



CB - Taupe




CW - Cherry Woodgrain



# Armstrong Aluminum Narrow Carts & Mini/Narrow Cart Accessories

ELECTRONIC TOUCHPAD LOCKING & KEY LOCKING  
Not MRI Compatible

## 11 Designer Colors



## Great For All Those Tight Spaces

### Features:

- Lightweight aluminum construction; weighs only 38.6kg
- Stabilizing frame with bumper
- Four 12.7cm swivel casters (two with brake and one tracking)
- Plastic top
- Eleven colors available

### Dimensions:

**Overall:** 111.8cm High x 64cm Deep x 62.7cm Wide (includes bumper and casters).

**Drawer configuration:** 76cm of total drawer height.

**Drawer size:** 43.2cm Deep x 38.1cm Wide.

**Drawer heights:** 7.3cm (3), 15.6cm (2), 23.8cm (1).

See below for Narrow Cart accessories.

**Glove  
Dispensers,  
Inside Back Cover**



AMC-2HG

### Aluminum 76cm Key Locking Narrow Cart

AMC-2Color Six-Drawer 76cm  
Key Locking Narrow Cart



AMC-3-BG

### Aluminum 76cm Electronic Touchpad Locking Narrow Cart

AMC-3Color Six-Drawer 76cm  
Electronic Touchpad Locking  
Narrow Cart

## Optional **MINI** & **NARROW** Cart Accessories

APH-1 Push Handle

MDM-1B Keyboard Drawer Weight  
capacity: 13.6kg

ATR-7N Three-rail Trellis System

ATR-8N Three-rail Trellis System with  
Four Hooks

AMT-4N Full Drawer Tray with Nine  
Ampule Dividers, Ampule Tray, Rail, and

Four Adjustable Dividers (Fits one per  
drawer), see page 38.

AMT-5N Full Drawer Tray with  
Two Rails and Four Adjustable Dividers  
(Fits one per drawer), see page 38.

AEF-9M Plastic Divider Set for 23.8cm

Mini Cart Drawer

MDM-1B Keyboard Drawer

AEF-6M Plastic Divider Set for 15.6cm

Mini Cart Drawer

AEF-3M Plastic Divider Set for 7.3cm  
Mini Cart Drawer

Also compatible with all accessories  
that fit on the sides of our Standard  
Mobile Workstations. See pages  
32-Inside Back Cover for more accessories.

# Aluminum Wide Cart and ScopeSafe™

## ELECTRONIC TOUCHPAD



AKE-B-6

## NEW Aluminum ELECTRONIC Touchpad WIDE Cart

### AKE-B-6 6-Drawer 76cm Electronic Touchpad Wide Cart, Beige

A-SMART® Wide Carts feature:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Extra-wide drawers: 7" wider than standard cart drawers
- Seamless, one-piece plastic top
- Standard stabilizing frame with bumper
- 5" all-swivel casters (two locking, one tracking)
- Beige drawer fronts in beige frame
- Carts weigh 53.5 kg.

#### Dimensions:

**Overall size:** 116.2cm H x 64.1cm D x 103.4cm W (includes bumper, handles, and casters).

**Actual Top Dimensions:** 82.5cm W x 43.2cm D

#### Available drawer configuration:

**Six-Drawer:** 7.3cm (3), 15.6cm (2), 23.8cm (1).

**NOTE:** No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.

A-SMART WIDE CART SYSTEMS ACCEPT ALL THE ACCESSORIES THAT FIT ON THE SIDES OF OUR STANDARD MOBILE WORKSTATIONS.

*\* Please note for the Armstrong Electronic Touchpad Carts, we suggest using alcohol only for the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.*



ATR-9N

### ATR-9N Wide Cart Trellis Rail System

### ATR-9X Additional Wide Cart Hanging Bar

### ANS-9 Narrow Shelf for Wide Cart

Dim.: 25.4cm x 87.6cm. Adds up to 53.3cm of height to cart.  
Wt. capacity: 23kg



ANS-9

## ScopeSafe™

### Key-Locking Storage for Delicate Scopes

#### ASH-4 ScopeSafe™ Locking Scope Holder

- Holds four scopes
- Has a padded interior for cushioned protection
- Contains multiple vents for ample ventilation
- Contains a protective plexiglas door front
- Attaches easily to wall or cart

- Built with a durable steel frame
  - Overall size: 122cm High x 15.2cm Deep x 40.6cm Wide
- Weight: 20 kg.

Takes two Accessory Panels to mount to a Premier Cart. Not compatible with Mini or Narrow Carts.

#### ASH-2 ScopeSafe™ Locking Scope Holder (not pictured)

- Holds two scopes
  - Same high-quality features as ASH-4
  - Compact: 122cm High x 15.2cm Deep x 22.9cm Wide
- Weight: 14.1 kg.

Takes two Accessory Panels to mount to a Premier Cart. Two ASH-2 ScopeSafes can be mounted to one side of a cart. Not compatible with Mini or Narrow Carts.

**Please note:** ASH-2 and ASH-4 will not hold colonoscopes.



ASH-4

## Isolation Carts



AKL-Y-400IP

### AKL-Y-400IP **STANDARD** Steel Key Locking Four-Drawer Isolation Cart Package

Includes AKL-Y-4 Four-Drawer Cart, PSR-1 Tubular Railing, APT-2 Plastic Top, AUH-1 Utility Hooks, and ACR-1 Coat Rack.

For Standard Cart features, see page 4.

For Standard Cart dimensions, see page 6.



PKL-Y-4

### PKL-Y-4 Basic **PREMIER™** Aluminum Key Locking Four-

Drawer 76cm Isolation Cart

### PBL-Y-4 Basic **PREMIER™** Aluminum Breakaway Locking

Four-Drawer 76cm Isolation Cart

Drawers: 7.3cm, 23.8cm (3)

For Premier™ Cart features, see page 5.

For Premier™ Cart dimensions, see page 7.



AMC-1Y

### AMC-1Y Five-Drawer 61cm Key Locking **STEEL** Mini Cart

Dimensions:

Overall: 88.3cm High x 52.4cm Deep x 54.6cm Wide.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm, 23.8cm.

Weight: 47kg



AMC-2Y

### AMC-2Y Six-Drawer 76cm Key Locking **ALUMINUM** Narrow Cart

Dimensions:

Overall: 111.8cm High x 62.9cm Deep x 62.2cm Wide (includes bumper & casters).

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm (2), 23.8cm (1).

Weight: 38.5kg



AKL-Y-4

### AKL-Y-4 Basic **STANDARD** Steel Key Locking Four-Drawer 76cm Isolation Cart

Drawers: 7.3cm, 23.8cm (3)

For Standard Cart features, see page 4.

For Standard Cart dimensions, see page 6.



ABL-Y-3

### ABL-Y-3 Basic **STANDARD** Steel Breakaway Locking Three-Drawer 61cm Isolation Cart

Drawers: 15.6cm, 23.8cm (2)

For Standard Cart features, see page 4.

For Standard Cart dimensions, see page 6.



## Cart Accessories



**ACG-1 Caster Guard, 4/set**

For 12.7cm dia. wheels.  
Durable with living hinge,  
locking tab, closure  
positioning/support pins.



**ATC-3 Tracking Swivel Caster**



**APH-1 STANDARD Push Handle**

Does not fit A-SMART  
Premier Carts.



**ATW-3 Soft-grip Twin Push Handles**



**ACC-24 Navy Cart Cover for 61cm A-SMART® Cart**  
**ACC-27 Navy Cart Cover for 69cm A-SMART Cart**  
**ACC-30 Navy Cart Cover for 76cm A-SMART Cart**

Nylon covers with Velcro tabs.



**AQC-1 Silent Knight Tablet Crusher**  
Includes fifty pouches.  
**AQC-2 Silent Knight Pouch, 1000/box**



**APS-5 Plastic Seals, 100/pkg.**  
Consecutively numbered.



**ADL-1 Drawer Liner**  
Black vinyl and latex free.  
43.2cm x 55.9cm. Fits  
standard size drawers.



**ALH-1 Label Holders and Cards (blank), 100/pkg.** 7.6cm x 2.5cm  
**ALH-2 Label Holders, 100/pkg.** 17.8cm x 5.1cm  
**ALH-12 Label Holders, 12/pkg.** 17.8cm x 5.1cm



**APS-1 Pull-tite Seals, 100/pkg.**  
25.1cm polypropylene seals  
open without using scissors.  
Consecutively numbered.



**AUH-1P PREMIER Utility Hooks, 2/pkg.**  
Mount in corner tracks.



**AUH-1 STANDARD Utility Hooks, 2/pkg.**  
Mount to sides or back  
of cart.

*All products are compatible with both Standard and Premier™ Carts unless otherwise stated.*

*To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.*



ACS-2

## \*ACS-2 Short Catheter Holder

- Has metal construction with neutral beige finish
- Requires no screws to mount
- 48.3cm High x 20.3cm Wide x 10.2cm Deep
- Divided into two compartments with non-removable divider

## \*ACS-3 Short Catheter Holder with Clear Front

- Clear front enables you to see inventory
- Front panel removes for easy cleaning
- Requires no screws to mount
- 48.3cm High x 10.2cm Wide x 10.2cm Deep



ACS-3



AOB-1

## \*ACL-2 Long Catheter Holder

- Fits A-SMART® Carts with 76cm drawer space
- Has metal construction with neutral beige finish
- Requires no screws to mount
- 71.1cm High x 20.3cm Wide x 10.2cm Deep
- Divided into two compartments with non-removable divider

## \*ACL-3 Long Catheter Holder with Clear Front

- Fits A-SMART Carts with 76cm drawer space
- Clear front enables you to see inventory
- Front panel removes for easy cleaning
- Requires no screws to mount
- 71.1cm High x 10.2cm Wide x 10.2cm Deep




ACL-2

## ADO-1 Divided Organizer with Lid

Dimensions: 8.3cm High x 18.1cm Wide x 18.1cm Deep.

## \*AOB-1 Oxygen Tank Bracket

Fits "D" and "E" sized tanks.



ADO-1

## AFT-1 **STANDARD** Fluid Tray System for Cart with 61cm Drawer Space, 3 trays/pkg.

## AFT-1P **PREMIER™** Fluid Tray System (all Cart sizes), 3 trays/pkg. (not pictured)

## AFT-2 **STANDARD** Fluid Tray System for Cart with 76cm Drawer Space, 3 trays/pkg.

Three height-adjustable trays on rail system. Inner tray dimensions: 43.5cm Long x 11.4cm Deep x 10.2cm High.



AFT-1



PSB-1

## \*PSB-1 **PREMIER™** Sharps Container Brackets

Fits Sharps containers that are 4.7 liters., 7.6 liters, or 11.3 liters

## ASB-1N **STANDARD** Sharps Container Bracket

Fits Sharps containers that are 4.7 liters., 7.6 liters, or 11.3 liters

*\*May require PHP-1 panels to mount to a Premier™ Cart.*



ASB-1N

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.



ASB-2

## ASB-2 Armstrong Security Box

Lockable box for portable use or inside a 7.3cm drawer. Designed to fit beside an AC-670S Drug Tray. The security box measures 27.9cm Wide x 40.6cm Deep x 6.8cm High.

## ATD-2 Tape Dispenser for Trellis Rail

## ATD-3 Tape Dispenser for Cart Top

5 tape brackets, each bracket width is 5.7cm (*not pictured*)

## AND-1 Armstrong Lockable Drawer Cover

Strong, secure cover can be fitted to any size A-SMART® Cart drawer (7.3cm, 15.6cm, 23.8cm, or 31.8cm drawer).

## AWB-1 Wire Basket\*

Plastic-coated wire basket. 38.1cm Long x 30.5cm Wide. 10.2cm High in back, 7.6cm High in front. Sloping top, flat bottom.

## AWB-3 Wire Chart Holder\*

Plastic-coated wire chart holder. 40.6cm Long x 10.2cm Deep x 25.4cm High.

## AWB-5 Sloping Wire Basket\*

Plastic-coated wire basket measures 30.5cm square. 20.3cm High in back, 10.2cm High in front. Sloping top, flat bottom.

## AWC-1 Waste Container\*

12.9 liters. 28.9cm x 21cm x 30.8cm.

## AWC-2 Fire Retardant Waste Container\*

13.2 liters. 31.1cm x 21cm x 29.2cm.

## AWC-3 Large Waste Container\*

26.5 liters. 36.5cm x 26cm x 38.1cm.

## ACV-1EL Electronic Push-Button ClearView Security Drug Box with Key Override (*manual locking*)

## CVS-1 ClearView Shelf

## PHP-1 **PREMIER™** Accessory Panels, 2/pkg.

40.6cm x 10.2cm.

## PHP-5 **PREMIER™** Accessory Panels for Back of Cart, 2/pkg.

58.4cm x 10.2cm.

## AAB-1 **STANDARD** Accessory Bar Bracket

Required for hanging some accessories on this page. Attaches to two locations on either side of cart. Size: 45.7cm x 2.5cm. (*For A-SMART Premier Carts, substitute **PHP-1** Accessory Panel*).


## PHP-1S **PREMIER™** Panel for Cart-mounted *new-DUET* (*not pictured*)

## AHP-3 **STANDARD** Panel for Cart-mounted *new-DUET* (*not pictured*)

\*Requires AAB-1 or PHP-1 to hang (*see above*)



AND-1



ATD-2



AWB-3



ACV-1EL with CVS-1



PHP-1



AWB-1



AWB-5



AWC-1



AAB-1

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see this page.





PSR-1


## PSR-1 Tubular Railing

Made of electro-polished stainless steel. Extends 8.3cm above the top of the cart.

## APT-2 **STANDARD** Plastic Top

Seamless ABS plastic protects cart metal and provides easy-to-clean work surface. Rounded corners mean infection has nowhere to hide.

Dimensions: 47cm x 68.6cm.



APT-2

## APT-5 **PREMIER™** Plastic Top (not pictured)

Dimensions: 47cm x 68.6cm.

## PIV-1 I.V. Pole with Brackets

Adjustable: extends 24.1cm to 100.3cm above cart top.



PIV-1

## ASN-3 **STANDARD** Side Shelf

### *Simple One-Hand Slow Shelf Retraction*

Adds 1.8 square feet of temporary work surface when extended and folds neatly away when not in use.

Dimensions: 45.1 cm Deep x 37.8 cm Width.

Weight Limit: 11.4 kg.



ASN-3

## PSN-3 **PREMIER™** Side Shelf,

### *Simple One-Hand Slow Shelf Retraction*

Same as ASN-3 (see above).



ACB-1

## ACB-1 **STANDARD** Plastic Cardiac Board with Brackets

Heavy-duty plastic board with handle that mounts to back of cart for easy access.

Dimensions: 39.4cm x 54.6cm x 1.9cm.



PCB-1

## PCB-1 **PREMIER™** Cardiac Board with Brackets

Same as ACB-1 (see above).

## ASU-1 Suction Unit Holder

Dimensions: 25.4cm High x 48.3cm Wide x 19.7cm Deep.



PSN-3



ASU-1

## AGB-1 Glove Box

Mounts to twenty different positions on a Standard Cart. Fasteners included. Beige. Gloves not included.

Dimensions: 25.4cm Long x 14.6cm Wide x 9.5cm Deep.



APC-1

## APC-1 Armstrong Wire Punch Card Tray, White

Will hold up to thirty punch cards per tray.

Dimensions: 17.8cm High x 35.6cm Deep x 17.8cm Wide. Great for fluid transfer.



AGB-1

*All products are compatible with both Standard and Premier™ Carts unless otherwise stated.*

*To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.*

# Shelving Units



ACD-1N PCD-1

**ACD-1N STANDARD** Defib Corner Shelf Unit  
Same features as PCD-1 (*see below*).

**PCD-1 PREMIER™** Defib Corner Shelf Unit  
Adds 24.1cm to height of cart. Swivel top rotates 360° – locks every 45°. Adjustable width: 33-46cm. Depth: 30.5cm. Strap included. Capacity: 11.3kg

**PNT-2 PREMIER™** Two-shelf Unit, Narrow, (*not pictured*)  
Dim.: 26cm x 70cm. Adds up to 53.3cm of height. Units include two uprights and two shelves.



PDS-1

**PNS-1 PREMIER™** One-shelf Unit, Narrow  
Dim.: 26cm x 70cm. Adds up to 53.3cm of height. Capacity: 25kg

**PDS-1 PREMIER™** One-shelf Unit, Deep  
Dim.: 37.1cm x 70cm. Adds up to 53.3cm of height. Capacity: 25kg

**PDT-2 PREMIER™** Two-shelf Unit, Deep (*not pictured*)  
Dim.: 37.1cm x 70cm. Adds up to 53.3cm of height. Units include two uprights and two shelves. Capacity: 25kg.



PNS-1

**ABS-1N STANDARD** One-shelf Unit, Deep  
Adds 50.8cm to height of cart. Outer dimensions: 37.1cm Deep x 70cm Wide. Capacity: 34kg

**ABS-2N STANDARD** One-shelf Unit, Narrow  
Dim.: 26cm x 70cm. Adds up to 50.8cm of height. Capacity: 25kg

**ATN-1N STANDARD** Two-shelf Unit, Deep  
Outer dimensions: 37.1cm Deep x 69.7cm Wide. Adds 58.4cm to height of cart. Capacity: 25kg

**ADS-1N STANDARD** Armstrong Defib One-shelf Unit with Swivel Top  
Outer dimensions: 35.8cm x 68.3cm. Adds 30.5cm to height of cart. Swivel top can rotate 360°. Includes strap to hold monitor. Capacity: 22.7kg

**PDS-2 PREMIER™** Armstrong Defib One-shelf Unit with Swivel Top  
Same as ADS-1N (*see above*).



ABS-2N



ADS-1N PDS-2

## Cart Trellis Systems

*See pages 26*

*All products are compatible with both Standard & Premier Carts unless otherwise stated.*

*To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.*

## Tilt Bins Organizers

Designed for modular assembly, Tilt Bins organize supplies to match your needs. Clear Tilt Bins leave stock visible for quick identification, yet are securely arranged.

Each bin has a slot for labeling and swings open for easy access. Bins are removable for cleaning or refilling.

**NOTE:** A Trellis Rail System or Shelving Unit (except ACD-1N, PCD-1, PLS-1, and ALS-1) is required to mount Tilt Bins to an A-SMART Cart. See page 26 for Trellis Rail Systems and page 36 for Shelving Units.

### ATB-6X Six-Bin Unit

Each bin measures: 11.4cm High x 13.3cm Deep

### ATB-5X Five-Bin Unit

Each bin measures: 16.5cm High x 13.3cm Deep

### ATB-4X Four-Bin Unit


Each bin measures: 20.6cm High x 16.8cm Deep

### ATB-3X Three-Bin Unit

Each bin measures: 24.1cm High x 19.7cm Deep

Each Tilt Bin Organizer has a total width of 57.2cm.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.



A combination of bin sizes can be used to accommodate multiple-sized supplies or equipment.

*Need More Organization  
For All Your Supplies?  
See pages 38-39.*



## Lockable Tilt Bins Organizers



**NOTE:** A Trellis Rail System or Shelving Unit (except ACD-1N, PCD-1, PLS-1, and ALS-1) is required to mount Tilt Bins to an A-SMART Cart. See page 26 for Trellis Rail Systems and page 36 for Shelving Units. Accessories (including Lockable Tilt Bins) can't be mounted directly above Lockable Tilt Bins.

*Tilt Bin dimensions shown above.*

### Key Lock Tilt Bins

ATB-5XK Lockable Five-Bin Unit

ATB-4XK Lockable Four-Bin Unit

ATB-3XK Lockable Three-Bin Unit

### Breakaway Lock Tilt Bins

ATB-5XL Lockable Five-Bin Unit

ATB-4XL Lockable Four-Bin Unit

ATB-3XL Lockable Three-Bin Unit



## A Place For Everything...



AMT-1N

### AMT-1N Modular Tray

Full drawer tray includes two vertical rails and eight adjustable dividers with card holder. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High.



AMT-6N

### AMT-2N Modular Tray

Same as the AMT-1 (see above), except with an ampule organizer, eleven adjustable dividers, and attached label holders. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High



AMT-2N

### AMT-3N Full Drawer Tray with Three Ampule Holders


Thirty-three dividers, three ampule trays, and four securing clips. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High



AMT-7N

### AMT-6N Modular Drawer Tray

Durable, plastic modular tray fits in a 15.6cm drawer. Includes two vertical rails and six dividers. Dimensions: 54.6cm Wide x 38.7cm Deep x 12.7cm High.



AMT-3N

### AMT-7N Modular Tray

Full drawer tray, eighteen ampule dividers, and seven long dividers. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High



AMT-8N

### AMT-8N Modular Tray


Full drawer tray, twenty-nine ampule dividers, and four long dividers. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

## Mini/Narrow Cart Trays

### AMT-4N Full Drawer Tray for Mini Cart

Full drawer tray with nine ampule dividers, ampule tray, rail, and four adjustable dividers. Fits one per drawer. Dimensions: 36.1cm Wide x 40.6cm Deep x 6.4cm High



AMT-4N

### AMT-5N Full Drawer Tray for Mini Cart

Full drawer tray with two rails and four adjustable dividers. Fits one per drawer. Dimensions: 36.1cm Wide x 40.6cm Deep x 6.4cm High



AMT-5N

## ...Everything In Its Place.

### ADV-1N Modular Drawer Divider Kit

Subdivides a 7.3cm drawer to better organize your supplies. Includes two 55.9cm rails, five 39.4cm rails, ten 16.5cm dividers, and three "L" clips.



**IMMEDIATE**, accurate access

**FASTER** emergency response

**EASIER**, better inventory control



### ADT-1 Crash Cart Drawer Tray with Lid

#### Control Your Drawer Inventory

Clear acrylic lid with thumb hole can be secured with locking seal.

Inside dim.: 53.3cm Wide x 38.1cm

Deep x 6.3cm High.


Outside dim.: 55.9cm Wide x 40.6cm

Deep x 7cm High.



### ADV-10N Modular Divider Kit for ADT-1 Drawer Tray

Fits inside ADT-1 Drawer Tray. Includes two 52.5cm Long grooved plastic rails, four 35.6cm Long rails, six 10.1cm Long dividers, and six 14cm Long dividers.



### ADV-6N Modular Drawer Divider Kit for 15.6cm, 23.8cm, and 30.5cm Drawers

12.7cm high modular kit fits 15.6cm, 23.8cm, or 30.5cm drawers. Includes two 55.9cm rails, four 39.4cm rails, and six 16.8cm dividers.



### APS-1 Pull Tite™ Seal, 100/pkg.

#### Easy-to-use Drawer Inventory Control

25.1cm polypropylene seals fit the ADT-1 Drawer Tray. Consecutively numbered.



### ADB-1 Security Bag, 20/pkg.

67.3cm x 47cm. Fits over 7.3cm trays (ADT-1, AMT-1, AMT-2, and AMT-3 Trays). Perforated center shows you when it's been used.

### ADB-2 Security Bag, 20/pkg.

82.5cm x 54.6cm. Fits over 12.7cm tray.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.



# Drawer Trays for Stock Exchange



## Drug Trays

### AC-670S Divisible Drug Tray, 7.3cm Deep

Sand color. Inside dimensions: 37.8cm Long x 23.5cm Wide. Holds up to eleven short dividers and seven long dividers (*not included*). Drawers hold a maximum of two trays, except for the Mini and Narrow Carts.

**AC-675S Short Divider for AC-670S**  
Max. slots: 11

**AC-680S Long Divider for AC-670S**  
Max. slots: 7

### AC-672S Divisible Drug Tray, 12.7cm Deep

Same as the AC-670S (*see above*), but made for deeper drawers. Drawers hold a maximum of two trays.

**AC-685S Short Divider for AC-672S**  
Max. slots: 11

**AC-690S Long Divider for AC-672S**  
Max. slots: 7

### ACC-655S Clear Acrylic Tray Cover, 2/pkg.

See-through cover fits AC-670S and AC-672S Divisible Drug Trays. Thumb notch ensures easy access.

### AC-665 Drug Tray, Small, 6.4cm Deep (*not pictured*)

Blue color. Inside dimensions: 23.5cm Long x 16.8cm Wide. Drawers hold a maximum of four trays.

**AC-675 Long Divider for AC-665**  
Max. slots: 5

**AC-677 Short Divider for AC-665**  
Max. slots: 7



## Armstrong Security Box

### ASB-2 Armstrong Security Box

#### Double-Lock Security

This Lockable Security Drug Box can be portable or fastened to the inside of a 7.3cm drawer. The Security Box is


designed to fit two per drawer, or next to our AC-670S 7.3cm Drug Tray.

Dimensions: 27.9cm Wide x 40.6cm Deep x 6.7cm High.

## Plastic Divider Sets

### AEF-3 Plastic Divider Set for 7.3cm Drawer

Two 55.9cm Long rails and six 41.3cm Long x 6.7cm High dividers.



### AEF-3L Long Plastic Divider Set for 7.3cm Drawer

Four 55.9cm Long rails and two 41.3cm Long x 6.7cm High dividers.

### AEF-6 Plastic Divider Set for 15.6cm Drawer

Two 55.9cm Long rails and six 41.3cm Long x 14.3cm High dividers.

### AEF-6L Long Plastic Divider Set for 15.6cm Drawer

Four 55.9cm Long rails and two 41.3cm Long x 14.3cm High dividers.

### AEF-9 Plastic Divider Set for 23.8cm Drawer

Two 55.9cm Long rails and six 41.3cm Long x 22.2cm High dividers.



### AEF-9L Long Plastic Divider Set for 23.8cm Drawer

Four 55.9cm Long rails and two 41.3cm Long x 22.2cm High dividers.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.



## Cart Extender

### ACE-1 Armstrong Cart Extender\*, Front-open

- Adds 1.5 cu. meters to cart
- Key lock to protect contents
- Includes three adjustable shelves
- Each shelf fits AC-670S or AC-672S drawer trays (one per shelf)
- Safety foot adds stability
- Dimensions: 87.3cm High x 45.7cm Deep x 34.3cm Wide

*\*Cart Extenders only fit A-Smart Carts with 76cm of drawer space. Does not fit on a Narrow or Mini Cart.*



## Looking For Extra Wide Drawers?

See page 30, Aluminum Electronic Touchpad WIDE Cart

## Glove Dispensers



AGB-1

### AGB-1 Cart-mount Glove Box

14.6cm Wide x 25.4cm High x 9.8cm Deep.

### AE-5359 Glove Dispenser, 6/case

Wall-mountable glove dispenser in a light beige. Adjustable shelf fits different sized boxes of gloves.

36.8cm x 15.2cm x 11.4cm.

### AGB-2 Steel, Baked-on Enamel Glove Box

Holds 2 boxes of gloves.

25.4cm Wide x 27.3cm High x 9.8cm Deep.

*Mounting hardware not included*

### AGB-3 Steel, Baked-on Enamel Glove Box

40.6cm Wide x 25.4cm High x 9.8cm Deep.



AE-5359

### AGB-4 Steel, Baked-on Enamel Glove Box

*(not shown)*

Holds 4 boxes of gloves.

57.2cm Wide x 25.4cm High x 9.8cm Deep.

*Mounting hardware not included*

### AGB-2P Plexiglas Glove Box

25.4cm Wide x 27.3cm High x 9.8cm Deep.

*Mounting hardware not included*

### AGB-3P Plexiglas Glove Box *(not shown)*

Holds 3 boxes of gloves.

40.6cm Wide x 25.4cm High x 9.8cm Deep.

*Mounting hardware not included*



AGB-2



AGB-3



AGB-2P



AGB-3V

### AGB-3V Steel, Baked-on Enamel Glove Box

Holds 3 boxes of gloves vertically. 14cm Wide x 76.2cm High x 10.2cm Deep.

*Mounting hardware not included*